IAEA Analytical Quality in Nuclear Applications Series No. 39

A Procedure for the Rapid Determination of ²²⁶Ra and ²²⁸Ra in Drinking Water by Liquid Scintillation Counting



A PROCEDURE FOR THE RAPID DETERMINATION OF ²²⁶Ra AND ²²⁸Ra IN DRINKING WATER BY LIQUID SCINTILLATION COUNTING The following States are Members of the International Atomic Energy Agency:

AFGHANISTAN ALBANIA ALGERIA ANGOLA ARGENTINA ARMENIA AUSTRALIA AUSTRIA AZERBALIAN BAHAMAS BAHRAIN BANGLADESH BELARUS BELGIUM BELIZE BENIN BOLIVIA BOSNIA AND HERZEGOVINA BOTSWANA BRAZIL **BRUNEI DARUSSALAM BULGARIA BURKINA FASO** BURUNDI CAMBODIA CAMEROON CANADA CENTRAL AFRICAN REPUBLIC CHAD CHILE CHINA COLOMBIA CONGO COSTA RICA CÔTE D'IVOIRE CROATIA CUBA **CYPRUS** CZECH REPUBLIC DEMOCRATIC REPUBLIC OF THE CONGO DENMARK DOMINICA DOMINICAN REPUBLIC ECUADOR EGYPT EL SALVADOR ERITREA **ESTONIA ETHIOPIA** FIJI FINLAND FRANCE GABON GEORGIA GERMANY

GHANA GREECE **GUATEMALA** HAITI HOLY SEE HONDURAS HUNGARY **ICELAND** INDIA **INDONESIA** IRAN, ISLAMIC REPUBLIC OF IRAO IRELAND ISRAEL ITALY JAMAICA JAPAN JORDAN **KAZAKHSTAN KENYA** KOREA, REPUBLIC OF **KUWAIT KYRGYZSTAN** LAO PEOPLE'S DEMOCRATIC REPUBLIC LATVIA LEBANON LESOTHO LIBERIA LIBYA LIECHTENSTEIN LITHUANIA LUXEMBOURG MADAGASCAR MALAWI MALAYSIA MALI MALTA MARSHALL ISLANDS MAURITANIA, ISLAMIC REPUBLIC OF MAURITIUS MEXICO MONACO MONGOLIA MONTENEGRO MOROCCO MOZAMBIQUE MYANMAR NAMIBIA NEPAL NETHERLANDS NEW ZEALAND NICARAGUA NIGER NIGERIA NORWAY

OMAN PAKISTAN PALAU PANAMA PAPUA NEW GUINEA PARAGUAY PERU PHILIPPINES POLAND PORTUGAL OATAR REPUBLIC OF MOLDOVA ROMANIA RUSSIAN FEDERATION RWANDA SAN MARINO SAUDI ARABIA SENEGAL SERBIA SEYCHELLES SIERRA LEONE SINGAPORE SLOVAKIA **SLOVENIA** SOUTH AFRICA SPAIN SRI LANKA **SUDAN SWAZILAND SWEDEN** SWITZERLAND SYRIAN ARAB REPUBLIC TAJIKISTAN THAILAND THE FORMER YUGOSLAV REPUBLIC OF MACEDONIA TOGO TRINIDAD AND TOBAGO TUNISIA TURKEY UGANDA UKRAINE UNITED ARAB EMIRATES UNITED KINGDOM OF GREAT BRITAIN AND NORTHERN IRELAND UNITED REPUBLIC OF TANZANIA UNITED STATES OF AMERICA URUGUAY UZBEKISTAN VENEZUELA, BOLIVARIAN REPUBLIC OF VIET NAM YEMEN ZAMBIA ZIMBABWE

The Agency's Statute was approved on 23 October 1956 by the Conference on the Statute of the IAEA held at United Nations Headquarters, New York; it entered into force on 29 July 1957. The Headquarters of the Agency are situated in Vienna. Its principal objective is "to accelerate and enlarge the contribution of atomic energy to peace, health and prosperity throughout the world".

IAEA Analytical Quality in Nuclear Applications No. 39

A PROCEDURE FOR THE RAPID DETERMINATION OF ²²⁶Ra AND ²²⁸Ra IN DRINKING WATER BY LIQUID SCINTILLATION COUNTING

INTERNATIONAL ATOMIC ENERGY AGENCY VIENNA, 2014

COPYRIGHT NOTICE

All IAEA scientific and technical publications are protected by the terms of the Universal Copyright Convention as adopted in 1952 (Berne) and as revised in 1972 (Paris). The copyright has since been extended by the World Intellectual Property Organization (Geneva) to include electronic and virtual intellectual property. Permission to use whole or parts of texts contained in IAEA publications in printed or electronic form must be obtained and is usually subject to royalty agreements. Proposals for non-commercial reproductions and translations are welcomed and considered on a case-by-case basis. Enquiries should be addressed to the IAEA Publishing Section at:

Marketing and Sales Unit, Publishing Section International Atomic Energy Agency Vienna International Centre PO Box 100 1400 Vienna, Austria fax: +43 1 2600 29302 tel.: +43 1 2600 22417 email: sales.publications@iaea.org http://www.iaea.org/books

For further information on this publication, please contact:

Terrestrial Environment Laboratory, Seibersdorf International Atomic Energy Agency 2444 Seibersdorf Austria Email: Official.Mail@iaea.org

A PROCEDURE FOR THE RAPID DETERMINATION OF 226Ra AND 228Ra IN DRINKING WATER BY LIQUID SCINTILLATION COUNTING

IAEA, VIENNA, 2014 IAEA/AQ/39 ISSN 2074–7659 © IAEA, 2014 Printed by the IAEA in Austria December 2014

FOREWORD

Since 2004, the environment programme of the IAEA has included activities aimed at the development of a set of procedures for the determination of radionuclides in terrestrial environmental samples. Reliable, comparable and 'fit for purpose' results are essential requirements for any decision based on analytical measurements. For the analyst, tested and validated analytical procedures are extremely important tools for the production of such data. For maximum utility, such procedures should be comprehensive, clearly formulated and readily available to both the analyst and the customer.

This publication describes a procedure for the rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water. The determination of radium in drinking water is important for protecting human health, since the consumption of drinking water containing radium may lead to an accumulation in the body, contributing to the radiological dose. The method is based on the separation of ²²⁶Ra and ²²⁸Ra from interfering elements using PbSO₄ and Ba(Ra)SO₄ co-precipitation steps. The isotopes ²²⁶Ra and ²²⁸Ra are then determined by liquid scintillation counting. The procedure is expected to be of general use to a wide range of laboratories, including the laboratories of the Analytical Laboratories for the Measurement of Environmental Radioactivity (ALMERA) network, both in emergency situations and for routine environmental monitoring purposes.

The method was established after an extensive review of papers from the scientific literature, and was tested and validated in terms of repeatability and trueness (relative bias) in accordance with International Organization for Standardization guidelines. Reproducibility tests were performed at expert laboratories. The calculation of massic activities, uncertainty budget, decision threshold and detection limit are also described.

The IAEA wishes to thank all the participants for their valuable contributions to develop the procedure and review this publication, especially C.-K. Kim for the design of the method development activity, R. Kleinschmidt and M. Cook (Australia), who assisted in drafting the manuscript and provided validation data, and the four expert laboratories (National Food Chain Safety Office, Hungary; Japan Chemical Analysis Center, Japan; Korea Institute of Nuclear Safety, Republic of Korea; and Jožef Stefan Institute, Slovenia) that contributed to the reproducibility testing. The IAEA officer responsible for this publication was A. Pitois of the IAEA Environment Laboratories.

EDITORIAL NOTE

This publication has been prepared from the original material as submitted by the contributors and has not been edited by the editorial staff of the IAEA. The views expressed remain the responsibility of the contributors and do not necessarily reflect those of the IAEA or the governments of its Member States.

This publication has not been edited by the editorial staff of the IAEA. It does not address questions of responsibility, legal or otherwise, for acts or omissions on the part of any person.

The use of particular designations of countries or territories does not imply any judgement by the publisher, the IAEA, as to the legal status of such countries or territories, of their authorities and institutions or of the delimitation of their boundaries.

The mention of names of specific companies or products (whether or not indicated as registered) does not imply any intention to infringe proprietary rights, nor should it be construed as an endorsement or recommendation on the part of the IAEA.

The contributors are responsible for having obtained the necessary permission for the IAEA to reproduce, translate or use material from sources already protected by copyrights.

The IAEA has no responsibility for the persistence or accuracy of URLs for external or third party Internet web sites referred to in this publication and does not guarantee that any content on such web sites is, or will remain, accurate or appropriate.

CONTENTS

1.	INTRODUCTION						
2.	SCOPE						
3.	NORMATIVE REFERENCES						
4.	RESPONSIBILITIES						
5.	TERM	S AND DEFINITIONS, SYMBOLS AND ABBREVIATIONS	3				
6.	PRINC	TPLE	5				
7.	EQUIF	PMENT AND CHEMICALS	6				
	7.1. 7.2.	EQUIPMENT CHEMICALS	6				
8.	PROC	EDURE	7				
	8.1. 8.2. 8.3. 8.4.	SEPARATION OF RADIUM BY PRECIPITATION PURIFICATION OF RADIUM SOURCE PREPARATION MEASUREMENT CONDITIONS SETTING AND DETERMINATION OF COUNTING EFFICIENCY 8.4.1. Source preparation for alpha/beta discrimination setting 8.4.2. Alpha/beta discrimination setting 8.4.3. Region of interest (ROI) setting 8.4.4. Counting efficiencies of ²²⁶ Ra and ²²⁸ Ra and chemical recovery	7 7 8 8 8 9 9 9 10				
9.	EXPRI	ESSION OF RESULTS	12				
	 9.1. 9.2. 9.3. 9.4. 	CALCULATION OF MASSIC ACTIVITIES OF ²²⁶ RA AND ²²⁸ RA AT THE SAMPLING DATE CALCULATION OF INDIVIDUAL UNCERTAINTY AND COMBINED UNCERTAINTY DECISION THRESHOLD DETECTION LIMIT	12 13 13 13				
10.	VALII	DATION OF THE PROCEDURE	14				
	10.1. 10.2. 10.3. 10.4. 10.5. 10.6.	LINEARITY, RANGE OF MEASUREMENT TRUENESS OF THE METHOD PRECISION OF THE METHOD REPEATABILITY REPRODUCIBILITY ACCEPTANCE CRITERIA	14 14 15 15 15 15				
11.	TEST REPORT16						
12.	CONCLUSION						

APPENDIX I. U AND TH DECAY SERIES	17
APPENDIX II. FLOW CHART FOR CHEMICAL SEPARATION OF RADIUM	19
APPENDIX III. EXPERIMENTAL CONDITIONS AS A FUNCTION OF MATRIX CONTENT	21
APPENDIX IV. ALPHA/BETA DISCRIMINATION SETTING	23
APPENDIX V. TYPICAL SPECTRA OF ²²⁶ RA AND ²²⁸ RA	25
APPENDIX VI. COUNTING EFFICIENCIES AND OVERALL EFFICIENCIES	27
APPENDIX VII. CALCULATION OF INDIVIDUAL UNCERTAINTY AND COMBINED UNCERTAINTY OF ACTIVITY CONCENTRATION OF ²²⁶ RA AND ²²⁸ RA IN WATER SAMPLES	29
APPENDIX VIII. VALIDATION RESULTS OF THE PROCEDURE	35
REFERENCES	41
CONTRIBUTORS TO DRAFTING AND REVIEW	43

1. INTRODUCTION

Radium is a radioactive element for which no stable isotope is known. The mass numbers of the known radioisotopes range from 206 to 230. Among them, two radium radioisotopes, ²²⁶Ra and ²²⁸Ra, are very significant from a radiological protection point of view due to their relatively long half-lives, presence in nature as radionuclides from the uranium and thorium decay series, and high dose conversion factors. The uranium and thorium decay series are shown in Figs. 1–2 in Appendix I. Radium radioisotopes are important because they can be easily incorporated into bones due to having similar properties to other elements from Group II (i.e. calcium) and producing short-lived radionuclides of high massic activity. For the general public, the consumption of food and water containing radium may lead to the accumulation of these radionuclides in the body, contributing to radiological dose.

Radium radioisotopes are often present in industrial waste products, referred to as Naturally Occurring Radioactive Material (NORM). It encompasses substances with enhanced concentrations of natural radionuclides, which occur at production sites and in the environment as a result of industrial or mining activities. Therefore, national and international regulators have shown an increasing interest in regulating radium levels in most industries and consumer products for the protection of human health and environment.

For radiological dose assessment, geochronological studies and environmental monitoring purposes, radium radioisotopes have been measured by several techniques, including gamma-ray spectrometry, alpha spectrometry, mass spectrometry and liquid scintillation counting techniques.

The radium radioisotopes ²²³Ra, ²²⁴Ra and ²²⁶Ra are alpha-particle emitters and can be directly measured by alpha spectrometry. The radium isotope ²²⁸Ra, which is a beta-particle emitter, can indirectly be measured by alpha spectrometry via its progeny ²²⁸Th. However, because of the long half-life of ²²⁸Th ($T_{1/2} = 1.9$ y), the in-growth period ranges between 3 and 18 months (depending on the ²²⁸Ra activity and required sensitivity), and hence alpha spectrometry is not an adequate measurement method for routine analysis of ²²⁸Ra.

Gamma-ray spectrometry is a useful non-destructive method that permits the simultaneous determination of many radionuclides in a bulk sample [1-2], without the need for complex and time consuming radiochemical separations as undertaken for alpha spectrometry. However, it is limited by the weak emission probabilities of many potentially useful emission lines, the relatively poor efficiency of the High Purity Germanium (HPGe) detectors over a wide range of energy, the difficult task of precisely calibrating the efficiency of the detector, and the need to evaluate self-absorption effects [3].

Mass spectrometric methods have been used for the determination of ²²⁶Ra and ²²⁸Ra, although mass spectrometric determination of ²²⁸Ra has only occasionally been reported due to its short half-life (5.75 y). The main advantages of mass spectrometric methods are high sensitivity and short analysis times, although radiochemical pre-concentration and separation procedures are usually unavoidable. The major problems in mass spectrometric analysis are the formation of isobaric and polyatomic interferences, the formation of multiple charged ions, the abundance sensitivity (or tailing) of neighbouring atoms, matrix effects and mass fractionation during detection [4].

Liquid scintillation counting (LSC) is often used for the determination of radium radioisotopes in environmental samples, especially in liquids [5-6]. Methods using LSC can provide high sample throughput, and acceptable detection limits for several radium radioisotopes [7-10]. Quenching, poor energy resolution and relatively high background, especially in comparison with alpha spectrometry, are however limitations of this technique [11-12].

Gross alpha/beta activity analysis has been widely applied to drinking water [13]. This method is however non radionuclide-specific and faces many sources of interferences leading to a large spread of results (up to two orders of magnitude), as observed in interlaboratory comparisons [14–15].

Taking into account the availability of LSC to simultaneously measure ²²⁶Ra and ²²⁸Ra at the same time with acceptable detection limit, this publication introduces a procedure for rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water, to be used for screening of activity levels in emergency situations and for routine environmental monitoring purposes. The procedure was validated in accordance with ISO guidelines.

2. SCOPE

This publication describes a method for the rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water using liquid scintillation counting following the chemical separation of radium. The method can be used for the fast detection of contamination of drinking water by radium in emergency situations or for routine environmental monitoring purposes.

The method is suitable for the analysis of 0.5 L of drinking water containing less than 100 mg kg^{-1} barium.

3. NORMATIVE REFERENCES

The principles of the normative references, which have been applied to this document, can be found in the list below.

- ISO 11929:2010. Determination of the characteristic limits (decision threshold, detection limit and limits of the confidence interval) for measurements of ionizing radiation Fundamentals and application;
- ISO/TS 21748:2004. Guidance for the use of repeatability, reproducibility and trueness estimates in measurement uncertainty estimation;
- ISO/IEC Guide 98-3:2008. Uncertainty of measurement Part 3: Guide to the expression of uncertainty in measurement;
- ISO 5725-1:1994. Accuracy (trueness and precision) of measurement methods and results Part 1: General principles and definitions;

- ISO 5725-2:1994. Accuracy (trueness and precision) of measurement methods and results Part 2: Basic method for the determination of repeatability and reproducibility of a standard measurement method;
- ISO 5725-4:1994. Accuracy (trueness and precision) of measurement methods and results Part 4: Basic methods for the determination of the trueness of a standard measurement method;
- ISO 7870-1:2007. Control charts Part 1: General guidelines;
- QUAM:2012.P1 EURACHEM/CITAC Guide, 2012. Quantifying Uncertainty in Analytical Measurement.

4. **RESPONSIBILITIES**

It is the responsibility of the analyst to follow this procedure and to establish appropriate safety and health practices. Compliance with national and/or internal chemical and radiological safety rules is obligatory. Although the method has been carefully tested each laboratory should perform an individual validation procedure.

5. TERMS AND DEFINITIONS, SYMBOLS AND ABBREVIATIONS

In the following, u(x) indicates the uncertainty in quantity x:

 A_x : Certified massic activity of the analyte on reference date (Bq·kg⁻¹).

 A_x^t : Massic activity of the analyte in the quality control sample on reference date (Bq·kg⁻¹).

 a_x : Massic activity of the specific analyte on the sampling date (Bq·kg⁻¹).

 C_x : Target massic activity of the specific analyte in the quality control sample prepared for the validation of the procedure (Bq·kg⁻¹).

 m_{s-x} : Mass of the certified standard solution taken for the analysis of the specific analyte (kg).

 m_{t-x} : Mass of the quality control sample taken for the analysis of the specific analyte (kg).

 m_s : Mass of sample (kg).

 n_x^s : Net count rates of the specific analyte in the certified standard solution (counts \cdot s⁻¹).

 n_x^t : Net count rates of the specific analyte in the quality control sample (counts \cdot s⁻¹).

 n_x : Net count rates of the specific analyte in the sample (counts \cdot s⁻¹).

PI: Precision index (%).

 $R_{226Ra:}$ Chemical recovery of ²²⁶Ra.

- R_{228Ra} : Chemical recovery of ²²⁸Ra.
- R_L : Reproducibility limit (Bq·kg⁻¹).
- r_L : Repeatability limit (Bq·kg⁻¹).
- r_{g-x}^{s} : Gross count rates of the specific analyte in the standard solution (counts s⁻¹).
- r_{g-x}^{t} : Gross count rates of the specific analyte in the quality control sample (counts \cdot s⁻¹).
- r_{g-x} : Gross count rates of the specific analyte in the sample (counts s⁻¹).
- r_{0-x} : Gross count rates of the specific analyte in the blank sample (counts s⁻¹).
- S_r : Standard deviation of repeatability (Bq·kg⁻¹).
- S_R : Standard deviation of reproducibility (Bq·kg⁻¹).
- T_{s-x} : Counting time of the specific analyte in the sample (s).
- t_{0-x} : Counting time of the specific analyte in the blank (s).
- t_{s-x} : Time interval between the measurement date and the reference date of the specific analyte in the certified standard solution (s).
- t_{t-x} : Time interval between the measurement date and the reference date of the specific analyte in the quality control sample (s).
- t_x : Time interval between the measurement date and the sampling date (s).
- ε_x : Counting efficiency of the specific analyte.
- ε_x^c : Overall efficiency of the specific analyte in the quality control sample.
- λ_x : Decay constant of the specific analyte (s⁻¹).
- $\overline{X_x}$: Mean of all measured values of the specific analyte in the quality control sample for the validation of the procedure (Bq·kg⁻¹).
- x_x^* : Decision threshold of the specific analyte (Bq·kg⁻¹).
- ξ_x : Detection limit of the specific analyte (Bq·kg⁻¹).
- δ : Relative bias of the method (%).
- u_{mb} : Uncertainty of the analytical balance (kg).
- AAS: Atomic Absorption Spectroscopy.
- ALMERA: Analytical laboratories for the measurement of environmental radioactivity.
- EDTA: Ethylenediaminetetraacetic acid.

FWHM: Full width at half-maximum.

HDPE: High-density polyethylene.

HPGe: High purity germanium

IAEA: International atomic energy agency.

ICP-OES: Inductively coupled plasma optical emission spectroscopy.

ICP-MS: Inductively coupled plasma mass spectrometry.

ISO: International organization for standardization.

LSC: Liquid scintillation counter.

MCA: Multichannel analyser.

MDL: Minimum detectable limit.

NORM: Naturally occurring radioactive material.

PP: Polypropylene.

ROI: Region of interest.

rpm: Revolutions per minute.

6. **PRINCIPLE**

Barium co-precipitation is used as a method of separation for radium due to the very similar chemical properties of barium and radium. The exploitation of the ability of barium to react with an excess of sulphate ions to produce a precipitate allows the quantitative analysis of environmental activity concentrations of radium in water.

The inclusion of a lead carrier allows the removal of ²¹⁰Pb from solution, which increases the accuracy of ²²⁸Ra measurement, as ²¹⁰Pb can cause spectral interference. The removal of ²¹⁰Pb is achieved by lowering the pH of the solution to re-precipitate barium sulphate using acetic acid in which lead sulphate is soluble. This allows ²¹⁰Pb to remain in solution and therefore be removed.

The source preparation is achieved by suspending the barium sulphate precipitate in the EDTA solution. Barium sulphate is insoluble in water, alkalis and acids, but EDTA increases the solubility due to the complexation of barium and the speciation effect. The EDTA molecule inhibits barium sulphate nucleation and added heat forces complexation. This enables the use of a naphthalene-based scintillation cocktail to gain better spectral resolution than with the use of a gel-forming cocktail.

The flow chart of the procedure is given in Appendix II.

Massic activities of ²²⁶Ra and ²²⁸Ra in the sample are calculated from net count rates of the sample, sample amount and the overall efficiency that can be obtained from spiked sample with known activities of ²²⁶Ra and ²²⁸Ra, and that shows the ability of the method to extract radium (chemical recovery) as well as the ability (counting efficiency) of the instrument to detect it.

7. EQUIPMENT AND CHEMICALS

7.1. EQUIPMENT

- Standard laboratory equipment;
- Hot plate with magnetic stirrer;
- Centrifuge;
- Analytical balance with accuracy of 0.1 mg;
- Water bath with temperature controller;
- Vortex;
- Low background liquid scintillation counter (LSC);
- Elemental analysis technique for Ba and Ca determination;
- pH meter or pH papers;
- 600 mL glass beaker;
- 50 mL HDPE or PP centrifuge tubes;
- 5 mL and 10 mL Eppendorf pipettes;
- Micropipette for 50 µL volume;
- Transfer pipettes;
- 20 mL low background plastic liquid scintillation vials.

7.2. CHEMICALS

- Pb carrier: 2.397 g Pb(NO₃)₂, 0.5 mL 14.4M HNO₃ made up to 100 mL with de-ionized water;
- Ba carrier: 2.836 g BaCl₂.2H₂O, 0.5 mL 14.4M HNO₃ made up to 100 mL with de-ionized water;
- 0.25M EDTA;
- 10M HCl (32%);
- 14.4M HNO₃ (65%);
- 9M HSO₄;
- 25% NH₃ solution;
- 96 % glacial acetic acid;
- Analytical grade ammonium sulphate (solid);
- 10% ammonium sulphate solution;
- Scintillation cocktail suitable for alpha/beta separation (Optiphase Hisafe 3)
- ²²⁶Ra and ²²⁸Ra certified standard solutions;
- Alpha emitter standard solution (241 Am or 210 Po or 242 Pu);
- Beta emitter standard solution 90 Sr/ 90 Y.

8. **PROCEDURE**

Prior to the start of the analysis of 226 Ra and 228 Ra in the water sample, it is recommended to determine the Ca and Ba contents of the sample using an elemental analysis technique (for example, AAS, ICP-OES or ICP-MS). If the Ba concentration is higher than the limit (100 mg·kg⁻¹), it is recommended to reduce the volume of the sample to be analysed so that the total content of Ba in the sample does not exceed 50 mg.

8.1. SEPARATION OF RADIUM BY PRECIPITATION

- (1) Acidify 0.5 L of the water sample to approximately pH 2 using drops of 10M HCl.
- (2) Add 2 mL of Pb carrier and 2 mL of Ba carrier to the solution. In case that the total content of Ba in the sample solution is 50 mg or higher, it is not necessary to add Ba carrier to the sample.
- (3) Add 4 mL of 9M HSO₄ and 5 g of ammonium sulphate accurately weighed.
- (4) Stir the solution to ensure that all solids are dissolved, allow the precipitate to form and then let the precipitate settle.
- (5) Decant the supernatant without disturbing the precipitate, leaving less than 30 mL of liquid in the glassware.
- (6) Quantitatively transfer the precipitate and the limited amount of liquid to a 50 mL HDPE or PP centrifuge tube rinsing the beaker with de-ionized water to avoid loss of precipitate.
- (7) Centrifuge the solution for five minutes at 3500 rpm.
- (8) Decant the excess supernatant carefully without disturbing the precipitate.
- 8.2. PURIFICATION OF RADIUM
- (1) Dissolve the precipitate in 10 mL of hot 0.25M EDTA solution (the solution volume is changeable, depending on total contents of Ca and Ba in the sample, and is given in Table 1 in Appendix III) and 3 mL of 25% NH₃ solution. Carefully agitate the solution to dissolve the precipitate.
- (2) Add 5 mL of 10% ammonium sulphate and adjust the pH to 4.2-4.5 using 96% glacial acetic acid. As the pH is lowered, the precipitate should begin to re-form.
- (3) Boil the solution in the HDPE or PP centrifuge tube in a water bath at 80°C for two minutes, cool it with cold tap water, then centrifuge for five minutes at 3500 rpm. The resulting precipitate is purified barium (radium) sulphate precipitate. Discard the supernatant.
- (4) Dissolve the precipitate in 10 mL of hot 0.25M EDTA solution (the solution volume is changeable, depending on total content of Ba in the sample, and is given in Table 2 in Appendix III) and agitate the solution carefully to dissolve the precipitate into solution. This may require the use of a vortex.
- (5) Add 3 mL of 10% ammonium sulphate solution to the solution and adjust the pH to 4.5 using the 96% glacial acetic acid. The precipitate should begin to re-form. Separate the precipitate by centrifugation for five minutes at 3500 rpm. (Precaution: don't leave the precipitate for overnight and don't allow it to dry. Prepare the source as soon as the precipitate has been prepared).

8.3. SOURCE PREPARATION

- (1) Wash the precipitate twice with 20 mL of de-ionized water, shake well, centrifuge and discard the supernatant to ensure removal of any residual solvents.
- (2) Add 3 mL of hot 0.25M EDTA to suspend the precipitate. Make sure to break up the precipitate at this point either through vortex or suction through a transfer pipette.
- (3) Heat the suspension in a water bath at a moderate temperature (30-50°C) for a minimum of 30 minutes.
- (4) Vortex the sample to ensure the precipitate gets well suspended and no cluster of precipitate is visible in the suspension. Steps (3) and (4) may need to be repeated for obtaining a homogeneous suspension and the absence of cluster of precipitate.
- (5) Quantitatively transfer the solution including partially dissolved BaSO₄ precipitate to a 20 mL plastic liquid scintillation vial using a transfer pipette
- (6) Rinse the HDPE or PP centrifuge tube with another 1 mL of hot 0.25M EDTA solution to ensure that no analyte remains in the tube.
- (7) Add 14 mL of Optiphase Hisafe 3 liquid scintillation cocktail to the source solution and vortex or shake well until the solution appears homogenous. (Precaution: the addition of the cocktail should be done all at once or in large portions to avoid any reaction with the source solution that could cause a cloudy, inhomogenous mixture).
- (8) Seal and shake the LSC source until the suspension appears homogenous.
- (9) Clean the vial with an alcohol wipe to remove any static interference.

Precaution: Sample stability is determined by the suspension of the precipitate. If the precipitate is not well reacted with the EDTA solution, it will fall out more easily and quickly. It is suggested to break up the precipitate in the EDTA solution either through vortex, shaking or continued suction and expulsion through a transfer pipette for stubborn precipitates.

Heat is also a major factor in the reaction between the precipitate and the EDTA solution. If a good suspension is difficult to achieve, heating in a water bath for a further 30 minutes followed by vortex/shaking will help force the reaction. This should be repeated until the sample can be held still without any settling occurring within a few minutes.

The grade of the ammonium sulphate used also has a bearing on the quality of the precipitate. Reagent grade ammonium sulphate produces a precipitate that is slightly chunkier and stickier and therefore harder to suspend during source preparation. Analytical grade ammonium sulphate produces a precipitate that is much easier to work with.

8.4. MEASUREMENT CONDITIONS SETTING AND DETERMINATION OF COUNTING EFFICIENCY

8.4.1. Source preparation for alpha/beta discrimination setting

- (1) Add 2 mL of barium carrier to two 50 mL volume HDPE or PP centrifuge tubes.
- (2) Add 3 mL of 10% ammonium sulphate solution and 1 mL of 25% NH₃ solution to each solution to obtain the BaSO₄ precipitates. Separate the precipitates by centrifuging for five minutes at 3500 rpm.

- (3) Dissolve the precipitates in 4 mL of hot 0.25M EDTA solutions and agitate the solutions carefully to dissolve and suspend the precipitates into solution. This may require the use of a vortex.
- (4) Follow section 8.3 steps (5) to (7).
- (5) Add 10–100 Bq of alpha and beta emitters, respectively (e.g. ²⁴¹Am or ²¹⁰Po or ²⁴²Pu and 90 Sr/⁹⁰Y) in 50 µL volume solutions.
- (6) Follow section 8.3 steps (8) to (9).

8.4.2. Alpha/beta discrimination setting

The discrimination factor integrates the tail of the pulse and compares it with the total charge in the pulse. The higher the discrimination factor setting is, the more pulses will be directed into the short pulse category (beta) and vice versa. This allows the simultaneous determination of alphas and betas.

The determination of an optimum discrimination factor requires two standards, one pure alpha and one pure beta emitter, ²⁴¹Am or ²¹⁰Po or ²⁴²Pu and ⁹⁰Sr/⁹⁰Y respectively. These radionuclides are used rather than ²²⁶Ra and ²²⁸Ra, as the latter are accompanied by progeny in-growth, which creates uncertainty in the determination of a discrimination factor. To arrive at the optimum setting, each standard is counted individually at a range of discrimination factor settings and the percent misclassification of alphas into the beta channel and vice versa are plotted against discrimination factor on the same graph. The point of minimum cross-over is selected as the optimum discrimination factor value.

As can be seen in Table 3 in Appendix IV, for testing of the method, a range of discrimination factors from 80 to 150 was chosen. Two blank samples were made following the above procedure, one spiked with ²⁴¹Am and one with ⁹⁰Sr/⁹⁰Y after source preparation. The activity concentrations of these spikes do not need to be documented as it is the raw count spillover that is being determined and not a final activity concentration. These samples were counted for 60 minutes each at each discrimination factor within the selected range. The region of interest used to analyse these samples utilised the full range of the instrument from channel 0 to channel 1024. The alpha and beta counts of each sample are then tabulated and the number of betas in the alpha MCA and vice versa is calculated, as can be seen in Table 3 in Appendix IV. These values are then graphed to determine the optimum discrimination factor for this laboratory was calculated to be 130. The optimum discrimination factor has to be determined at each laboratory performing this procedure.

8.4.3. Region of interest (ROI) setting

The setting of multichannel analysers (MCA) and windows for the measurement of radium is performed using the standard Alpha/Beta configuration. This involves the first MCA analysing all beta events in the first half and all alpha events in the second half, as determined by the discrimination factor. The second MCA then analyses all guard anticoincidence events in the first half and all guard coincidence events in the second half. These guard events are not considered as true sample events. Any liquid scintillation analyser with the ability to discriminate between alpha and beta events will have a standard Alpha/Beta configuration for use with this method.

The lower and upper limits of the analysis windows are set using the known emission energies of ²²⁶Ra and ²²⁸Ra, approximately 4500–5000 keV and 0–50 keV, respectively. When using channels in the MCA, this converts to approximately 600–640 and 50–300, respectively. Each liquid scintillation counter should be calibrated for radium analysis, including setting of regions, prior to any analysis.

Regions of interest are set specifically to contain only the analyte of interest. The range of each region has been set following the width of the FWHM of each peak. This enables a greater accuracy when calculating specific activity. It also negates the use of spillover calculations to correct for interfering peaks in the ROI as the difference in the number of counts lost and gained in each region is negligible as shown in Fig. 4 to 7 in Appendix V.

8.4.4. Counting efficiencies of ²²⁶Ra and ²²⁸Ra and chemical recovery

8.4.4.1. Counting efficiencies of 226 Ra and 228 Ra determined by direct measurement of 226 Ra and 228 Ra standard sources

The counting efficiencies of the liquid scintillation counter for ²²⁶Ra and ²²⁸Ra were determined using two blank samples, resulting from the BaSO₄ precipitates in de-ionized water. The samples were spiked after preparation of BaSO₄ precipitates, one with ²²⁶Ra and one with ²²⁸Ra. Spiking of the samples after preparation eliminates the chemical recovery variable. For the determination, each spiked sample was counted seven times. This provided a reasonable data set for calculations. The counting efficiencies of ²²⁶Ra and ²²⁸Ra were calculated using equations (1) and (2):

$$\mathcal{E}_{226Ra} = \frac{n_{226Ra}^{s}}{A_{226Ra} \times m_{s-226Ra} \times e^{-\lambda_{226Ra} \times t_{s-226Ra}}}$$
(1)

$$\mathcal{E}_{228Ra} = \frac{n_{228Ra}^{s}}{A_{228Ra} \times m_{s-228Ra} \times e^{-\lambda_{228Ra} \times t_{s-228Ra}}}$$
(2)

The counting efficiencies were determined to be 1.03 and 0.646 with 6% and 10% of relative standard deviation for ²²⁶Ra and ²²⁸Ra respectively, as shown in Fig. 8 in Appendix VI.

8.4.4.2. Overall efficiencies of ²²⁶Ra and ²²⁸Ra determined by analysis of quality control samples with known activities of ²²⁶Ra and ²²⁸Ra

The overall efficiency describes the ability of the method to extract radium as well as the ability of the instrument to detect it. These parameters are known as chemical recovery and counting efficiency. This ensures that all factors of sample preparation have been taken into account. As part of a thorough validation process, these parameters have been considered separately.

A quality control sample (spiked sample with known activities of ²²⁶Ra and ²²⁸Ra) was used for the determination of the overall efficiency of ²²⁶Ra and ²²⁸Ra. Overall efficiencies for both ²²⁶Ra and ²²⁸Ra were calculated using equations (3) and (4).

$$\varepsilon_{226Ra}^{c} = \frac{n_{226Ra}^{t}}{A_{226Ra}^{t} \times m_{t-226Ra} \times e^{-\lambda_{226Ra} \times t_{t-226Ra}}}$$
(3)

$$\varepsilon_{228Ra}^{c} = \frac{n_{228Ra}^{t}}{A_{228Ra}^{t} \times m_{t-228Ra} \times e^{-\lambda_{228Ra} \times t_{t-228Ra}}}$$
(4)

The control chart for overall efficiencies of 226 Ra and 228 Ra is shown in Fig. 9 in Appendix VI. The overall efficiencies of 226 Ra and 228 Ra were determined to be 0.908 and 0.516 with 6.6 % and 18.9 % of relative standard deviation for 226 Ra and 228 Ra, respectively.

8.4.4.3. Chemical recovery

Using counting efficiency and overall efficiency data, the chemical recovery was determined using equations (5) and (6).

$$R_{226Ra} = \frac{\varepsilon_{226Ra}^c}{\varepsilon_{226Ra}}$$
(5)

$$R_{228Ra} = \frac{\varepsilon_{228Ra}^c}{\varepsilon_{228Ra}} \tag{6}$$

The chemical recovery was calculated indirectly using this method instead of the conventional gravimetric or tracer techniques due to the chemistry of the procedure and the technique of source preparation used. The addition of ¹³³Ba as a tracer has previously shown inconsistent recovery, possibly due to the chemistry used during the procedure, and it is therefore not suitable. Gravimetric determination requires drying of the precipitate, which causes inconsistent source preparation and is therefore not suitable. In this procedure, the chemical recovery was not used for the calculation of specific activities of ²²⁶Ra and ²²⁸Ra in the sample, but used for quality control purpose only. The chemical recoveries of ²²⁶Ra and ²²⁸Ra ranged from 78–90%.

8.4.4.4. Measurement

Samples are counted using a liquid scintillation counter with pulse shape analysis for 60 minutes each including quality control and blank samples. The spectra should be inspected for reasonable alpha peak resolution and any visible quenching (movement of spectra outside regions of interest). Those are shown in Appendix V. For quality assurance purpose, a quality control sample and a blank were analysed for every set of samples. Samples should be counted immediately after source preparation to avoid extensive in-growth of progeny and the degradation of alpha spectrum due to the formation of flakes of precipitation.

As illustrated in Figs. 10–11 in Appendix VI, samples have a period of approximately 4 hours of stability within 2% of the variation of the overall efficiency of 226 Ra in a shaken source and 12% in an unshaken source due to in-growth of progeny and sinking of the BaSO₄ precipitate.

High selectivity of the method is provided by chemical separation and the ability of pulse shape analysis afforded by liquid scintillation counting. However, possible interferences cannot be excluded. If the solid precipitate is not completely suspended in EDTA, alpha peak broadening can occur due to the inability of alpha particles contained within unsuspended precipitate to fluoresce efficiently.

9. EXPRESSION OF RESULTS

9.1. CALCULATION OF MASSIC ACTIVITIES OF ²²⁶RA AND ²²⁸RA AT THE SAMPLING DATE

Activities of ²²⁶Ra and ²²⁸Ra in water sample were expressed in the procedure as units of massic activity ($Bq\cdot kg^{-1}$). In case the activities of ²²⁶Ra and ²²⁸Ra are expressed as units of activity concentration ($Bq\cdot L^{-1}$), the density of the sample should be measured to convert the unit of massic activity to activity concentration. All nuclear data used for the calculation of massic activities of ²²⁶Ra and ²²⁸Ra were taken from the recommended data of the DDEP database from the Laboratoire National Henri Becquerel [16] and Nudat 2.5 [17].

Massic activities of both ²²⁶Ra and ²²⁸Ra are calculated using the counts per minute as given by the spectral analysis software, separate efficiency factors as determined from the validation process and the background count rate as specified by the blank sample data.

The peak resolution using this method is good and therefore the use of spillover calculations is not necessary. This is due to the fact that the ROI is very specific and the difference between counts lost and gained in each region is negligible. Any background counts produced due to the environment or chemiluminescence is also corrected for, since blank/background samples are always run simultaneously.

In this procedure, massic activities of 226 Ra and 228 Ra in water sample are calculated based on the overall efficiencies of 226 Ra and 228 Ra obtained from the analysis of the quality control sample with known activities of 226 Ra and 228 Ra as shown in equations (7) and (8).

$$a_{226Ra} = \frac{n_{226Ra}}{\varepsilon_{226Ra}^{c} \times m_{S}} \times e^{\lambda_{226Ra} \times t_{226Ra}}$$
(7)

$$a_{228Ra} = \frac{n_{228Ra}}{\varepsilon_{228Ra}^c \times m_S} \times e^{\lambda_{228Ra} \times t_{228Ra}}$$
(8)

It is therefore recommended that real samples should be analysed with minimum one quality control sample and one blank sample at the same time. The analyst should pay attention to the quality control of overall efficiencies of ²²⁶Ra and ²²⁸Ra, because they are changeable; testing showed changes within a relative standard deviation of 7% for ²²⁶Ra and 19% for ²²⁸Ra, as shown in Fig. 9 in Appendix VI. It is advised that the analyst prepares control charts for counting efficiency and overall efficiency of ²²⁶Ra and ²²⁸Ra, similar to the ones shown in Figs. 8–9 in Appendix VI.

9.2. CALCULATION OF INDIVIDUAL UNCERTAINTY AND COMBINED UNCERTAINTY

The combined uncertainty is calculated according to the guide to the expression of uncertainty in measurement [18], and is described in Appendix VII.

An example of spreadsheets for the measurement of ²²⁶Ra and ²²⁸Ra is given in Table 4 in Appendix VII. It can be noticed that the overall efficiencies of ²²⁶Ra and ²²⁸Ra are a dominant contributor (90%) to the combined uncertainty of the measurement result.

9.3. DECISION THRESHOLD

In accordance with ISO 11929:2010, the decision threshold, x_x^* , is obtained from the following equation (9):

$$x_x^* = k_{1-\alpha} \times u(0_x) \tag{9}$$

 $u(0_x)$ is the uncertainty of the measurement when the activity of the analyte is equal to zero, and $k_{1-\alpha} = 1.65$ is chosen referring to the probability that error of the first kind (wrongly rejecting the hypothesis that the analyte is present) is not greater than 5%. $u(0_x)$ can be quantified as:

$$u(0_x) = \frac{1}{\varepsilon_x \times m_s} \times \sqrt{\frac{r_{0-x}}{T_{s-x}} + \frac{r_{0-x}}{t_{0-x}}}$$
(10)

9.4. DETECTION LIMIT

In accordance with ISO 11929:2010, the detection limit, ξ_x^* , is calculated by the following equations (11) and (12), assuming $\alpha = \beta$ then $k_{I-\alpha} = k_{I-\beta} = k$,

$$\xi_x = 2 \times q_x \tag{11}$$

$$q_{x} = k_{1-\alpha} \times u(0_{x}) + \frac{1}{2} \times \frac{k_{1-\alpha}^{2}}{a_{x}} \times [u^{2}(a_{x}) - u^{2}(0_{x})]$$
(12)

 a_x is the massic activity of the analyte and $u(0_x)$ is the uncertainty of the measurement when the activity of the analyte is equal to zero and $k_{I-\alpha} = 1.65$ is chosen.

The detection limit of ²²⁶Ra and ²²⁸Ra in the procedure are 0.01 Bq·kg⁻¹ and 0.06 Bq·kg⁻¹, respectively for a 1 hour counting time in a low background liquid scintillation spectrometer with 0.5 kg of sample.

10. VALIDATION OF THE PROCEDURE

For the validation of the procedure, four types of spiked water samples with different activity ratios of ²²⁸Ra/²²⁶Ra, used as quality control samples, and a blank were prepared. One of the quality control samples was used for testing the reproducibility and the other three quality control samples were used for repeatability and trueness testing. The repeatability, reproducibility and trueness of the procedure were tested according to ISO-5725 and ISO/TS-21748.

The following parameters were established via the validation process:

- Linear response range;
- Trueness;
- Precision;
- Repeatability;
- Reproducibility;
- Detection limit.

Quality control samples are run to ensure the accuracy of both the method and the instrumentation. Blank samples are run to monitor reagent radiological purity and ensure that there is no cross-contamination during the procedure. This also includes the monitoring of the background to ensure that there is no excess background activity being counted by the instrumentation that could affect the results.

10.1. LINEARITY, RANGE OF MEASUREMENT

The linear range of measurement of the method was tested by measuring the ²²⁶Ra and ²²⁸Ra activities at various activity concentrations in de-ionized water. Results of measurements as a function of reference values have been tabulated and graphed in Fig. 12 in Appendix VIII. Linearity of the method was proven in the range 0–100 Bq·kg⁻¹ for both ²²⁶Ra and ²²⁸Ra.

It is very likely that the method can be used in a wider range, from detection limit to approximately 1000 Bq·kg⁻¹. This makes this method suitable for ²²⁶Ra and ²²⁸Ra massic activities determination ranging from detection limit to higher than environmental levels.

10.2. TRUENESS OF THE METHOD

The trueness of the procedure was tested using twelve replicates of three types of quality control samples and calculated using equation (13).

$$\delta(\%) = \frac{C_x - \overline{X_x}}{C_x} \times 100 \tag{13}$$

Trueness of the method for the three types of quality control samples is given in Tables 5–7 in Appendix VIII.

10.3. PRECISION OF THE METHOD

The precision index was determined using twelve replicates of three types of quality control samples and calculated using the equation (14):

$$PI\% = 100 \times \sqrt{\frac{u(C_x)^2}{C_x^2} + \frac{u(a_x)^2}{a_x^2}}$$
(14)

Precision of the method for the three types of quality control samples is given in Tables 5 to 7 in Appendix VIII.

10.4. REPEATABILITY

The repeatability limit was estimated by analysing twelve replicates of three types of quality control samples with known activity concentrations of ²²⁶Ra and ²²⁸Ra in repeatability conditions. The repeatability limit is calculated using equation (15):

$$r_L = S_r \times 2.8 \tag{15}$$

The results of repeatability testing in one laboratory are shown in Tables 5 to 7 in Appendix VIII. The procedure provides 7% of relative repeatability standard deviation (S_r) for 1–2 Bq·kg⁻¹ of ²²⁶Ra and ²²⁸Ra activity level, as shown in Table 7 in Appendix VIII.

10.5. REPRODUCIBILITY

The reproducibility limit was estimated by analysing ten replicates of a quality control sample with known activity concentrations of ²²⁶Ra and ²²⁸Ra in five laboratories. The reproducibility limit is calculated using the equation (16):

$$R_L = S_R \times 2.8 \tag{16}$$

The results of reproducibility testing obtained from an inter-laboratory exercise with five laboratories are shown in Table 8 in Appendix VIII. The relative reproducibility standard deviation is determined to be 17% for ²²⁶Ra and 15% for ²²⁸Ra at 0.5 Bq·kg⁻¹ activity level of both isotopes in inter-laboratory comparison exercise.

10.6. ACCEPTANCE CRITERIA

According to the criteria of ISO/TS 21748:2004, the test results meet the requirements if:

- the absolute mean bias $|\delta|$ in the repeatability test is smaller than the repeatability limit (r_L) : $|\delta| < r_{L_1}$
- the absolute mean bias $|\delta|$ in the reproducibility test is smaller than the reproducibility limit (R_L): $|\delta| < R_{L_i}$
- PI < 25%.

According to these criteria, the procedure for the rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water samples performs well.

11. TEST REPORT

The test report shall conform to ISO 17025 requirements and shall contain the following information:

- Sample code;
- Sample description;
- Sampling date (if known);
- Analysis request code;
- Responsible person;
- Reference to the measurement and evaluation procedure used;
- Test result;
- Units in which the test result is expressed;
- Date to which the test result has been corrected.

Unless otherwise required, the test result should be stated together with the expanded uncertainty [19] calculated using a coverage factor k = 2.

The following form is recommended: $a_x \pm u(a_x)$ with the associated k value.

12. CONCLUSION

A procedure was developed and tested for the rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water using liquid scintillation counting. Radium is radiochemically separated from interfering elements using co-precipitation steps. The analysis time for chemical separation and source preparation of radium for a set of samples (3 samples, one blank and one quality control sample) is about 8 hours.

According to the repeatability and reproducibility tests, the measurement process is considered to perform adequately.

The method can be used for the rapid determination of ²²⁶Ra and ²²⁸Ra in drinking water, both for screening of activity levels in emergency situations and for routine environmental monitoring purposes.

APPENDIX I. U AND TH DECAY SERIES



FIG. 1. ²²⁶Ra is a member of the ²³⁸U decay series [17].



FIG. 2. ²²⁴Ra and ²²⁸Ra are members of the ²³²Th decay series [17].

APPENDIX II. FLOW CHART FOR CHEMICAL SEPARATION OF RADIUM



APPENDIX III. EXPERIMENTAL CONDITIONS AS A FUNCTION OF MATRIX CONTENT

TABLE 1. VOLUME OF 0.25M EDTA SOLUTION REQUIRED DEPENDING ON CONTENTS OF CALCIUM AND BARIUM IN THE SAMPLE AT THE FIRST DISSOLUTION STEP OF BaSO₄ PRECIPITATE (8.2. STEP 1)

Contents of Ca and E	Volume of 0.25 M EDTA	
Са	Ba	solution (mL)
0	50	15
0	75	26
0	100	35
300	30	52
400	30	65
500	30	75

TABLE 2. VOLUME OF 0.25M EDTA SOLUTION REQUIRED DEPENDING ON THE CONTENT OF BARIUM IN THE SAMPLE SOLUTION AT THE SECOND DISSOLUTION STEP OF BaSO₄ PRECIPITATE (8.2. STEP 4)

Content of Ba in the sample	Volume of 0.25M EDTA
(mg)	solution (mL)
50	10
75	15
100	25

APPENDIX IV. ALPHA/BETA DISCRIMINATION SETTING

Discrimination	Alpha	counts	Beta c	counts				
factor					Total alpha	Total beta	Alpha	Beta
	Alpha MCA	Beta MCA	Alpha MCA	Beta MCA	counts	counts	spinover	spinover
80	4596	7.49	10193	26500	4603	36693	0.277783	0.001628
90	4627	6.43	7823	30505	4633	38328	0.204105	0.001388
100	4622	6.94	5032	34832	4628	39864	0.126225	0.001499
110	4608	8.86	2514	39945	4617	42459	0.059216	0.001919
120	4591	10.1	909	45474	4601	46384	0.019608	0.002201
130	4609	11.8	231	49241	4621	49472	0.004677	0.002553
140	4606	25.8	45.6	50674	4632	50720	0.000899	0.005564
150	4546	69.6	11.3	50783	4616	50795	0.000223	0.015072

TABLE 3. RESULTS OF OPTIMUM ALPHA/BETA DISCRIMINATION SETTING



FIG. 3. Correlation between spillover and discrimination factor.

APPENDIX V. TYPICAL SPECTRA OF ²²⁶RA AND ²²⁸RA



FIG. 4. Negligible spectral overlap of (a) ²²⁶Ra spectra and (b) ²²⁸Ra spectra.



FIG. 5. ²²⁶Ra spectra and counting window as shown in spectral analysis software.



Channels FIG. 6. ²²⁸*Ra spectra and counting window as shown in spectral analysis software, alpha peaks are* ²²⁸*Ra decay progeny.*



FIG. 7. Combined spectra of ²²⁶Ra and ²²⁸Ra as shown in spectral analysis software.



FIG. 8. Control chart for counting efficiencies of ²²⁶Ra and ²²⁸Ra.



FIG. 9. Control chart for overall efficiencies of ²²⁶Ra and ²²⁸Ra.



FIG. 10. Variation of overall efficiency for ²²⁶Ra with elapsed time after source preparation.



FIG. 11. Variation of overall efficiency for ²²⁸Ra with elapsed time after source preparation.

APPENDIX VII. CALCULATION OF INDIVIDUAL UNCERTAINTY AND COMBINED UNCERTAINTY OF ACTIVITY CONCENTRATION OF ²²⁶RA AND ²²⁸RA IN WATER SAMPLES

VII.1. Uncertainty in weighing the sample mass

The uncertainty associated with the mass of sample is estimated using the data from the calibration certificate of the analytical balance and the manufacturer's recommendations on the uncertainty estimation of the balance. If it is assumed that a balance with u_{mb} mass uncertainty is used for weighing the sample and the linearity distribution is assumed to show a triangular distribution, the standard uncertainty in the weighing of the sample mass can be calculated as:

$$u_{\sin gle_weighing} = \frac{u_{mb}}{\sqrt{3}}$$
(17)

The contribution for the linearity has to be accounted for twice, once for the tare and once for the gross mass (tare + sample), leading to an uncertainty $u(m_s)$ according to equation (18):

$$u(m_s) = \sqrt{2 \times \left(\frac{u_{mb}}{\sqrt{3}}\right)^2}$$
(18)

VII.2. Uncertainty of net counts of ²²⁶Ra and ²²⁸Ra in standard solution, spiked sample and real sample

The uncertainty of net count rates at ROIs of ²²⁶Ra and ²²⁸Ra in LSC spectra can be calculated as follows:

$$n_{226Ra}^{s} = r_{g-226Ra}^{s} - r_{0-226Ra} \quad \text{and} \quad u(n_{226Ra}^{s}) = \sqrt{\left(\frac{r_{g-226Ra}^{s}}{T_{s-226Ra}}\right) + \left(\frac{r_{0-226Ra}}{t_{0-226Ra}}\right)}$$
(19)

$$n_{228Ra}^{s} = r_{g-228Ra}^{s} - r_{0-228Ra} \quad \text{and} \quad u(n_{228Ra}^{s}) = \sqrt{\left(\frac{r_{g-228Ra}^{s}}{T_{s-228Ra}}\right) + \left(\frac{r_{0-228Ra}}{t_{0-228Ra}}\right)}$$
(20)

 n_{226Ra}^{t} , $u(n_{226Ra}^{t})$, n_{228Ra}^{t} and $u(n_{228Ra}^{t})$ for spiked sample, and n_{226Ra} , $u(n_{226Ra})$, n_{228Ra} and $u(n_{228Ra})$ for real sample can be calculated in the same way with n_{226Ra}^{s} and $u(n_{226Ra}^{s})$ for standard solution.

VII.3. Uncertainties of counting efficiency and overall efficiency of ²²⁶Ra and ²²⁸Ra

The uncertainty of the counting efficiency $u(\varepsilon_{226Ra})$ can be calculated based on equation (1). Differentiating equation (1), each partial derivative can be calculated as follows:

$$\frac{\partial \varepsilon_{226Ra}}{\partial n_{226Ra}^s} = \frac{\varepsilon_{226Ra}}{n_{226Ra}^s}$$
(21)

$$\frac{\partial \varepsilon_{226Ra}}{\partial A_{226Ra}} = \frac{\varepsilon_{226Ra}}{A_{226Ra}}$$
(22)

$$\frac{\partial \varepsilon_{226Ra}}{\partial m_{s-226Ra}} = \frac{\varepsilon_{226Ra}}{m_{s-226Ra}}$$
(23)

$$\frac{\partial \varepsilon_{226Ra}}{\partial \lambda_{226Ra}} = -\varepsilon_{226Ra} \times t_{s-226Ra}$$
(24)

$$\frac{\partial \varepsilon_{226Ra}}{\partial t_{s-226Ra}} = \varepsilon_{226Ra} \times \lambda_{226Ra}$$
(25)

The uncertainty of counting efficiency $(u(\varepsilon_{226Ra}))$ can be calculated from the law of error propagation as follows:

$$u(\varepsilon_{226Ra})^{2} = \left(\frac{\varepsilon_{226Ra}}{n_{226Ra}^{s}}\right)^{2} \times u(n_{226Ra}^{s})^{2} + \left(\frac{\varepsilon_{226Ra}}{A_{226Ra}}\right)^{2} \times u(A_{226Ra})^{2} + \left(\frac{\varepsilon_{226Ra}}{m_{s-226Ra}}\right)^{2} \times u(m_{s-226Ra})^{2} + (\varepsilon_{226Ra} \times t_{s-226Ra})^{2} \times u(\lambda_{226Ra})^{2} + (\varepsilon_{226Ra} \times \lambda_{226Ra})^{2} \times u(t_{s-226Ra})^{2}$$
(26)

The uncertainty of the counting efficiency $u(\varepsilon_{228Ra})$ can be calculated based on equation (2). Differentiating equation (2), each partial derivative can be calculated as follows:

$$\frac{\partial \varepsilon_{228\,Ra}}{\partial n_{228\,Ra}^s} = \frac{\varepsilon_{228\,Ra}}{n_{228\,Ra}^s}$$
(27)

$$\frac{\partial \varepsilon_{228Ra}}{\partial A_{228Ra}} = -\frac{\varepsilon_{228Ra}}{A_{228Ra}}$$
(28)

$$\frac{\partial \varepsilon_{228Ra}}{\partial m_{s-228Ra}} = -\frac{\varepsilon_{228Ra}}{m_{s-228Ra}}$$
(29)

$$\frac{\partial \varepsilon_{228Ra}}{\partial \lambda_{228Ra}} = \varepsilon_{228Ra} \times t_{s-228Ra}$$
(30)

$$\frac{\partial \varepsilon_{228Ra}}{\partial t_{s-228Ra}} = \varepsilon_{228Ra} \times \lambda_{228Ra}$$
(31)

The uncertainty of counting efficiency $(u(\varepsilon_{228Ra}))$ can be calculated from the law of error propagation as follows:

$$u(\varepsilon_{228Ra})^{2} = \left(\frac{\varepsilon_{228Ra}}{n_{228Ra}^{s}}\right)^{2} \times u(n_{228Ra}^{s})^{2} + \left(-\frac{\varepsilon_{228Ra}}{A_{228Ra}}\right)^{2} \times u(A_{228Ra})^{2} + \left(-\frac{\varepsilon_{228Ra}}{m_{s-228Ra}}\right)^{2} \times u(m_{s-228Ra})^{2} + (\varepsilon_{228Ra} \times t_{s-228Ra})^{2} \times u(\lambda_{228Ra})^{2} + (\varepsilon_{228Ra} \times \lambda_{228Ra})^{2} \times u(t_{s-228Ra})^{2}$$
(32)

Overall efficiency uncertainties of ²²⁶Ra and ²²⁸Ra, $u(\varepsilon_{226Ra}^{c})$ and $u(\varepsilon_{228Ra}^{c})$ can be calculated in the same way with calculation of counting efficiency uncertainties, $u(\varepsilon_{226Ra})$ and $u(\varepsilon_{228Ra})$. The uncertainties of $t_{s-226Ra}$ and $t_{s-228Ra}$ are regarded negligible.

VII.4. Uncertainty of chemical recovery

Uncertainties of chemical recovery of ²²⁶Ra and ²²⁸Ra are calculated as follows:

$$u(R_{226Ra}) = \frac{\varepsilon_{226Ra}^{c}}{\varepsilon_{226Ra}} \times \sqrt{\left(\frac{u(\varepsilon_{226Ra}^{c})}{\varepsilon_{226Ra}^{c}}\right)^{2} + \left(\frac{u(\varepsilon_{226Ra})}{\varepsilon_{226Ra}}\right)^{2}}$$
(33)

$$u(R_{228Ra}) = \frac{\varepsilon_{228Ra}^{c}}{\varepsilon_{228Ra}} \times \sqrt{\left(\frac{u(\varepsilon_{228Ra}^{c})}{\varepsilon_{228Ra}^{c}}\right)^{2} + \left(\frac{u(\varepsilon_{228Ra})}{\varepsilon_{228Ra}}\right)^{2}}$$
(34)

VII.5. Combined uncertainty in the determination of the analyte at the sampling date

The combined uncertainty of the specific activity of 226 Ra in the sample is calculated based on equation (7) using the general formula given elsewhere [18, 19]. Differentiating equation (7), each partial derivative can be calculated as follows:

$$\frac{\partial a_{226Ra}}{\partial n_{226Ra}} = \frac{a_{226Ra}}{n_{226Ra}}$$
(35)

$$\frac{\partial a_{226Ra}}{\partial \varepsilon_{226Ra}} = -\frac{a_{226Ra}}{\varepsilon_{226Ra}}$$
(36)

$$\frac{\partial a_{226Ra}}{\partial m_s} = -\frac{a_{226Ra}}{m_s}$$
(37)

$$\frac{\partial a_{226Ra}}{\partial \lambda_{226Ra}} = a_{226Ra} \times t_{226Ra}$$
(38)

$$\frac{\partial a_{226Ra}}{\partial t_{226Ra}} = a_{226Ra} \times \lambda_{226Ra}$$
(39)

The uncertainty of counting efficiency $(u(\varepsilon_{226Ra}))$ can be calculated from the law of error propagation as follows:

$$u(a_{226Ra})^{2} = \left(\frac{a_{226Ra}}{n_{226Ra}}\right)^{2} \times u(n_{226Ra})^{2} + \left(-\frac{a_{226Ra}}{\varepsilon_{226Ra}}\right)^{2} \times u(\varepsilon_{226Ra})^{2} + \left(-\frac{a_{226Ra}}{m_{s}}\right)^{2} \times u(m_{s})^{2} + \left(a_{226Ra} \times t_{226Ra}\right)^{2} \times u(\lambda_{226Ra})^{2} + \left(a_{226Ra} \times \lambda_{226Ra}\right)^{2} \times u(t_{226Ra})^{2}$$
(40)

The uncertainty of $t_{226Ra} \times u(t_{226Ra})$ is often regarded negligible.

The combined uncertainty of the specific activity of 228 Ra on the sampling date can be calculated by the same way with that of 226 Ra.

Components	Value of variable	Standard uncertainty	Percent contribution to combined standard uncertainty
⁻²²⁶ Ra			
r _{g-226Ra} (cps)	0.651	0.014	10
r _{0-226Ra} (cps)	0.00113	0.00057	
E _{226Ra}	0.908	0.060	90
Half-life of ²²⁶ Ra (y)	1600	7	negligible
m _s (kg)	0.508	0.00008	negligible
$a_{226Ra}, (Bq \cdot kg^{-1})$	1.41	0.09	
²²⁸ Ra			
r _{g-228Ra} (cps)	0.1986	0.0075	8
$r_{0-228Ra}$ (cps)	0.040	0.0034	2
E _{228Ra}	0.516	0.097	90
Half-life of ²²⁸ Ra (y)	5.75	0.04	negligible
m _s (kg)	0.508	0.00008	negligible
a _{228Ra} (Bq·kg ⁻¹)	0.599	0.100	

TABLE4.EXAMPLEOFSPREADSHEETSFORTHEMEASUREMENTUNCERTAINTIES OF226 Ra AND228 Ra IN SPIKED WATER SAMPLE





	²²⁶ Ra			²²⁸ Ra			
Number of run	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index	
1	1.39±0.09	-9.7	6.5	$0.69{\pm}0.07$	-5.5	10.5	
2	1.52 ± 0.09	-1.3	6.0	$0.79{\pm}0.07$	8.2	9.3	
3	1.52 ± 0.09	-1.3	6.0	0.81 ± 0.07	11.0	9.1	
4	1.54 ± 0.09	0.0	5.9	0.81 ± 0.07	11.0	2.7	
5	1.52 ± 0.09	-1.3	6.0	$0.80{\pm}0.07$	9.6	9.2	
6	1.47 ± 0.09	-4.5	6.2	$0.84{\pm}0.07$	15.1	8.8	
7	1.47 ± 0.09	-4.5	6.2	0.81 ± 0.07	11.0	9.1	
8	1.35 ± 0.08	-12.3	6.0	0.87 ± 0.07	19.2	8.5	
9	1.38 ± 0.08	-10.4	5.8	0.87 ± 0.07	19.2	8.5	
10	1.30 ± 0.08	-15.6	6.2	$0.89{\pm}0.07$	21.9	8.3	
11	1.33±0.08	-13.6	6.0	$0.84{\pm}0.07$	15.1	9.9	
12	1.40±0.10	-9.1	7.2	0.80 ± 0.07	9.6	12.8	
Mean	1.43			0.82			
Absolute mean bias	0.11			0.09			
Repeatability standard deviation (S _r)	0.08			0.05			
Relative repeatability standard deviation Sr (%)	5.59			6.10			
Repeatability limit (r _L)	0.22			0.20			
Assigned value at reference date	1.54±0.03*			0.73±0.02*			
*: coverage factor (k =1)							

TABLE 5. REPEATABILITY TESTING RESULTS FOR $^{226}\mathrm{Ra}$ and $^{228}\mathrm{Ra}$ in quality control sample 1

		²²⁶ Ra		²²⁸ Ra		
Number of run	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index
1	1.01 ± 0.07	0.0	7.0	1.56±0.09	6.1	6.7
2	$1.00{\pm}0.07$	-1.0	7.1	1.62 ± 0.09	10.2	6.5
3	1.01 ± 0.07	0.0	7.0	1.60 ± 0.10	8.8	7.1
4	1.03 ± 0.07	2.0	6.9	1.50 ± 0.09	2.0	3.4
5	0.98 ± 0.07	-3.0	7.2	1.54 ± 0.09	4.8	6.8
6	1.01 ± 0.07	0.0	7.0	1.56 ± 0.09	6.1	6.7
7	0.99 ± 0.07	-2.0	7.1	1.54 ± 0.09	4.8	6.8
8	0.97 ± 0.07	-4.0	7.3	1.60 ± 0.10	8.8	7.1
9	0.95 ± 0.07	-5.9	7.4	1.60±0.10	8.8	7.1
10	0.98 ± 0.07	-3.0	7.2	1.50±0.09	2.0	6.9
11	1.02 ± 0.07	1.0	6.9	1.53±0.09	4.1	6.8
12	1.00±0.1	-1.0	10.0	1.60±0.20	8.8	13.0
Mean	1.00			1.56		
Absolute mean bias	0.01			0.09		
Repeatability standard deviation (S _r)	0.02			0.04		
Relative repeatability standard deviation S _r (%)	2.00			2.56		
Repeatability limit (r _L)	0.06			0.11		
Assigned value at reference date	1.01±0.01*			1.47±0.05*		
*: coverage factor (k=1)					

TABLE 6. REPEATABILITY TESTING RESULTS FOR $^{\rm 226}\rm{Ra}$ and $^{\rm 228}\rm{Ra}$ in quality control sample 2

		²²⁶ Ra		²²⁸ Ra		
Number of run	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index	massic activity (Bq·kg ⁻¹)	relative bias (%)	precision index
1	2.0±0.1	2.0	5.0	2.5±0.1	8.2	5.0
2	2.1 ± 0.1	2.9	4.8	2.4±0.1	3.9	5.2
3	2.0 ± 0.1	2.0	5.0	2.3±0.1	0.4	5.3
4	2.0 ± 0.1	2.0	5.0	2.5±0.1	8.2	3.0
5	2.0 ± 0.1	2.0	5.0	2.5±0.1	8.2	5.0
6	1.7 ± 0.1	16.7	5.9	2.3±0.1	0.4	5.3
7	1.8 ± 0.1	11.8	5.6	2.5±0.1	8.2	5.0
8	1.8 ± 0.1	11.8	5.6	2.8±0.1	21.2	4.7
9	1.8 ± 0.1	11.8	5.6	2.7±0.1	16.9	4.8
10	1.7 ± 0.1	16.7	5.9	2.7±0.1	16.9	4.8
Mean	1.89			2.52		
Absolute mean bias	0.15			0.21		
Repeatability standard deviation (S _r)	0.14			0.17		
Relative repeatability standard deviation Sr (%)	7.41			6.75		
Repeatability limit (r _L)	0.39			0.48		
Assigned value at reference date	2.04±0.01*			2.31±0.07*		

TABLE 7. REPEATABILITY TESTING RESULTS FOR 226 Ra AND 228 Ra IN QUALITYCONTROL SAMPLE 3

*: coverage factor (k =1)

TABLE 8. REPRODUCIBILITY TEST RESULTS FOR THE DETERMINATION OF ^{226}Ra and ^{228}Ra in a quality control sample by inter-laboratory comparison exercise with 5 laboratories (unit: Bq·kg⁻¹)

Factors	²²⁶ Ra	²²⁸ Ra
Mean of measurement results at target date	0.55	0.54
Assigned value at target date	$0.40{\pm}0.02^{*}$	$0.52{\pm}0.02^{*}$
Absolute mean bias	0.15	0.02
Repeatability variance (S_r^2)	0.002	0.001
Between laboratory variance (S_L^2)	0.002	0.001
Reproducibility variance (S_R^2)	0.004	0.003
Repeatability limit (r _L)	0.13	0.11
Reproducibility limit (R _L)	0.17	0.15

*: coverage factor (k =1)

REFERENCES

- [1] RIHS, R., CONDOMINES, M., An improved method for Ra isotope (²²⁶Ra, ²²⁸Ra, ²²⁴Ra) measurements by gamma spectrometry in natural waters: Application to CO₂-rich thermal waters from the French Massif Central, Chem. Geol. **182** (2002) 409–421.
- [2] SEMKOW, T.M., et al., Low-background gamma spectrometry for environmental radioactivity, Appl. Radiat. Isot. **57** (2002) 213–223.
- [3] JODLOWSKI, P., Self-absorption correction in gamma-ray spectrometry of environmental samples An overview of methods and correction values obtained for the selected geometries, Nukleonika **51** (2006) 21–25.
- [4] ZORIY, M.V., VARGA, Z., PICKHARDT, C., OSTAPCZUK, P., HILLE, R., HALICZ, L., SEGAL, I., BECKER, J.S., Determination of ²²⁶Ra at ultratrace level in mineral water samples by sector field inductively coupled plasma mass spectrometry, J. Environ. Monit. 7 (2005) 514–518.
- [5] VRSKOVA, M., et.al., Comparison of determining of Ra-228 via LSC and gamma spectrometry in mineral waters, in Proceedings of Advances in Liquid Scintillation Spectrometry, 2005, Tucson, Arizona, USA.
- [6] WALLNER, G., et.al., Determination of natural radionuclides in drinking water from Waldviertel, Austria, in Proceedings of Advances in Liquid Scintillation Spectrometry, 2008, Tucson, Arizona, USA.
- [7] CHAŁUPNIK, S., LEBECKA J., Determination of ²²⁶Ra, ²²⁸Ra and ²²⁴Ra in water and aqueous solutions by liquid scintillation counting, in Proceedings of Advances in Liquid Scintillation Spectrometry, 1993, Tucson, Arizona, USA.
- [8] REPINC, U., BENEDIK, L., Natural radioactivity in drinking and mineral water, in Proceedings of Advances in Liquid Scintillation Spectrometry, 2001, Tucson, Arizona, USA.
- [9] SCHOENHOFER, F. et.al., Natural radionuclides in mineral waters sold in Upper Silesia, Poland: their measurement, doses and compliance with regulations, in Proceedings of Advances in Liquid Scintillation Spectrometry, 2008, Tucson, Arizona, USA.
- [10] KIM, Y.-J., et al., Simultaneous determination of Ra-226 and Pb-210 in ground water and soil samples using the liquid scintillation counter suspension gel method. Appl. Radiat. Isot. **54** (2001) 275–281.
- [11] KÖHLER, M., et al., Comparison of methods for the analysis of ²²⁶Ra in water samples, Appl. Radiat. Isot. **56** (1–2) (2002) 387–392.
- [12] BURNETT, W.C., TAI W.C., Determination of radium in natural waters by liquid scintillation, Anal. Chem. **64** (1992) 1691–1697.
- [13] JOBBAGY, V., WAETJEN, U., MERESOVA, J., Current status of gross alpha/beta activity analysis in water samples: A short overview of methods, J. Radioanal. Nucl. Chem. **286** (2010) 393–399.
- [14] JOBBAGY, V., MERESOVA, J., WAETJEN, U., Critical remarks on gross alpha/beta activity analysis in drinking waters: Conclusions from a European interlaboratory comparison, Appl. Radiat. Isot. **87** (2014) 429–434.
- [15] JOBBAGY, V., DIRICAN, A., WAETJEN, U., Radiochemical characterization of mineral waters for a European interlaboratory comparison, Microchem. J. 110 (2013) 675–680.

- [16] Monographie BIPM-5 Vol. 1 Monographie BIPM-5 Vol. 5, Decay Data Evaluation Project (DDEP) accessed on 20 March 2009. <u>http://www.nucleide.org/DDEP_WG/DDEPdata.htm</u>
- [17] NuDat 2.5, NuDat database version 2.5 accessed on 20 March 2009. http://www.nndc.bnl.gov/nudat2/
- [18] INTERNATIONAL ORGANIZATION FOR STANDARDIZATION, ISO/IEC Guide 98:2008, Guide to the expression of uncertainty in measurement, International Organization for Standardization, Geneva (2008).
- [19] ELLISON, S.L.R., ROSSLEIN, M., WILLIAMS, A. EURACHEM/CITAC Guide: Quantifying Uncertainty in Analytical Measurement, J. Anal. Chem, **58** (2002) 191–198.

CONTRIBUTORS TO DRAFTING AND REVIEW

Chmielewska, I.	Central Mining Institute, Katowice, Poland						
Cook, M.	Queensland Health, Brisbane, Australia						
Gaigl, F.A.	International Atomic Energy Agency						
Kleinschmidt, R.	Queensland Health, Brisbane, Australia						
Kim, C.K.	International Atomic Energy Agency						
Kim, Y.J.	Korea Institute of Nuclear Safety, Republic of Korea						
Kis-Benedek, G.	International Atomic Energy Agency						
Kozar Logar, J.	Jožef Stefan Institute, Ljubljana, Slovenia						
Martin, P.	International Atomic Energy Agency						
Maurizio, F.	Agenzia Regionale per la Protezione dell'Ambiente della Lombardia, Italy						
Pitois, A.	International Atomic Energy Agency						
Rusconi, R.	Lombardia, Agenzia Regionale per la Protezione dell'Ambiente della Lombardia, Italy						
Sansone, U.	International Atomic Energy Agency						
Sanada, T.	Japan Chemical Analysis Center, Japan						
Tarjan, S.	International Atomic Energy Agency						
Varga, B.	National Food Chain Safety Office, Food and Feed Safety Directorate, Hungary						



ORDERING LOCALLY

In the following countries, IAEA priced publications may be purchased from the sources listed below or from major local booksellers.

Orders for unpriced publications should be made directly to the IAEA. The contact details are given at the end of this list.

AUSTRALIA

DA Information Services

648 Whitehorse Road, Mitcham, VIC 3132, AUSTRALIA Telephone: +61 3 9210 7777 • Fax: +61 3 9210 7788 Email: books@dadirect.com.au • Web site: http://www.dadirect.com.au

BELGIUM

Jean de Lannoy Avenue du Roi 202, 1190 Brussels, BELGIUM Telephone: +32 2 5384 308 • Fax: +32 2 5380 841 Email: jean.de.lannoy@euronet.be • Web site: http://www.jean-de-lannoy.be

CANADA

Renouf Publishing Co. Ltd.

5369 Canotek Road, Ottawa, ON K1J 9J3, CANADA Telephone: +1 613 745 2665 • Fax: +1 643 745 7660 Email: order@renoufbooks.com • Web site: http://www.renoufbooks.com

Bernan Associates

4501 Forbes Blvd., Suite 200, Lanham, MD 20706-4391, USA Telephone: +1 800 865 3457 • Fax: +1 800 865 3450 Email: orders@bernan.com • Web site: http://www.bernan.com

CZECH REPUBLIC

Suweco CZ, spol. S.r.o. Klecakova 347, 180 21 Prague 9, CZECH REPUBLIC Telephone: +420 242 459 202 • Fax: +420 242 459 203 Email: nakup@suweco.cz • Web site: http://www.suweco.cz

FINLAND

Akateeminen Kirjakauppa PO Box 128 (Keskuskatu 1), 00101 Helsinki, FINLAND Telephone: +358 9 121 41 • Fax: +358 9 121 4450 Email: akatilaus@akateeminen.com • Web site: http://www.akateeminen.com

FRANCE

Form-Edit

5 rue Janssen, PO Box 25, 75921 Paris CEDEX, FRANCE Telephone: +33 1 42 01 49 49 • Fax: +33 1 42 01 90 90 Email: fabien.boucard@formedit.fr • Web site: http://www.formedit.fr

Lavoisier SAS

14 rue de Provigny, 94236 Cachan CEDEX, FRANCE Telephone: +33 1 47 40 67 00 • Fax: +33 1 47 40 67 02 Email: livres@lavoisier.fr • Web site: http://www.lavoisier.fr

L'Appel du livre

99 rue de Charonne, 75011 Paris, FRANCE Telephone: +33 1 43 07 50 80 • Fax: +33 1 43 07 50 80 Email: livres@appeldulivre.fr • Web site: http://www.appeldulivre.fr

GERMANY

Goethe Buchhandlung Teubig GmbH

Schweitzer Fachinformationen Willstätterstrasse 15, 40549 Düsseldorf, GERMANY Telephone: +49 (0) 211 49 8740 • Fax: +49 (0) 211 49 87428 Email: s.dehaan@schweitzer-online.de • Web site: http://www.goethebuch.de

HUNGARY

Librotade Ltd., Book Import PF 126, 1656 Budapest, HUNGARY Telephone: +36 1 257 7777 • Fax: +36 1 257 7472 Email: books@librotade.hu • Web site: http://www.librotade.hu

INDIA

Allied Publishers

1st Floor, Dubash House, 15, J.N. Heredi Marg, Ballard Estate, Mumbai 400001, INDIA Telephone: +91 22 2261 7926/27 • Fax: +91 22 2261 7928 Email: alliedpl@vsnl.com • Web site: http://www.alliedpublishers.com

Bookwell

3/79 Nirankari, Delhi 110009, INDIA Telephone: +91 11 2760 1283/4536 Email: bkwell@nde.vsnl.net.in • Web site: http://www.bookwellindia.com

ITALY

Libreria Scientifica "AEIOU"

Via Vincenzo Maria Coronelli 6, 20146 Milan, ITALY Telephone: +39 02 48 95 45 52 • Fax: +39 02 48 95 45 48 Email: info@libreriaaeiou.eu • Web site: http://www.libreriaaeiou.eu

JAPAN

Maruzen Co., Ltd.

1-9-18 Kaigan, Minato-ku, Tokyo 105-0022, JAPAN Telephone: +81 3 6367 6047 • Fax: +81 3 6367 6160 Email: journal@maruzen.co.jp • Web site: http://maruzen.co.jp

NETHERLANDS

Martinus Nijhoff International Koraalrood 50, Postbus 1853, 2700 CZ Zoetermeer, NETHERLANDS Telephone: +31 793 684 400 • Fax: +31 793 615 698 Email: info@nijhoff.nl • Web site: http://www.nijhoff.nl

Swets Information Services Ltd.

PO Box 26, 2300 AA Leiden Dellaertweg 9b, 2316 WZ Leiden, NETHERLANDS Telephone: +31 88 4679 387 • Fax: +31 88 4679 388 Email: tbeysens@nl.swets.com • Web site: http://www.swets.com

SLOVENIA

Cankarjeva Zalozba dd Kopitarjeva 2, 1515 Ljubljana, SLOVENIA Telephone: +386 1 432 31 44 • Fax: +386 1 230 14 35 Email: import.books@cankarjeva-z.si • Web site: http://www.mladinska.com/cankarjeva_zalozba

SPAIN

Diaz de Santos, S.A. Librerias Bookshop • Departamento de pedidos Calle Albasanz 2, esquina Hermanos Garcia Noblejas 21, 28037 Madrid, SPAIN Telephone: +34 917 43 48 90 • Fax: +34 917 43 4023 Email: compras@diazdesantos.es • Web site: http://www.diazdesantos.es

UNITED KINGDOM

The Stationery Office Ltd. (TSO) PO Box 29, Norwich, Norfolk, NR3 1PD, UNITED KINGDOM Telephone: +44 870 600 5552 Email (orders): books.orders@tso.co.uk • (enquiries): book.enquiries@tso.co.uk • Web site: http://www.tso.co.uk

UNITED STATES OF AMERICA

Bernan Associates 4501 Forbes Blvd., Suite 200, Lanham, MD 20706-4391, USA Telephone: +1 800 865 3457 • Fax: +1 800 865 3450 Email: orders@bernan.com • Web site: http://www.bernan.com

Renouf Publishing Co. Ltd.

812 Proctor Avenue, Ogdensburg, NY 13669, USA Telephone: +1 888 551 7470 • Fax: +1 888 551 7471 Email: orders@renoufbooks.com • Web site: http://www.renoufbooks.com

United Nations

300 East 42nd Street, IN-919J, New York, NY 1001, USA Telephone: +1 212 963 8302 • Fax: 1 212 963 3489 Email: publications@un.org • Web site: http://www.unp.un.org

Orders for both priced and unpriced publications may be addressed directly to:

IAEA Publishing Section, Marketing and Sales Unit, International Atomic Energy Agency Vienna International Centre, PO Box 100, 1400 Vienna, Austria Telephone: +43 1 2600 22529 or 22488 • Fax: +43 1 2600 29302 Email: sales.publications@iaea.org • Web site: http://www.iaea.org/books

14-46971