



# Minimizing matrix effects in multi-residue pesticide analysis using a new sensitivity-enhanced triple quadrupole instrument in combination with extensive sample dilution

László Tölgyesi  
LC/MS scientist  
EFSC Demo Group

*International Symposium on Food Safety and Quality:  
Applications of Nuclear and Related Techniques*

# Pesticides screening – regulation and residue situation

- **EU regulation**

- Maximum residue levels regulated in commission regulation (EC) 396/2005
  - appendix II – IV regulates >170 000 maximum residue limits (matrix-pesticide combinations)
- European and national monitoring programs are in place also regulated by (EC) 396/2005
- 121 notifications related to pesticide residues in tea have been entered in the Rapid Alert System for Food and Feed (RASFF)



# Why do you want a 6495?

Utmost sensitivity is required for ...

- the analysis of a large number of relevant compounds with a single injection
- dilution of samples for reduced matrix effects, better accuracy and improved robustness
- improved LOQs for pesticides in complex matrices
- highest precision

## Experimental setup

- Small validation with tomato, orange, black tea
- QuEChERS extraction
- Dilution of black tea samples with acetonitrile prior injection
- Quantitation with solvent calibration – one calibration for all matrices?




# UHPLC/MS/MS Setup

## 1290 Infinity UHPLC


- column: Zorbax RRHD Eclipse Plus C18
- A: H<sub>2</sub>O / 0.1% HCOOH / 5mM NH<sub>4</sub>COOH
- B: MeOH / 0.1% HCOOH / 5mM NH<sub>4</sub>COOH

## 6495 LC/MS system


- single run, fast polarity switching, dynamic MRM



# Proven iFunnel Technology




## Agilent Jet Stream




- Thermal confinement of ESI plume
- Efficient desolvation
- Creates an ion rich zone

## Hexabore Capillary



- Six capillary bores
- Samples x6 times more ion rich gas from the source

## Dual Ion Funnel




- Removes the gas but captures the ions
- Removes neutral noise






# 6495 QQQ Technologies Continued Development



## Proven iFunnel Technology

- Agilent Jet Stream
- Hexabore Capillary
- Dual Ion Funnel
- Increased ion generation
- Enhanced ion sampling



1

- New Enhanced Q1 Ion Optics
- Improved ion transmission

2

- New Tapered Hexapole Collision Cell
- Effective ion collection and transmission

3

- New Detector with High Energy Conversion Dynode (up to 20 kV)
- Improved NEG ion detection with low noise



Agilent Technologies


IAEA 2014 @ Vienna

11/10/2014

6

# New Enhanced Q1 Ion Optics


*Improved transmission of ions and system robustness*



- New optimized MS 1 prefilter geometry for improved precursor **ion transmission**
  - Improved peak area response and peak area %RSD, more **sensitive** and **precise**
- New optical lens elements for reduced contamination and easier autotune
  - More **reliable** and **robust** performance

# New Curved and Tapered Hexapole Collision Cell

*Effective collection and transmission of ions*




- Curved and Tapered Hexapole Assembly for efficient collection and transmission of product ions
- Consistent collision cell pressure for **higher quality MS/MS spectra**
- Designed for consistent collision energies across all QQQ platforms




# New Detector with 20 kV High Energy Conversion Dynode

*More efficient detection of ions with low noise characteristics*




- Improved **ion detection** efficiency with High Energy Dynode (HED) voltage up to 20 kV
  - Improved peak area response and peak area %RSD in positive & negative ion mode
  - Improved **sensitivity** and **precision** for a wide m/z range of product ions
- Low noise level at 20 kV
  - Improved signal to noise

# 6495 QQQ Technologies Continued Development



## Proven iFunnel Technology

- Agilent Jet Stream
- Hexabore Capillary
- Dual Ion Funnel
- Increased ion generation
- Enhanced ion sampling



1

- New Enhanced Q1 Ion Optics
- Improved ion transmission

2

- New Tapered Hexapole Collision Cell
- Effective ion collection and transmission

3

- New Detector with High Energy Conversion Dynode
- Improved NEG ion detection with low noise

## Quantitative Applications

- Enhanced peak area response
- Improved peak area %RSD
- More sensitive and precise
- Lower Limits of Detection (IDL) and Quantitation (LLOQ)
- More reliable and robust



Agilent Technologies

IAEA 2014 @ Vienna

11/10/2014


10

# Bentazon: Superb Precision & Sub-fg IDL in Black Tea

Amount measured	Replicates	%RSD	t (99%)	IDL
1 pg/mL, 2.0 fg (LLOQ)	n = 7 injections	11.4	3.143	0.36 pg/mL, 0.72 fg


$$\text{IDL} = t \times (\% \text{RSD} / 100) \times \text{Amount} = 3.143 \times (11.4 / 100) \times 2.0 \text{ fg} = 0.72 \text{ fg}$$

*Injection volume = 2 µL*



- Excellent peak area **precision** (%RSD) are observed at the lowest level in black tea matrix

# Signal Improvements for High Relevance Pesticides




- Improved **sensitivity** (peak area gain of 3x) are observed on the new **6495** vs. the **6490**

# Improved Precision & Lower LLOQs in Black Tea




## Comparison of Area %RSD for 50 Pesticides at MRL



## # of pesticides detected at MRL

6495




- 67% of pesticides were easily detected with a %RSD < 20 with 1:100 dilution
- This means LLOQs < 20 pg/mL (ppt)
- %RSD < 20 meets SANCO guidelines

- Improved **precision** (%RSD) are observed, particularly at the lowest levels (LLOQs)
- The enhanced peak area response and improved precision (%RSD) means pesticides achieve **lower LLOQs** ( $\leq$  MRL level) using the new 6495.




# Dilution allows More Efficient Ionization of Pesticides

## Alanycarb, 10 µg/kg in Black Tea



## Oxamyl, 10 µg/kg in Black Tea





# Recovery (%) in Black Tea with Different Dilutions


Pesticides	No dilution (n = 5)	Dilution 1:5 (n = 5)	Dilution 1:10 (n = 5)	Dilution 1:20 (n = 5)	Dilution 1:50 (n = 5)	Dilution 1:100 (n = 5)
Acetamiprid	29.4 ± 0.8	57.3 ± 1.4	67.5 ± 3.7	<u>79.9 ± 2.9</u>	<u>91.8 ± 5.2</u>	<u>109.5 ± 3.4</u>
Alanycarb	10.4 ± 1.3	<u>73.9 ± 2.2</u>	<u>81.5 ± 14.3</u>	<u>85.7 ± 11.1</u>	<u>87.6 ± 4.7</u>	<u>121.7 ± 10.8</u>
Aldicarb	36.9 ± 1.0	69.9 ± 1.4	<u>78.0 ± 3.5</u>	<u>91.0 ± 4.2</u>	<u>95.2 ± 8.8</u>	<u>104.9 ± 14.1</u>
Carbaryl	56.9 ± 1.8	<u>80.1 ± 3.8</u>	<u>80.8 ± 4.1</u>	<u>96.1 ± 7.2</u>	<u>102.6 ± 6.6</u>	<u>116.4 ± 9.6</u>
Dimethoate	33.9 ± 1.7	68.6 ± 2.4	<u>84.1 ± 5.4</u>	<u>89.0 ± 7.9</u>	<u>88.2 ± 8.8</u>	<u>84.7 ± 7.5</u>
Diuron	<u>79.7 ± 4.0</u>	<u>90.4 ± 7.0</u>	<u>91.7 ± 4.9</u>	<u>94.9 ± 7.2</u>	<u>89.2 ± 7.3</u>	<u>100.9 ± 13.5</u>
Flufenoxuron	<u>95.4 ± 1.1</u>	<u>88.8 ± 1.6</u>	<u>89.4 ± 3.8</u>	<u>93.3 ± 5.8</u>	<u>100.0 ± 6.1</u>	<u>119.2 ± 13.9</u>
Monocrotophos	4.6 ± 0.3	13.9 ± 0.3	21.8 ± 0.8	33.8 ± 1.1	58.5 ± 2.0	<u>95.1 ± 5.7</u>
Oxamyl	20.8 ± 0.7	52.6 ± 1.9	65.0 ± 2.0	<u>79.7 ± 3.0</u>	<u>91.2 ± 4.6</u>	<u>110.6 ± 5.2</u>
Thiamethoxam	40.0 ± 1.4	45.9 ± 0.9	46.6 ± 3.8	52.2 ± 1.7	<u>70.9 ± 2.9</u>	<u>97.3 ± 2.0</u>

- Cells shaded in green shows %**recovery** of 70 – 120, which are in compliance with SANCO requirements.
- Pesticides show full recovery and no signal suppression with 1:100 dilution
- Neat solvent calibration curve can be used for pesticides quantitation with 1:100 dilution

# Summary

- The method benefits from the highly sensitive Agilent 6495 Triple Quadrupole LC/MS System.
- It takes full advantage of the low delay volumes of the Agilent 1290 Infinity LC System.
- Enhanced sensitivity allowed for more flexibility in the degree of sample dilution.
- Extensive dilution of sample extracts was applied to minimize matrix effects
- With any dilution, a lower matrix amount is introduced into the LC/MS system

# Thank you for your attention!



Agilent Technologies

IAEA 2014 @ Vienna

11/10/2014

20