

IAEA HUMAN HEALTH SERIES

No. 2

Quality Assurance Programme for Screen Film Mammography



IAEA

International Atomic Energy Agency

IAEA HUMAN HEALTH SERIES PUBLICATIONS

The mandate of the IAEA human health programme originates from Article II of its Statute, which states that the “Agency shall seek to accelerate and enlarge the contribution of atomic energy to peace, health and prosperity throughout the world”. The main objective of the human health programme is to enhance the capabilities of IAEA Member States in addressing issues related to the prevention, diagnosis and treatment of health problems through the development and application of nuclear techniques, within a framework of quality assurance.

Publications in the IAEA Human Health Series provide information in the areas of: radiation medicine, including diagnostic radiology, diagnostic and therapeutic nuclear medicine, and radiation therapy; dosimetry and medical radiation physics; and stable isotope techniques and other nuclear applications in nutrition. The publications have a broad readership and are aimed at medical practitioners, researchers and other professionals. International experts assist the IAEA Secretariat in drafting and reviewing these publications. Some of the publications in this series may also be endorsed or co-sponsored by international organizations and professional societies active in the relevant fields.

There are two categories of publications in this series:

IAEA HUMAN HEALTH SERIES

Publications in this category present analyses or provide information of an advisory nature, for example guidelines, codes and standards of practice, and quality assurance manuals. Monographs and high level educational material, such as graduate texts, are also published in this series.

IAEA HUMAN HEALTH REPORTS

Human Health Reports complement information published in the IAEA Human Health Series in areas of radiation medicine, dosimetry and medical radiation physics, and nutrition. These publications include reports of technical meetings, the results of IAEA coordinated research projects, interim reports on IAEA projects, and educational material compiled for IAEA training courses dealing with human health related subjects. In some cases, these reports may provide supporting material relating to publications issued in the IAEA Human Health Series.

All of these publications can be downloaded cost free from the IAEA web site:

<http://www.iaea.org/Publications/index.html>

Further information is available from:

Sales and Distribution Unit
International Atomic Energy Agency
Vienna International Centre
PO Box 100
1400 Vienna, Austria

Readers are invited to provide their impressions on these publications. Information may be provided via the IAEA web site, by mail at the address given above, or by email to:

Official.Mail@iaea.org.

QUALITY ASSURANCE PROGRAMME
FOR SCREEN FILM MAMMOGRAPHY

The following States are Members of the International Atomic Energy Agency:

AFGHANISTAN	GHANA	NIGERIA
ALBANIA	GREECE	NORWAY
ALGERIA	GUATEMALA	OMAN
ANGOLA	HAITI	PAKISTAN
ARGENTINA	HOLY SEE	PALAU
ARMENIA	HONDURAS	PANAMA
AUSTRALIA	HUNGARY	PARAGUAY
AUSTRIA	ICELAND	PERU
AZERBAIJAN	INDIA	PHILIPPINES
BAHRAIN	INDONESIA	POLAND
BANGLADESH	IRAN, ISLAMIC REPUBLIC OF	PORTUGAL
BELARUS	IRAQ	QATAR
BELGIUM	IRELAND	REPUBLIC OF MOLDOVA
BELIZE	ISRAEL	ROMANIA
BENIN	ITALY	RUSSIAN FEDERATION
BOLIVIA	JAMAICA	SAUDI ARABIA
BOSNIA AND HERZEGOVINA	JAPAN	SENEGAL
BOTSWANA	JORDAN	SERBIA
BRAZIL	KAZAKHSTAN	SEYCHELLES
BULGARIA	KENYA	SIERRA LEONE
BURKINA FASO	KOREA, REPUBLIC OF	SINGAPORE
BURUNDI	KUWAIT	SLOVAKIA
CAMEROON	KYRGYZSTAN	SLOVENIA
CANADA	LATVIA	SOUTH AFRICA
CENTRAL AFRICAN REPUBLIC	LEBANON	SPAIN
CHAD	LESOTHO	SRI LANKA
CHILE	LIBERIA	SUDAN
CHINA	LIBYAN ARAB JAMAHIRIYA	SWEDEN
COLOMBIA	LIECHTENSTEIN	SWITZERLAND
CONGO	LITHUANIA	SYRIAN ARAB REPUBLIC
COSTA RICA	LUXEMBOURG	TAJIKISTAN
CÔTE D'IVOIRE	MADAGASCAR	THAILAND
CROATIA	MALAWI	THE FORMER YUGOSLAV REPUBLIC OF MACEDONIA
CUBA	MALAYSIA	TUNISIA
CYPRUS	MALI	TURKEY
CZECH REPUBLIC	MALTA	UGANDA
DEMOCRATIC REPUBLIC OF THE CONGO	MARSHALL ISLANDS	UKRAINE
DENMARK	MAURITANIA	UNITED ARAB EMIRATES
DOMINICAN REPUBLIC	MAURITIUS	UNITED KINGDOM OF GREAT BRITAIN AND NORTHERN IRELAND
ECUADOR	MEXICO	UNITED REPUBLIC OF TANZANIA
EGYPT	MONACO	UNITED STATES OF AMERICA
EL SALVADOR	MONGOLIA	URUGUAY
ERITREA	MONTENEGRO	UZBEKISTAN
ESTONIA	MOROCCO	VENEZUELA
ETHIOPIA	MOZAMBIQUE	VIETNAM
FINLAND	MYANMAR	YEMEN
FRANCE	NAMIBIA	ZAMBIA
GABON	NEPAL	ZIMBABWE
GEORGIA	NETHERLANDS	
GERMANY	NEW ZEALAND	
	NICARAGUA	
	NIGER	

The Agency's Statute was approved on 23 October 1956 by the Conference on the Statute of the IAEA held at United Nations Headquarters, New York; it entered into force on 29 July 1957. The Headquarters of the Agency are situated in Vienna. Its principal objective is "to accelerate and enlarge the contribution of atomic energy to peace, health and prosperity throughout the world".

IAEA HUMAN HEALTH SERIES No. 2

QUALITY ASSURANCE PROGRAMME FOR SCREEN FILM MAMMOGRAPHY

INTERNATIONAL ATOMIC ENERGY AGENCY
VIENNA, 2009

COPYRIGHT NOTICE

All IAEA scientific and technical publications are protected by the terms of the Universal Copyright Convention as adopted in 1952 (Berne) and as revised in 1972 (Paris). The copyright has since been extended by the World Intellectual Property Organization (Geneva) to include electronic and virtual intellectual property. Permission to use whole or parts of texts contained in IAEA publications in printed or electronic form must be obtained and is usually subject to royalty agreements. Proposals for non-commercial reproductions and translations are welcomed and considered on a case-by-case basis. Enquiries should be addressed to the IAEA Publishing Section at:

Sales and Promotion, Publishing Section
International Atomic Energy Agency
Vienna International Centre
PO Box 100
1400 Vienna, Austria
fax: +43 1 2600 29302
tel.: +43 1 2600 22417
email: sales.publications@iaea.org
<http://www.iaea.org/books>

© IAEA, 2009

Printed by the IAEA in Austria
October 2009
STI/PUB/1381

IAEA Library Cataloguing in Publication Data

Quality assurance programme for screen film mammography. — Vienna :
International Atomic Energy Agency, 2009.
p. ; 29 cm. — (IAEA human health series, ISSN 2075-3772 ; no. 2)
STI/PUB/1381
ISBN 978-92-0-101609-6
Includes bibliographical references.

1. Breast — Radiography — Quality control. 2. Breast — Cancer —
Diagnosis — Quality control. I. International Atomic Energy Agency.
II. Series.

IAEAL

09-00592

FOREWORD

The application of radiation in human health, for both the diagnosis and treatment of disease, is an important component of the work of the IAEA. In the area of diagnostic radiology, this work is focused on quality assurance methods to both the promotion of the effective use of radiation for diagnostic outcome, through achieving and maintaining appropriate image quality, and also on dose determination to allow the monitoring and reduction of dose to the patient.

In response to heightened awareness of the importance of patient dose contributed by radiology procedures, the IAEA published *Dosimetry in Diagnostic Radiology: An International Code of Practice* (Technical Reports Series No. 457) in 2007, to form a basis for patient dose determination for the Member States. Further to this, it is recognized that for complex diagnostic procedures, such as mammography, a detailed guidance document is required to give the professionals in the clinical centre the knowledge necessary to assess the patient dose, as well as to ensure that the procedure gives the maximal patient benefit possible. It is well documented that without the implementation of a quality culture and a systematic quality assurance programme with appropriate education, the detection of breast cancer cannot be made at an early enough stage to allow effective curative treatment to be undertaken.

Currently there are a number of established quality assurance protocols in mammography from national and regional institutions, however, many of these protocols are distinctive and so a harmonized approach is required. This will allow the Member States to facilitate quality assurance in mammography in a standardized way which will also facilitate the introduction of national quality assurance programmes that are needed to underpin effective population screening programmes for breast cancer.

Development of a quality assurance document for screen film mammography was started in 2005 with the appointment of a drafting committee of international experts. The current publication is endorsed by the European Federation of Organisations for Medical Physics and the Asia-Oceania Federation of Organizations of Medical Physics.

The IAEA acknowledges the special contribution of the drafting committee chaired by M. Yaffe (Canada), with M. Chevalier (Spain), J.C. Heggie (Australia), P. Mora (Costa Rica) and K. Young (United Kingdom), and the American College of Radiology (ACR), which gave permission for the use of material from its Quality Assurance Manual. PAHO involvement is also acknowledged. The IAEA officers responsible for this publication were I.D. McLean (Division of Human Health), F. Pernička (Division of Human Health) and P. Ortiz López (Division of Radiation, Transport and Waste Safety).

EDITORIAL NOTE

Although great care has been taken to maintain the accuracy of information contained in this publication, neither the IAEA nor its Member States assume any responsibility for consequences which may arise from its use.

The use of particular designations of countries or territories does not imply any judgement by the publisher, the IAEA, as to the legal status of such countries or territories, of their authorities and institutions or of the delimitation of their boundaries.

The mention of names of specific companies or products (whether or not indicated as registered) does not imply any intention to infringe proprietary rights, nor should it be construed as an endorsement or recommendation on the part of the IAEA.

CONTENTS

1.	INTRODUCTION	1
1.1.	Why high quality in mammography is necessary	1
1.2.	Purpose	1
1.3.	Philosophy	2
2.	ELEMENTS OF HIGH QUALITY MAMMOGRAPHY	3
2.1.	Personnel	3
2.2.	Equipment	3
2.2.1.	Mammography unit	3
2.2.2.	Image receptor system	3
2.2.3.	Processing	4
2.2.4.	Viewing conditions	4
2.2.5.	Quality assurance	4
2.2.6.	Regular maintenance	4
3.	BASIC PRINCIPLES OF QUALITY ASSURANCE IN MAMMOGRAPHY	5
3.1.	Quality assurance activities	5
3.2.	Roles and responsibilities	7
3.2.1.	The owner	7
3.2.2.	Lead interpreting physician (radiologist)	7
3.2.3.	Radiographer (technologist)	7
3.2.4.	Medical physicist	8
4.	CLINICAL IMAGE QUALITY	9
4.1.	Patient positioning and compression	9
4.1.1.	Introduction	9
4.1.2.	Labelling of mammograms	11
4.1.3.	Breast compression	16
4.1.4.	Patient positioning	18
4.1.5.	Additional views	35
4.1.6.	Special circumstances	52
4.1.7.	Post-mastectomy imaging	59
4.2.	Clinical image evaluation	59
4.2.1.	Introduction	59
4.2.2.	Mediolateral oblique view positioning	60
4.2.3.	Craniocaudal view positioning	66
4.2.4.	Compression	70
4.2.5.	Exposure	72
4.2.6.	Contrast	77
4.2.7.	Sharpness	79
4.2.8.	Noise	81
4.2.9.	Artefacts	82
4.2.10.	Collimation	85
4.2.11.	Labelling	85
4.2.12.	Conclusion	86
	Bibliography to Chapter 4	87

5.	OUTLINE OF QC TESTS	89
6.	RADIOGRAPHER'S QUALITY CONTROL TESTS	91
6.1.	Visual inspection	93
6.1.1.	Radiological equipment	93
6.2.	Film storage	94
6.2.1.	Temperature and humidity	94
6.3.	Darkroom and film processing	95
6.3.1.	Darkroom cleanliness	95
6.3.2.	Evaluation of temperature, humidity and ventilation conditions	96
6.3.3.	White light leakage and safe lights	97
6.3.4.	Automatic processor developer temperature and other parameters	99
6.3.5.	Sensitometry	100
6.3.6.	Artefact detection	104
6.3.7.	Transition between film emulsion numbers	105
6.4.	Imaging system	106
6.4.1.	Cleaning of intensifying screens	106
6.4.2.	Screen film contact and light tightness of cassettes	107
6.4.3.	Uniformity of cassette sensitivity, attenuation and artefacts	109
6.5.	Automatic exposure control	111
6.5.1.	Test of system constancy	111
6.5.2.	Test of AEC thickness compensation	114
6.6.	Image quality	116
6.6.1.	Evaluation of image quality	116
6.7.	Film rejection rate	118
6.7.1.	Film rejection rate	118
7.	MEDICAL PHYSICIST'S TESTS	119
7.1.	Unit assembly	121
7.1.1.	Unit assembly evaluation	121
7.2.	Sensitometry and darkroom	122
7.2.1.	Sensitometry and darkroom conditions	122
7.3.	Radiological equipment	123
7.3.1.	Radiation leakage	123
7.3.2.	Accuracy and repeatability of kVp	125
7.3.3.	Half-value layer (HVL)	127
7.3.4.	Output repeatability and linearity	129
7.3.5.	Compression	131
7.3.6.	Evaluation of the automatic exposure control	133
7.4.	Collimation system	136
7.4.1.	Collimation system	136
7.5.	Image viewing conditions	138
7.5.1.	Luminance and homogeneity of the viewboxes	138
7.5.2.	Interpretation room ambient illumination	140
7.6.	Image quality	141
7.6.1.	Evaluation of image quality	141
7.6.2.	System spatial resolution	144
7.7.	Dosimetry	145
7.7.1.	Incident air kerma at the entrance surface of the phantom	145
7.7.2.	Determination of the mean glandular dose (D_G)	148

APPENDIX I:	GUIDELINES FOR THE TRAINING OF A MEDICAL PHYSICIST SPECIALIZING IN MAMMOGRAPHY	151
APPENDIX II:	MAMMOGRAPHY ROOM DESIGN	152
APPENDIX III:	DARKROOM DESIGN	155
APPENDIX IV:	TEST EQUIPMENT SPECIFICATIONS	160
APPENDIX V:	AUTOMATIC PROCESSOR DEVELOPER TIME, SPECIFIC GRAVITY, REPLENISHMENT AND pH	161
APPENDIX VI:	EXAMPLE OF A CALCULATION FOR FILM CROSSOVER PROCEDURE ...	163
APPENDIX VII:	CALCULATION OF AVERAGE FILM GRADIENT	164
APPENDIX VIII:	THE STANDARD BREAST	166
REFERENCES	167
BIBLIOGRAPHY	169
ANNEX I:	RADIOGRAPHER'S DATA COLLECTION SHEETS	171
ANNEX II:	MEDICAL PHYSICIST'S DATA COLLECTION SHEETS	195
GLOSSARY	217
CONTRIBUTORS TO DRAFTING AND REVIEW	221

1. INTRODUCTION

1.1. WHY HIGH QUALITY IN MAMMOGRAPHY IS NECESSARY

Breast cancer is the most common cancer among women worldwide and is a leading cause of cancer mortality in women. Breast cancer incidence increased 30–40% from the 1970s to the 1990s in most countries, with the most marked increases among women aged 50 years and older, although the incidence for women under 50 years is increasing. Overall, North American and northern European countries have the highest incidence rates of breast cancer; intermediate levels have been reported in western Europe, Oceania, Scandinavia and Israel; the lowest levels are observed in eastern Europe, Central America and South America, and Asia. Breast cancer incidence and mortality rates vary fourfold by geographic location between countries with the highest and lowest rates.

Mammography is an X ray examination of the breast. Its principal purpose is to facilitate the detection of breast cancer at a point earlier in its natural history than is possible by clinical examination. It has been demonstrated that routine screening with high quality mammography is effective in reducing mortality from breast cancer in women aged 40–69 years [1, 2]. In countries with such mammography screening programmes, there has been a marked decrease in breast cancer mortality over the past two decades [2]. Mammography is also useful in refining the diagnosis of breast cancer (assessment or workup) after a suspicious area in the breast has been detected and for localizing a lesion for therapy.

The radiological signs of breast cancer include mass densities which are typically slightly more attenuating of X rays than the surrounding normal tissue, small microcalcifications, asymmetry between the two breasts and architectural distortion of tissue patterns. In order to detect breast cancer accurately and at the earliest possible stage, the image must have excellent contrast to reveal mass densities and spiculated fibrous structures radiating from them that are characteristic of cancer. In addition, the spatial resolution must be excellent to reveal the calcifications, their number and their shape. The imaging system must have adequate latitude to provide this contrast and resolution over the entire breast effectively. The geometrical characteristics of the X ray unit and the positioning of the breast by the radiographer must be such that as much breast tissue as possible is included in the mammogram. Finally, the noise (signal fluctuation) of the image must be sufficiently low to reveal the subtle structures in a reliable manner and the X ray dose must be as low as reasonably achievable while being compatible with these image quality requirements.

1.2. PURPOSE

It has been well established that to achieve high quality mammography, the following elements are essential:

- (1) Well trained and experienced personnel (radiologist, radiographer, medical physicist);
- (2) Modern, well designed equipment;
- (3) Equipment in good working order;
- (4) Proper positioning and technical factors for exposure;
- (5) Appropriate image viewing conditions.

An effective quality assurance (QA) programme is necessary to ensure that all of these elements remain in place over time. The part of this programme that is concerned with the technical aspects is referred to as quality control (QC).

This publication is intended to provide a standardized framework for QC for mammography which can be used in Member States. It is intended to provide practical tests to help ensure high quality of screen film mammography. In order to be feasible in areas where resources may be limited, the tests are designed to be carried out with the simplest test equipment possible.

1.3. PHILOSOPHY

Several well established QC programmes for mammography currently exist in different jurisdictions [3, 4]. These are comprehensive and reflect resources available in those countries. The IAEA recognizes the different resources and needs of Member States and has developed specific programmes for individual areas. An example is a QC programme for mammography developed for the Latin American countries, implemented in the framework of the Regional Cooperative Agreement for the Advancement of Nuclear Science and Technology in Latin America and the Caribbean (known as ARCAL) [5] and associated national protocols [6]. The present report attempts to incorporate the most important components of these programmes in a harmonized manner to be useful to the broad range of Member States. It has been developed with the philosophy that if mammography is to be performed, it must be of high quality in order to allow the earliest detection of cancers. In some areas, both human and technological resources are limited and, therefore, this publication was developed with the concept of practicality in mind.

Recent publications indicate that digital mammography can provide equal or superior accuracy to screen film mammography [7]. Digital mammography also has potential for increased efficiency in image archiving and retrieval and the possibility to avoid the costs, complexity and waste disposal problems associated with chemical processing of film. These factors have stimulated interest in the acquisition of digital systems. This presents both opportunities and challenges to those involved in delivering mammography services. Like screen film mammography, it is apparent that the clinical performance of digital mammography depends on the proper design of the equipment and its use in an optimized manner. Therefore, one of the important challenges is to have in place in a timely fashion an appropriate framework of QA for digital mammography systems. The present publication addresses only QA issues relevant to screen film mammography. It is intended that an additional publication will address the specific issues associated with the use of digital mammography.

2. ELEMENTS OF HIGH QUALITY MAMMOGRAPHY

2.1. PERSONNEL

There are many elements that contribute to the mammography process. The experience of personnel directly and indirectly involved in the mammographic task is crucial to the final outcome. In that respect, it is essential that:

- The mammography images be acquired by experienced radiographers trained specifically in mammography (see Section 4.1, for example, on the importance of correct positioning and the use of compression);
- The images be interpreted by an appropriately trained and experienced radiologist;
- A medical physicist be available as a consultant to the facility. This may be either on a full time or part time basis, according to the needs of the facility in QA and radiation protection. (See Section 3.2.4 and Appendix I for the specific requirements that a medical physicist should meet.) The availability and qualifications of the medical physicist must be in compliance with local regulations.

2.2. EQUIPMENT

2.2.1. Mammography unit

The X ray unit must be specifically designed for mammography and include the following key features:

- X ray tube with a nominal focal spot of 0.3 mm [8];
- If magnification mammography is performed (this capability should be present on systems that are used for diagnostic mammography and not exclusively for screening), a magnification stand and a second, smaller focal spot of nominal size ≤ 0.15 mm;
- Molybdenum target. Supplementary targets composed of materials such as tungsten or rhodium may also be available;
- Tube current ≥ 80 mA for a Mo target for contact mammography and ≥ 20 mA for magnification mammography;
- Beryllium exit window;
- Beam filter of molybdenum. An additional filter composed of rhodium is highly desirable;
- Motorized compression device;
- Readout of compression thickness and force is highly desirable;
- Automatic exposure control (AEC) with a sensor whose position is adjustable;
- Fine control of optical density on AEC;
- Moving grid designed for mammography;
- Focus–film distance ≥ 60 cm;
- Buckys that can accommodate film of sizes 18 cm \times 24 cm and 24 cm \times 30 cm are desirable.

The room in which the mammography unit is sited should have a stable temperature and humidity for satisfactory operation. This may require appropriate air conditioning. More complete details for siting a mammography unit are provided in Appendix II.

2.2.2. Image receptor system

An acceptable image receptor system is characterized by having:

- Cassettes designed for mammography;
- A high resolution single mammography screen in cassette;
- Appropriate mammography film matched to be used with the selected screen.

2.2.3. Processing

Acceptable film processing requires:

- Automatic processor with digital readout of developer temperature.
- Appropriate replenishment system for mammography, matched to the film volume processed at the facility.
- Use of processing chemicals matched to the film that is used.
- Correct processing time, temperature and replenishment setting.
- Proper water quality and temperature control.
- If a darkroom is used, it must be free of light leaks and internal light sources that could fog the film. The level of ionizing radiation in the darkroom must also be within acceptable tolerances.
- Darkroom, cassettes and screens must be free from dust that could cause artefacts in the mammograms.
- Proper ventilation of processing area.
- Appropriate storage of films and chemicals.
- Management of film and chemical inventory to ensure that fresh materials are used for imaging.

See Appendix III for more details on some of the requirements mentioned.

2.2.4. Viewing conditions

Successful viewing conditions require:

- A viewbox designed for mammography with luminance ≥ 3000 cd/m²;
- Lamps in the viewbox matched for brightness and colour;
- The ability to mask edges of mammograms;
- Low ambient light in the room.

2.2.5. Quality assurance

To ensure high quality, all of the above factors are necessary; however, in addition, it is essential that a comprehensive QA programme be in place. It is also imperative that time be allocated to allow the necessary QC tests to be performed regularly, that results be carefully recorded and that corrective action be taken promptly when indicated. The basic elements of a QA programme in mammography are outlined in Section 3.

2.2.6. Regular maintenance

In addition to regular quality assurance, it is also essential that all mammographic units and associated film processors undergo regular maintenance consistent with best practice or recommendations from the manufacturers.

3. BASIC PRINCIPLES OF QUALITY ASSURANCE IN MAMMOGRAPHY

3.1. QUALITY ASSURANCE ACTIVITIES

A QA programme in diagnostic radiology, as defined by WHO [9], is an organized effort by the staff operating a facility to ensure that the diagnostic images produced are of sufficiently high quality so that they consistently provide adequate diagnostic information at the lowest possible cost and with the least possible exposure of the patient to radiation. Registrants and licensees shall establish a comprehensive QA programme for medical diagnosis with the participation of appropriate medical physicists, taking into account the principles established by WHO [9].

QA programmes for medical exposures shall include:

- (1) Measurements of the physical parameters of the radiation generators and imaging devices at the time of commissioning and periodically thereafter.
- (2) Verification of the appropriate physical and clinical factors used in patient diagnosis (or treatment).
- (3) Written records of relevant procedures and results. This includes a manual that defines clear lines of responsibility, outlines the individual QC tests performed, gives the test frequencies, is useful for staff training, facilitates audit of a service and helps to keep information within the service.
- (4) Verification of the appropriate calibration and conditions of operation of the dosimetry and monitoring equipment.
- (5) Optimization of clinical protocols and equipment operation to achieve the aims of QA as stated previously.
- (6) Regular and independent quality audit reviews of the QA programme.

QA programmes are designed to ensure that the radiology equipment and staff procedures yield the desired information. They include:

- (1) Administrative procedures or management actions designed to verify that:
 - QC tests are performed properly and according to a planned timetable;
 - Results of these tests are evaluated promptly and accurately;
 - The necessary corrective measures are taken in response to these results;
 - There is an appropriate assignment of responsibility for QA actions;
 - Standards of quality are established for equipment in the facility;
 - Adequate training is provided;
 - Appropriate equipment for each examination is selected, including the writing of adequate equipment specifications.
- (2) Acceptance testing and commissioning (see Fig. 1):
 - Acceptance tests are those performed to verify that the purchase specifications have been met by the vendor. These tests are often performed by the company installing the equipment, under the supervision of the medical physicist [10] or, preferably, independently by the medical physicist.
 - Commissioning tests are those undertaken at the time the equipment is put into service; they are used to establish baseline levels of performance and are performed by the medical physicist.
 - To a large extent, these tests overlap. This publication primarily describes tests that form a comprehensive, ongoing QC programme for mammography, but it is recognized that it is necessary to ensure that the equipment as delivered conforms to specified standards and that appropriate initial baseline values are established and used to ensure the maintenance for the quality of the equipment throughout its service life. These acceptance and commissioning tests are included in this publication and are indicated as such. During acceptance testing, a qualified person should check the electrical and mechanical safety of any new installation.




FIG. 1. Life cycle of a piece of equipment.

- (3) QC tests (also classified as either constancy or status tests by the IEC) are used to test the components of the radiological system and verify that the equipment is operating satisfactorily.
- (4) Verification of QC equipment and material.
- (5) Follow-up of any corrective actions proposed:
 - It is important that routine QC testing be properly performed in the mammography facility and that results be documented thoroughly and carefully. It is equally important that problems and potential problems be clearly documented and communicated to the facility in a timely manner and that the medical physicist be assured that the receiving party has received and understood the supplied information. This is especially the case when safety concerns are raised.
 - The reporting structure in the facility should be understood by the medical physicist, who should ideally report problems to an individual who is empowered to call in service personnel and, if necessary, who can ensure that the equipment is not used until the problems are corrected. The medical physicist may be asked to explain the problems to service personnel and share test results with them. The medical physicist and the representative from the facility should work together to ensure that the problems have been appropriately corrected.
- (6) Education and training of staff, including the radiologist, radiographer and medical physicist: each must meet a minimum level of qualification.
- (7) Continuing education: Each team member must undertake sufficient continuing education to ensure that they are up to date on new techniques and that they are refreshed relative to their basic knowledge, for example, of radiation safety.
- (8) Experience: To ensure proficiency on an annual basis, the radiologist must read a sufficient number of cases; the radiographer must do a minimum number of cases; and the medical physicist must perform a sufficient number of acceptance tests and carry out routine QC testing on a sufficient number of mammography units.

3.2. ROLES AND RESPONSIBILITIES

3.2.1. The owner

The owner has specific responsibilities to ensure that all regulatory and/or licensing requirements are met. Further, the owner must ensure that all radiologists, radiographers, medical physicists and other personnel who work at the facility are appropriately qualified and trained, and meet all continuing education and experience requirements.

It is the responsibility of the owner to ensure that a QA programme is in place encompassing all aspects of the mammography imaging process. The specific tasks within that programme may be delegated to appropriate staff who may have more expertise to carry out those tasks. Notwithstanding the above delegation, it remains the ultimate responsibility of the owner that the elements of the QA programme are fulfilled.

A lead interpreting physician (LIP), usually a mammography radiologist, should be identified by the mammography facility to have the specific responsibility of ensuring that all required QA activities are performed.

3.2.2. Lead interpreting physician (radiologist)

Although it is recognized that the LIP will delegate many of the following tasks, the LIP still has the responsibilities for:

- (1) Ensuring that the technical personnel and/or radiographers have adequate training and continuous education courses in mammography.
- (2) Motivating, supervising and managing all the aspects related to the QC programme in the area of mammography.
- (3) Providing an orientation programme for radiographers based on a carefully established procedures manual.
- (4) Selecting a single radiographer to be the primary QC radiographer to perform the prescribed QC tests and oversee those that have been delegated to other individuals.
- (5) Ensuring the availability of the equipment and necessary materials for implementation of the QC tests.
- (6) Arranging staffing and scheduling so that adequate time is available to carry out the QC tests and to record and interpret the results.
- (7) Ensuring that a medical physicist is available to oversee the equipment related QC programme and to perform the medical physicist's tests.
- (8) Reviewing the radiographer's test results at least every three months, or more frequently if consistency has not yet been achieved; and reviewing the medical physicist's test results annually, or more frequently when needed.
- (9) Overseeing, or designating an individual to oversee, the radiation protection programme for employees, patients and other individuals in the surrounding area.
- (10) Ensuring that records concerning employee qualifications, mammography technique and procedures, infection control procedures, QC, safety and protection are properly maintained and updated.
- (11) Providing feedback continually, both positive and negative, to the technical personnel and/or radiographers, on the image quality and the QC procedures.
- (12) Verifying the percentage of rejected films performed by the radiographers and ensuring that if it exceeds the specified limit, appropriate corrective action is implemented.

3.2.3. Radiographer (technologist)

The responsibilities of the radiographer include:

- (1) Ensuring that the QC tests are performed, interpreted and recorded appropriately. This is best achieved when one radiographer assumes overall responsibility for QC matters and is able to train others to assist in QC activities.

- (2) Recording imaging problems.
- (3) Undertaking additional continuous education courses in mammography.

3.2.4. Medical physicist

The medical physicist is a person trained in medical physics and certified as a medical physicist according to the applicable programme in the Member State, if such a programme exists. Guidelines for the training of a medical physicist are given in Appendix I. If it is not possible to have a medical physicist carry out the specified medical physicist's tests in this document, then such tests may be delegated to a trained radiographer or to a service person employed by the vendor of the mammography equipment.

The minimum requirements for an individual who is delegated to carry out the annual medical physicist's QC tests described in this publication are:

- (1) Training in radiation safety;
- (2) Training in the physics of mammography;
- (3) Practical training in the testing of mammography equipment.

The responsibilities of the medical physicist include:

- (1) Advising the facility on the safe and effective use of X rays for mammography. This includes image quality and radiation protection of the patient and personnel.
- (2) Advising the facility on equipment for mammography.
- (3) Conducting tests to ensure the safety and proper performance of equipment used in mammography. These tests include acceptance, commissioning and routine QC tests.
- (4) Providing oversight and advice to the radiographer who carries out the radiographer's component of the QC programme.

4. CLINICAL IMAGE QUALITY

For copyright reasons, the material in Chapter 4 is available in the hard copy version only of this publication.

5. OUTLINE OF QC TESTS

QC tests are intended to verify the stability in the operation of the equipment or elements used to acquire the mammogram. The tests have been classified into two types: essential and desirable, with respect to their importance in influencing image quality and dose. The performance of the first category of tests is considered indispensable; however, it is recommended that the tests in the second category be carried out if adequate human resources and equipment can be made available.

Test Priority

Essential refers to tests that must be done in a facility.
Desirable describes the test procedures that should be performed if feasible.

Many of the tests need to be performed very frequently (weekly and daily). Therefore, it is recommended that these tests be performed by local personnel who are present daily in the installation (technical personnel, normally radiographers). The lower frequency tests have been assigned in the majority of the cases to medical physicists and radiologists. Tolerance values for the tests are indicated. Again, these are classified into two categories: acceptable and achievable. In some cases, only the acceptable level has been defined.

Performance Standards

Acceptable indicates that performance must be within these tolerances and if it is not, the equipment should not be used.
Achievable indicates the level of performance that should be attained under favourable circumstances; this is the level at which a facility should work if it is feasible.

A facility should strive to ensure that equipment operates at the achievable level of performance, as this will produce the highest image quality and the most appropriate dose performance. It is recognized, however, that limited resources and other factors may occasionally prevent the achievable levels from being obtained. In no case should the facility continue to perform mammography if the equipment does not meet the acceptable standard of operation, because below this level the value of the procedure and/or its safety is considered unacceptable. Each test in the QC programme has a specified tolerance level for achievable and acceptable results as applicable. Should the results of a test fall outside the specified tolerance, the test should usually be repeated to confirm the result before action is taken.

Suitable minimum specifications for test equipment are provided in Appendix IV.

Table 2 in Chapter 6 and Table 4 in Chapter 7 list all the tests to be carried out by the radiographer and the medical physicist, respectively. In some cases, either could perform a particular test, so the tests are listed in both tables. In such cases, a decision must be made as to who will actually do the test and that person should consistently do the test thereafter.

6. RADIOGRAPHER'S QUALITY CONTROL TESTS

A brief description of the methodology to be undertaken when performing the radiographer's QC tests is provided in this section (Table 2). The order in which tests appear in this publication does not necessarily indicate that in which they should be performed. The preferred order will depend on various factors relating to the mammography facility as well as the evaluator's preferences, always having in mind that tests exist whose results affect the execution of others. Data collection sheets are found in Annex I and are available in electronic format.

TABLE 2. LIST OF RADIOGRAPHER'S QUALITY CONTROL TESTS

Test	Priority (E ^a , D ^b)	Suggested frequency	Tolerances
<i>Visual inspection</i>			
Visual inspection and evaluation of the mammography unit	E	Monthly	See Section 6.1
<i>Film storage</i>			
Temperature	E	Monthly	15–21°C
Humidity			40–60%
Position of film boxes and cassettes	E	Monthly	See Section 6.2.1
Film inventory	D	Monthly	Time period for inventory updating: <3 months
<i>Darkroom and film processing</i>			
Darkroom cleanliness	E	Daily	—
Temperature	E	Monthly	15–21°C
Humidity	E	Monthly	30–70%
Ventilation conditions	D	Monthly	See Sections 6.2.1 and 6.3.2
White light leakage	E	Annually	See Section 6.3.3
Safe lights	E	Annually	Rating ≥ 15 W FSL ^c <0.05 OD ^d in 2 min
Developer temperature	E	Daily	Achievable: $\pm 0.5^\circ\text{C}$ Acceptable: $\pm 1.0^\circ\text{C}$ of the manufacturer recommended value
Sensitometry	E	Daily	See Table 3 in Section 6.3.5.5
Development time, specific gravity, pH and replenishment rate		Only when problems are detected	See Appendix V
Artefact detection during processing	E	Weekly	Acceptable: no clinically significant artefacts
<i>Imaging system</i>			
Screen cleanliness	E	Weekly	See Section 6.4.1
Screen film contact	E	Semi-annually	Acceptable: spots ≤ 5 mm
Light tightness of cassettes	E	Semi-annually	Acceptable: blackening ≤ 2 mm chest wall edge, ≤ 5 mm other edges

TABLE 2. LIST OF RADIOGRAPHER'S QUALITY CONTROL TESTS (cont.)

Test	Priority (E ^a , D ^b)	Suggested frequency	Tolerances
Matching of cassette sensitivity ^c	E	Semi-annually	Achievable: maximum deviation ≤0.20 OD Acceptable: maximum deviation ≤0.30 OD
Cassettes uniformity	D	Semi-annually	Acceptable: maximum deviation ≤5% mAs
Artefacts from each cassette	E	Semi-annually	Acceptable: no clinically significant artefacts
<i>Automatic exposure control</i>			
Test of system constancy	E	Daily	Achievable: OD = OD _{target} ± 0.15 Acceptable: OD = OD _{target} ± 0.20 Acceptable: mAs within ±10% of mAs that produces OD _{target} Acceptable: no clinically significant artefacts
Compensation of the AEC ^f for different thickness	E	Monthly	Achievable: OD = OD _{target} ± 0.15 Acceptable: OD = OD _{target} ± 0.20 Acceptable: ±10% of baseline mAs
<i>Image quality</i>			
ACR ^g phantom score	D	Weekly	Acceptable: Fibres: ≥4 Microcalcifications: ≥3 Masses: ≥3
Optical density difference between disc and background	D	Weekly	Achievable: ≥0.55 OD Acceptable: ≥0.40 OD
<i>Reject analysis</i>			
Reject films analysis	E	Quarterly	Achievable: ≤3% Acceptable: ≤8%

^a E: Essential, basic requirement.

^b D: Desirable, recommended.

^c FSL: Fog due to the safety light.

^d OD: Optical density.

^e This includes speed of screens and cassette attenuation.

^f AEC: Automatic exposure control.

^g ACR: American College of Radiology.

6.1. VISUAL INSPECTION

6.1.1. Radiological equipment

6.1.1.1. Scope

- Objective: To verify the mechanical and electrical operation of the mammography unit;
- References [3, 11];
- Frequency: Monthly.

6.1.1.2. Instrumentation

- (1) Spirit level.
- (2) Tape measure.

6.1.1.3. Methodology

Execute the visual verification of the operation of the X ray machine following the checklist shown on the data collection sheet.

6.1.1.4. Recommendations and corrective actions

- (1) If some of the movements or control scales are not functioning correctly, the technical service should be called for immediate repair.
- (2) If the generator or the compression is not functioning, mammography should not be performed until the problem has been corrected and proper performance of the system has been verified.
- (3) If some of the components related to unit electrical safety show problems of electric discharges or the high voltage electrical wiring is damaged, the technical service should be called for immediate repair.

6.2. FILM STORAGE

6.2.1. Temperature and humidity

6.2.1.1. Scope

- Objectives: To verify temperature and humidity of the storage place for films and chemicals; to check the positioning and organization of film boxes, cassettes and chemicals; to verify proper storage and utilization order of the film and chemicals;
- References [3, 5, 11, 12];
- Minimum frequency: Monthly and after changes such as renovation and relocation.

6.2.1.2. Instrumentation

- (1) Thermometer.
- (2) Hygrometer.

6.2.1.3. Methodology

- (1) Measure the temperature and humidity level at the storage place.
- (2) Observe the positioning and organization of film boxes.
- (3) Review the film inventory.

6.2.1.4. Interpretation of results and conclusions

Tolerances:²

- (1) Temperature: Manufacturers' recommendations (15–21°C).
- (2) Humidity: 30–70%.
- (3) Cassettes should be stored vertically.
- (4) Film boxes should be located in a vertical position (never horizontally in order to avoid films being damaged or suffering artefacts due to the weight of film boxes placed one on top of another) and organized in chronological order in accordance with their expiration date. This placement system facilitates control of the inventory and ensures that the oldest boxes are used first.
- (5) Time period for inventory updating: <3 months.

6.2.1.5. Recommendations and corrective actions

- (1) If positioning and organization of film boxes, cassettes and chemicals is not appropriate, communicate this to the responsible person for the service.
- (2) If the environmental conditions do not fulfil the specifications of the film and chemical manufacturers, investigate the causes and inform the person responsible for undertaking corrective measures.

² For films and chemicals, the storage temperatures should be those recommended by the manufacturers.

6.3. DARKROOM AND FILM PROCESSING

6.3.1. Darkroom cleanliness

6.3.1.1. Scope

- Objective: To maintain cleanliness in the darkroom in order to minimize possible artefacts on the X ray films;
- References [3, 5, 12];
- Minimum frequency: Every day before beginning the working day.

6.3.1.2. Instrumentation

- (1) Bucket.
- (2) Lint-free cloths.

6.3.1.3. Methodology

- (1) Once a week, clean or vacuum the grills of the air ventilation channels, safety lights and walls.
- (2) Clean all the work surfaces with a wet cloth. Clean the tray that receives the film in the processor.
- (3) Ensure that the darkroom floor is clean.

6.3.1.4. Interpretation of results and conclusions

The dirt present in the darkroom affects the quality of the mammographic images, as it is introduced into the cassettes when they are handled. Such dirt becomes affixed in the form of powder specks on the intensifying screens, producing artefacts (spots of variable size) which can be observed after developing the films. Such spots can interfere with the detail of the mammographic image of diagnostic interest. Accordingly, part of the evaluation of darkroom cleanliness is to observe whether these artefacts exist in the developed films.

6.3.1.5. Recommendations and corrective actions

- (1) It is advisable for darkroom walls to be covered with a paint that does not produce reflections, preferably a dull oil paint.
- (2) The ingestion of food and beverages, and smoking must not be allowed in the darkroom.
- (3) It is necessary to ensure that hands are clean and dry at all times.
- (4) It is advisable to eliminate any object that contributes to dust accumulations.
- (5) False ceilings (e.g. suspended soft tiles) should be avoided, as they allow the buildup of dust and dirt which can easily be dislodged and fall into the processor, into open cassettes or onto films.
- (6) It is recommended that the air supplied to the darkroom be filtered.
- (7) Any of these actions should be documented in the comments section on the data collection sheet.

6.3.2. Evaluation of temperature, humidity and ventilation conditions

6.3.2.1. Scope

- Objective: To verify temperature, humidity and ventilation conditions in the darkroom;
- References [3, 5, 6, 11, 12];
- Minimal frequency: Monthly and after changes.

6.3.2.2. Instrumentation

- (1) Thermometer.
- (2) Hygrometer (if available).

6.3.2.3. Methodology

- (1) Measure the temperature and the humidity level of the darkroom. Note the values on the data collection sheet.
- (2) Check whether a strong odour of processing chemicals is present. Confirm the existence and operation of the air circulation system (pay special attention to the air extractor if there is one).
- (3) Evaluate the exhaust fan operation. 'Flutter test' (see also Appendix III.7, item (8) in the list of darkroom ventilation conditions:

“...with the processor *turned off*, a piece of tissue is to be held near the gap in the exhaust duct. The tissue should be seen to be drawn towards the gap. If there is not noticeable movement in the tissue, the flexible duct is to be removed and the negative pressure measured 25 cm into the opening of the duct, using a manometer.”).

6.3.2.4. Interpretation of results and conclusions

Tolerances:

- (1) Temperature: 15–21°C.
- (2) Humidity: 30–70%.
- (3) There should be no perceived odour of developing chemicals in the darkroom (this is a sign of deficient air circulation).
- (4) The 'flutter test' fails if there is inadequate suction, or if the tissue blows away from the gap.

6.3.2.5. Recommendations and corrective actions

If the environmental and ventilation conditions are outside of the tolerable range, as seen from the 'flutter test', the causes should be investigated and the service responsible informed about the actions that should be carried out.

6.3.3. White light leakage and safe lights

6.3.3.1. Scope

- Objective: To confirm that safe lights and possible leakage of white light in the darkroom do not produce fog on the mammographic films;
- References [3, 11–13];
- Minimum frequency: Annually unless problems are found. The test should be repeated when the safe light bulbs or filters are changed.

6.3.3.2. Instrumentation

- (1) Mammographic unit.
- (2) Phantom (45 mm of PMMA or ACR phantom).
- (3) Densitometer.
- (4) Mammographic film.
- (5) Opaque paper.
- (6) Cassettes.
- (7) Stopwatch.

6.3.3.3. Methodology

- (1) Ensure that filters, safe light bulb power, and the distance between filter and work surfaces are the recommended ones, and that the lamp is directed to the ceiling.
- (2) Turn off all the lights in the darkroom and wait for 5 min in order to adapt to the darkness.
- (3) Observe if there is any white light leaking around the doors, the processor, film exchange boxes, extractors or the ceiling. If there is any white light leakage, correct it before continuing.
- (4) In complete darkness, load a cassette with a film.³
- (5) Place the cassette in the cassette holder of the mammography unit.
- (6) Place the phantom on the image receptor, aligned with the chest wall side and centred between the two lateral sides.
- (7) Lower the compression paddle until it is in contact with the phantom.
- (8) Confirm the position of the automatic exposure control device: it should be placed in the centre of the phantom.
- (9) Make an exposure with the automatic exposure control or select those factors used to obtain the image of an average breast (of 45 mm of thickness). If a phantom is not available, expose the film in order to obtain an optical density in the range of the target optical density 1.5–1.9 optical density (OD).
- (10) With the darkroom in complete darkness, take out the film from the cassette and place it on the work surface with the emulsion side upwards. Cover half of the film with the opaque paper, placing it perpendicularly to the corresponding side to the chest wall edge of the film.
- (11) Turn on the safety lights and wait 2 min.
- (12) Process the film.
- (13) Measure the optical density of the image in the side that has been covered with the opaque paper (if using the ACR phantom, avoid measuring the density in the location of the phantom details, such as specks, masses and fibres).
- (14) Measure the optical density at a point adjacent to the previous one, but in an area that is in the part directly exposed to the safety lights.
- (15) Take note of both optical densities. This difference is used to evaluate the fog produced by the safe lights (FSL).
- (16) Record the FSL value on the data collection sheet.

³ If more than one type of film is used, this test should be repeated for each type of film.

- (17) Document the execution of this test on the data collection sheet.
- (18) If not previously done, use a permanent marker to write the installation date on the safe light filters.

6.3.3.4. *Interpretation of the results and conclusions*

Tolerances:

- (1) Bulb power: ≤ 15 W;
- (2) FSL: ≤ 0.05 in 2 min at work surface (i.e. where the cassettes are handled);
- (3) White light leakage should not be observed.

6.3.3.5. *Recommendations and corrective actions*

If FSL exceeds the tolerance, corrective measures should be implemented immediately. If the problem is not solved, the cassettes should be loaded only in complete darkness. The leading causes of fog are:

- Incorrect or expired filters;
- Cracks in the filters or in the lamp housings;
- Safe lights too near to the work surfaces;
- Incorrect bulb power (higher than 15 W);
- Indicator lights in the processor, clocks, etc.;
- Entry of light around the doors, processor or film exchange boxes;
- Light leakage in the ceiling (especially tile ceilings) or through ventilation ducts.

See also Appendix III for more complete advice about darkroom design.

6.3.4. Automatic processor developer temperature and other parameters

6.3.4.1. Scope

- Objective: To verify the developer temperature;⁴
- References [3, 4, 12, 14];
- Minimal frequency: Daily.

6.3.4.2. Instrumentation

Digital thermometer: $\pm 0.1^{\circ}\text{C}$ (never use mercury). If the processor is equipped with built-in digital temperature readout, this can be used provided that its accuracy has been confirmed in the past year.

6.3.4.3. Methodology

- (1) After the processor is turned on, wait enough time for the developer temperature to become stabilized.
- (2) If an external thermometer is used for this test, introduce the thermometer into the developer tank. Record the reading and the time on the data collection sheet.

6.3.4.4. Interpretation of results and conclusions

Tolerances:

Achievable: $\pm 0.5^{\circ}\text{C}$;

Acceptable: $\pm 1.0^{\circ}\text{C}$ or the value indicated by the film manufacturer. Check that the temperature indicator is working appropriately.

6.3.4.5. Recommendations and corrective actions

If the temperature is found to be outside of tolerances, contact the service representative.

⁴ Note: For other processor related measurement tests, such as development time, specific gravity, replenishment rate and pH, the services of a specialized technician are recommended. If this is not possible, some test procedures are suggested in Appendix III.

6.3.5. Sensitometry

6.3.5.1. Scope

- Objective: To verify that the processor is working in a stable manner;
- References [3–5, 11–17];
- Minimal frequency: Daily and after changes.

6.3.5.2. Instrumentation

- (1) Sensitometer.
- (2) Mammography film.
- (3) Densitometer. If a densitometer/computer system capable of automatically determining sensitometry information is available, such as film gradient and speed (S) at a set density, this is highly desirable and should be used. However, the main procedures described here assume that this is not available and use speed and contrast indices to establish processor performance.
- (4) Viewbox.
- (5) Thermometer.

6.3.5.3. Methodology

Establishing the initial operating levels

It may be useful for the medical physicist to work with the radiographer in establishing initial operating levels at the time the processor is installed. These levels may have to be revised when changes are introduced in:

- Film type (see Section 6.3.7);
- Chemical products brand;
- Developer and fixer temperatures;
- Replenishing rate of the developer and fixer;
- Development time;
- Sensitometer and densitometer;
- Significant processor components;
- Film box (see Section 6.3.7).

Steps:

- (1) Confirm that the darkroom meets the adequate conditions in accordance with the tolerances defined in Section 6.3.1.
- (2) Clean as required the processor, mixing and or replenishment tanks and chemistry feed lines. (Be careful not to use cleaning agents that will leave chemical residues which will change the required pH levels.) Mix developer and fixer solutions for the replenishment tanks, ensuring that the specific gravity and pH are as specified by the manufacturer. Chemistry should be added to the processor. Note that starter solution may be required to be added to the replenisher developer to ‘restrain’ the developer.
- (3) Confirm that the displayed developer temperature is consistent with that recommended by the manufacturer for the films in use.
- (4) Confirm that the replenishment rate for the chemicals is as specified by the manufacturer.

- (5) Expose the film on the emulsion side using the sensitometer⁵ with the control set to green light (for orthochromatic films).
- (6) Use the densitometer to measure the optical densities of each step on the test film. Make the measurements in the centre of each step.
- (7) Repeat the previous action for five consecutive days. The sensitometry should be carried out at the same hour of the day and the film should be placed on the input tray of the processor in the same orientation (i.e. short axis of the film parallel to the feed direction and emulsion side uppermost).
- (8) Measure with the densitometer the densities of each step of the five test films obtained. Carry out the measurement in the centre of each step.
- (9) Obtain the average values for each step using the reading of the five initial patterns.
- (10) Determine the step number with average density closest to (but not lower than) 1.20 OD. Identify this step as the speed index step. Note the number of the step and the value of its average density on the data collection sheet. This OD is the initial operating level (IOL) for S. (In daily control tests, this step number will be used to measure the value of S.⁶)
- (11) Determine the step number with average density above 2.20 OD but closest to 2.20 OD. Identify it as the step of high density (HD) and note the step number and OD on the data collection sheet. The difference between the average densities of HD and S is called density difference (DD). Record DD on the data collection sheet. This difference is the IOL for DD. (In daily QC tests, these steps and their difference will be used to obtain the value of DD.) If an automatic densitometer is used, it will also provide a graph of the gradient versus OD, step number or log exposure level.
- (12) Identify as density of base plus fog (B+F) the average density value of the first step of the pattern of densities (or of any unexposed area). Note the value on the data collection sheet. This is the IOL for B+F. (In daily QC tests, this step will be used to obtain the value of B+F.)
- (13) Write the IOL values obtained in steps (10)–(12) near the central lines of the three corresponding graphs and on the data collection sheet (see example in Fig. 42).
- (14) Establish in each part of the processor control sheet, the higher and lower tolerances for S, DD and B+F (see section on results interpretation and conclusions).
- (15) Note that the IOL should not be revised on a daily basis but should be revised only when new conditions in the processor or film emulsion cause the daily DD values to consistently move away from the IOL and approach the control limit ($IOL \pm 0.15$).

6.3.5.4. *Quality control daily tests*

Sensitometry should be carried out every day at the beginning of the workday, after the processor has reached its proper operating temperature and before processing the first film:

- (1) Expose the film using the sensitometer and process it⁷ before processing any daily mammography.
- (2) Measure the densities in the steps obtained upon establishing the IOLs. Note the values of S, B+F, and DD on the data collection sheet.
- (3) Plot the values of S, DD and B+F on the processor QC sheet.
- (4) Determine if any point is outside the tolerance values. If none is, proceed to point (7).
- (5) If there are points outside the tolerances, repeat the test. If after that they are still outside the tolerances, investigate the possible causes and correct the problem. Repeat the test again in order to verify that the

⁵ Verify that the sensitometer emits light from a single side, since they usually emit light from both sides for double emulsion films used for general radiography. Also verify that the sensitometer produces light of the appropriate colour to match the sensitivity of the film.




⁶ If there are changes on the processor or the films, the S and DD values as well as the step numbers should be determined again.

⁷ Once the film has been exposed with the sensitometer, it should be processed immediately. The way of feeding the film into the processor can affect the densities of the steps. This means that it should always be introduced in the entrance tray on the same side (right or left), in the same orientation and with the emulsion side consistently either upwards or downwards.

X-RAY PROCESSING CONTROL CHART
Department of Diagnostic Radiology

Processor: Mammo Film: ABC Emulsion # 12345 Month: April Year: 2007

Day:	2	3	4	5	6	9	10	11	12	13	16	17	18	19	20	23	24	25	26	27	30	
Month:	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4



Replenishment Rate			Temperature (C)								
Date	Developer	Fixer	Date	Developer	Fixer	Date	Developer	Fixer	Date	Developer	Fixer

FIG. 42. Sample film processor QC chart. Note that the measurement on 6 April 2007 was repeated because the initial value of DD was outside the action level.

problem has been solved and record the values in the tables for processor control on the data collection sheets.

- (6) Note the cause of the problem and type of corrective action carried out on the data collection sheet.
- (7) Observe if the graphed values are indicative of any trend (three or more values that move in the same direction). If a trend exists but the points are not outside the control limits, the mammography studies can be carried out. However, it is important to determine the cause of this behaviour and not to wait until the values go outside the permitted limits.

6.3.5.5. Interpretation of results and conclusions

Tolerances: Table 3 outlines tolerances for fundamental film properties.

TABLE 3. TOLERANCES FOR FUNDAMENTAL FILM PROPERTIES

Parameter	Acceptable	Achievable
IOL B+F	≤ 0.25	≤ 0.21
B+F	$\leq \text{IOL} + 0.03$	$\leq \text{IOL} + 0.02$
S	$\text{IOL} \pm 0.15 \text{ OD}$	$\text{IOL} \pm 0.10$
DD	$\text{IOL} \pm 0.15 \text{ OD}$	$\text{IOL} \pm 0.10$

6.3.5.6. Recommendations and corrective actions

- (1) If differences for DD and S with respect to the IOLs are found below the tolerance of $\pm 0.15 \text{ OD}$ but are higher than $\pm 0.10 \text{ OD}$, the test should be repeated. If the result is the same, clinical films can be processed, but the operation of the processor should be monitored more closely; this means measuring the temperature, pH and processing time. If the cause is not any of these factors, the replenishing rate should be monitored.
- (2) If the differences are higher than the limit of $\pm 0.15 \text{ OD}$, the source of the problem should be corrected immediately and clinical films should not be processed.
- (3) If the value of B+F is higher than the IOL by more than 0.03 OD , the problem should be corrected immediately.
- (4) The causes and corrective actions carried out in each case should be noted on the data collection sheet.

6.3.6. Artefact detection

6.3.6.1. Scope

- Objective: To determine if the processor introduces artefacts;
- References [3–6, 12, 14, 18];
- Frequency: Weekly.

6.3.6.2. Instrumentation

- (1) PMMA phantom of 45 mm thickness, free of any imperfections or artefacts.
- (2) Mammography cassettes and films.
- (3) Viewbox.

6.3.6.3. Methodology

- (1) Select two cassettes that are in optimal condition and load them with mammography films.
- (2) Put a lead number on the PMMA phantom in the right upper position of the cassettes.
- (3) Expose each cassette with the automatic exposure control or choose the exposure factors to produce a phantom image with optical density slightly higher than 1.20 OD.
- (4) Develop the films under the same conditions used for clinical films. The second film should be introduced into the processor feed tray perpendicularly with respect to the first one.
- (5) Observe both films in the viewbox and identify the origin of any spots or marks.

6.3.6.4. Interpretation of results and conclusions

Tolerance: No clinically significant artefacts are present, such as grid lines, blotches, streaks, high or low density striations [14, 19].

6.3.6.5. Recommendations and corrective actions

If artefacts are located on the part of the film where the breast image would appear, the processor maintenance person should be contacted.

6.3.7. Transition between film emulsion numbers

6.3.7.1. Scope

- Objective: To determine the transitional change in processor QC parameters when the film emulsion batch is changed or when the brand or manufacturer of films changes;
- References [3–5, 14, 18];
- Frequency: As required.

6.3.7.2. Instrumentation

- (1) Sensitometer.
- (2) Densitometer.
- (3) Film from two different boxes.

6.3.7.3. Methodology

- (1) The transition should be carried out when the box of films for QC needs to be changed or when the brand or manufacturer of films changes.
- (2) The transition should be done when at least five films remain in the box.
- (3) It is necessary to be sure that the processor is operating correctly.
- (4) At the same hour of the day, expose films from the box in use and films from the new box with the sensitometer.
- (5) Calculate the averages of S, B+F and DD for the films of both boxes and note the values on the data collection sheet.
- (6) Calculate the differences among the average values of S, B+F and DD obtained from the films from the new and old boxes. Note them on the data collection sheet.
- (7) Adjust the new levels in accordance with the differences found in the previous point. This is carried out by adding the differences, with their sign (positive or negative), to the values corresponding to the operating levels of S, B+F and DD of the old films (see Appendix VI for an example of this procedure).
- (8) Record the new operating levels and their respective limits on the data collection sheet and in the graphs on the processor control sheet.
- (9) Note the lot number of the box and the day in which the transition was carried out.

6.4. IMAGING SYSTEM

6.4.1. Cleaning of intensifying screens

6.4.1.1. Scope

- Objective: To ensure that the intensifying screens of the mammography cassettes are not damaged, are dust-free and without particles or other types of dirt that can degrade or interfere in the detail of the image of diagnostic interest;
- References [3, 5, 6];
- Minimum frequency: Weekly and after changes.

6.4.1.2. Instrumentation

- (1) Non-alkaline soap solution composition and concentration as approved by the screen manufacturer.
- (2) Lint-free cloth.
- (3) Ultraviolet lamp.
- (4) Brush.

6.4.1.3. Methodology

- (1) Confirm by inspection that the exterior part of the cassette is free of dirt before opening the cassette.
- (2) Choose a clean area of the darkroom to work in.
- (3) Dry clean the screens following the recommendations of the manufacturer.
- (4) With the aid of the ultraviolet lamp, confirm that the screens are free of dust, dirt, lint, streaks and pencil marks, nail polish residue or other spots. When cleaning is indicated, the intensifying screens should be cleaned with the soap solution and dried with the cloth, leaving the cassette partially open in a vertical position.
- (5) After loading the cassettes and before use, wait at least 15 min or for a period of time that is recommended by the manufacturer.
- (6) Document these actions on the data collection sheet.

6.4.1.4. Recommendations and corrective actions

Review the cleaning of the darkroom whenever dirt on the screens is detected.

6.4.2. Screen film contact and light tightness of cassettes

6.4.2.1. Scope

The air trapped between the screen and the film is a common cause of poor contact in single-screen cassettes. Thus, before commencing mammography with patients, it is important to wait 15 min after loading the cassette. It is recommended that the institution have enough cassettes in order for the quality of the provided service not to be affected.

- Objective: To confirm that the screen film contact is uniform; confirm the light tightness of the cassettes;
- References [3–6, 15–17];
- Minimum frequency: Semi-annually and after changes.

6.4.2.2. Instrumentation

- (1) Screen film contact test tool.
- (2) Mammography film.
- (3) Screens and cassettes to be inspected.
- (4) Densitometer with an aperture of at least 2.0 mm diameter.

6.4.2.3. Methodology

This test should be carried out specifically with a screen film contact test tool designed for QC in mammography and never with one used for general radiology.

- (1) Check the cleanliness of the cassettes and the screens to be inspected. If it is necessary to perform moist cleaning of the screens, allow them to dry in air at least 15 min before continuing the test.
- (2) Load all the cassettes to be checked with films and wait 15 min to allow the release of trapped air.
- (3) Number each cassette in order to facilitate its identification.⁸
- (4) Place the numbered cassette on the breast support plate of the mammography unit.
- (5) Place the screen film contact test tool on the cassette.
- (6) Make an exposure using a manual radiographic technique (25–28 kVp) to obtain a density between 1.5–2.0 OD. This density range provides an adequate film gradient to facilitate the detection of spots due to poor contact.
- (7) Develop the film under the same conditions as clinical mammography studies.
- (8) Expose the other cassettes under the same conditions.
- (9) Observe the films in a (previously controlled) viewbox, at a distance of 1 m, looking for dark or clear areas (bad contact). If necessary, mask the film in order to avoid any excess of light that can dazzle and impede the visualization of poor contact areas.
- (10) Check for areas of blackening, most specifically near margins. Measure the size of these areas in each margin of the film. Note the results on the data collection sheet for light tightness.
- (11) Re-clean the cassettes which do not pass the test and repeat the test.
- (12) Place in the viewbox the two images obtained of the cassette that did not pass the test the first time.
- (13) Observe and compare the poor contact areas in order to see if these areas are found in the same position (if they are not in the same position, this is probably due to pieces of dust or grit and the cassettes should be cleaned again).
- (14) Document this test on the data collection sheet.

⁸ The screen should be identified in a permanent manner by marking it with an optically opaque marker in an area which receives radiation but is outside the breast. This enables detection and correction of artefacts and other problems seen on clinical images related to a specific screen or cassette.

6.4.2.4. *Interpretation of results and conclusions*

Tolerance:

- (1) Acceptable: Small dark spots indicating the areas of poor contact ≤ 5 mm in diameter; size of the blackening areas ≤ 5 mm in the margins of each edge of the film, except at the chest wall, where 2 mm is the maximum.
- (2) The appearance of darkened areas towards the centre of the films is not acceptable.

6.4.2.5. *Recommendations and corrective actions*

- (1) The cassettes that show areas of bad contact are not suitable for clinical use and should be replaced, and matching of the speed should be verified.
- (2) The cassettes, which are not light-tight, should be repaired or discarded as indicated.

6.4.3. Uniformity of cassette sensitivity, attenuation and artefacts

6.4.3.1. Scope

- Objective: To confirm uniformity in the speed and attenuation of the intensifying screens of the cassettes used in mammography; to evaluate screen and cassette related artefacts;
- References [3–6, 15, 17];
- Minimum frequency: Semi-annually and after changes.

6.4.3.2. Instrumentation

- (1) Cassettes used routinely in the mammography service.
- (2) Phantom ACR accreditation or PMMA with 45 mm of uniform thickness. For artefact evaluation, it is preferable to use a slab of material that is large enough to cover the entire active area of the image receptor.
- (3) Densitometer.
- (4) Mammography films of the same box.
- (5) Radiographic markers (numbers).

6.4.3.3. Methodology

- (1) Select all the cassettes that are to be evaluated, number them and load them with film.

Perform steps (2)–(9) for each cassette:

- (2) Record on the data collection sheet the cassette number, screen type (indicated near the edge of the screen), condition of cassette latches and the reference data of the film used for the test (manufacturer, type and lot number).
- (3) Place the phantom on the breast support of the mammography unit, ensuring that it covers the sensors of the automatic exposure control (AEC). With the ACR accreditation phantom, it is necessary to ensure that the sensors are under the central part of the phantom.
- (4) It is preferable that the exposure be carried out with the AEC. If the equipment does not have AEC, select the technique factors manually (and use the same for all cassettes). These factors should be those clinically used to obtain images of a breast of 50 mm thickness under compression. Make the exposure, and on the data collection sheet, record the factors used (e.g. density control setting, anode, filter, kVp, mAs) to make the exposure of each cassette.
- (5) Develop all the films in the processor normally used for processing the mammography studies. Ensure that the processor is working optimally.
- (6) Measure the optical density of the images at a point located 40 mm from the chest wall and laterally centred. For the ACR phantom, measure the optical density in the geometric centre of the image. Note the result on the data collection sheet.
- (7) Determine the values of mAs maximum (mAs_{max}) and minimum (mAs_{min}), and find the difference ($mAs_{max} - mAs_{min}$) between them. Note the result on the data collection sheet.
- (8) Determine the values of maximum optical density (OD_{max}) and minimum (OD_{min}), and calculate the difference between them. Note the result on the data collection sheet.
- (9) Inspect the film for each cassette for artefacts.

6.4.3.4. Interpretation of results and conclusions

Tolerances:

- (1) Maximum deviation in OD = $OD_{\max} - OD_{\min}$:
Achievable: ≤ 0.2 OD;
Acceptable: ≤ 0.3 OD.
Record the values on the data collection sheet.
- (2) Maximum deviation in mAs = $100 (mAs_{\max} - mAs_{\min}) / mAs_{\text{mean}}$: acceptable $\leq 5\%$.
- (3) There should be no significant artefacts due to damage to individual screens.

6.4.3.5. Recommendations and corrective actions

- (1) Those cassettes that do not perform within acceptable tolerances should be removed from clinical use since they will generate excessive dose values.
- (2) Screens with significant artefacts that cannot be removed by cleaning should be replaced. Note that modern screens have a limited life expectancy and should be regularly replaced.

Note: This test assumes that both the film processor and the AEC are stable and operating properly. The processor performance can be verified with sensitometry. AEC performance can be demonstrated by performing this test several times with the same cassette (see Sections 6.3.5 and 6.5.1).

6.5. AUTOMATIC EXPOSURE CONTROL

6.5.1. Test of system constancy

6.5.1.1. Scope

- Objective: To confirm that the mammographic X ray unit is producing a consistent optical density and output. The following values must first be established for each unit in conjunction with the medical physics service:
 - Target optical density (OD_{target}): This is the optical density that the unit is aiming for and it should be selected by taking into account local factors (such as the type of film in use) but should be in the range 1.5–1.9 OD.
 - Standard measurement position: This is the location on the film image at which densitometer measurements are to be made. It is important to use the same position each time and it should be chosen as a point located 40 mm from the chest wall and laterally centred.
 - Baseline mAs: This is the mAs needed to achieve the OD_{target} when carrying out the test according to the protocol below.
- References [3–6, 15, 17];
- Minimum frequency: Daily and under circumstances such as following service or repair, after moving a mobile unit and when malfunction is suspected.

6.5.1.2. Instrumentation

- (1) PMMA slab(s) with a total thickness of 45 mm. The slabs should be labelled and the same ones used each time.
- (2) A test cassette used for all of the measurements. This is a regular cassette used for clinical mammography; however, this cassette is marked so that it can be used for QC testing.
- (3) Densitometer.
- (4) Film in current use.

6.5.1.3. Methodology



- (1) Place the PMMA centrally on the breast support table against and slightly overlapping the chest wall edge.
- (2) Ensure that the AEC chamber is at the chest wall position.
- (3) Place the loaded test cassette in the cassette holder.
- (4) Select all operating parameters as used clinically for this breast thickness. These will include AEC density control setting, AEC mode, collimator and compression paddle. If kVp/target/filter are selected manually, always use the same settings for this test.
- (5) Apply sufficient compression to activate the AEC (if required). Different types of mammography machine have different AEC systems and it is recommended that local testing protocols specify the level of compression applied for this test.
- (6) Make an exposure and record all post-exposure factors, such as mAs, target material, filter material and kVp.
- (7) If processing facilities are available, process the film after carrying out sensitometry. If processing facilities are not available, the film should be returned with the others for processing and reading.
- (8) Measure the resultant optical density at the standard position and record it on the data collection sheet. It is useful to record the optical density graphically in order to identify trends.
- (9) Check the film for artefacts and other faults involving, for example, beam alignment, grid lines, scratches and processing.
- (10) Record the optical density and exposure factors (target, filter kVp and mAs) used for each exposure. Plot graphs of OD and mAs versus the date as illustrated in Fig. 43.
- (11) Record on the data collection sheet if the performance falls outside of tolerances.
- (12) Comment about the presence of any significant artefacts on the data collection sheet.

X-RAY SYSTEM AEC CONSTANCY CHART

Department of Diagnostic Radiology

Room: 1 Cassette: QC6 Month: April Year: 2007
 AEC Mode: Auto kV Density setting: "0" AEC detector position: chest
 kVp: 27 Anode/Filter: Mo/Mo Phantom: 45 mm PMMA

Day:	2	3	4	5	6	9	10	11	12	13	16	17	18	19	20	23	24	25	26	27	30
Month:	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4	4
Initials:																					



Comments

Date	Comment
12th April	Unit switched to 26 kV giving high mAs

Date	Comment

FIG. 43. Sample X ray AEC constancy QC chart.

6.5.1.4. *Interpretation of results and conclusions*

Tolerances:

- (1) Constancy of optical density:
Achievable: $OD = OD_{\text{target}} \pm 0.15$;
Acceptable: $OD = OD_{\text{target}} \pm 0.2$;
- (2) Constancy of mAs: Acceptable: $\pm 10\%$ of baseline mAs;
- (3) No clinically significant artefacts.

6.5.1.5. *Recommendations and corrective actions*

In case of non-compliance, contact the person in charge of the maintenance service.

6.5.2. Test of AEC thickness compensation

6.5.2.1. Scope

- Objective: To confirm that the AEC produces the target optical density and compensates for different thicknesses;
- References [3, 4, 6, 15, 17, 20, 21];
- Minimum frequency: Monthly.

6.5.2.2. Instrumentation

- (1) PMMA slabs to be able to produce three thicknesses of approximately 20 mm, 45 mm and 70 mm.
- (2) Radio-opaque objects for film identification.
- (3) A test cassette used for all of the measurements. This is a regular cassette used for clinical mammography; however, this cassette is marked so that it can be used for QC testing.
- (4) Densitometer.

6.5.2.3. Methodology

- (1) Select the clinically relevant automatic exposure mode. Select position '0' or 'Normal' or 'N' on the density control selector.
- (2) Use the same (test) cassette to carry out all the tests.
- (3) Place the loaded cassette in the Bucky.
- (4) Place a radio-opaque object to identify the film.
- (5) Place the thinnest PMMA slab on the breast support, with its edge aligned with the chest wall margin of the breast support and laterally centred. Confirm that the PMMA covers the AEC sensors completely. It is important that the AEC sensor be placed in the same position every time the test is performed (e.g. at the position closest to the chest wall edge).
- (6) If the mammography system is fully automatic (i.e. it automatically selects the target/filter combination and kVp), operate it in this manner. Otherwise, select the kVp and target/filter combination used clinically for a breast corresponding to this absorber. Note the mAs for this exposure.⁹ Process the film and measure the OD at the standard position.
- (7) For each of the other two thicknesses of PMMA, expose a film using the clinically relevant kVp and target/filter combination for that thickness. If automatic mode is used clinically then this mode should be utilized. A technique chart is essential in the absence of automatic mode.
- (8) Record the optical densities and exposure factors (target, filter kVp and mAs) used for each exposure. Plot graphs of OD and mAs for the three thicknesses versus the date, as illustrated in Fig. 44.
- (9) Record on the data collection sheet if the performance falls outside of tolerances.

6.5.2.4. Interpretation of results and conclusions

Tolerances:

- (1) Optical density for each thickness:
Achievable: $OD = OD_{\text{target}} \pm 0.15$;
Acceptable: $OD = OD_{\text{target}} \pm 0.2$;
- (2) Exposure (mAs) for each thickness: Acceptable: $\pm 10\%$ of baseline mAs.

⁹ If it is noticed that the automatic system causes the kVp or target/filter combination to change for different exposures, change the thickness of PMMA slightly (by 5 mm) to force the system to select a single, consistent spectrum. Then repeat the test.

6.5.2.5. Recommendation and corrective action

In case of non-compliance, contact the person in charge of the maintenance service.




FIG. 44. Sample X ray AEC compensation QC chart.

6.6. IMAGE QUALITY

6.6.1. Evaluation of image quality

6.6.1.1. Scope

- Objective: To ensure that image quality is acceptable;
- Reference [3];
- Minimum frequency: Weekly.

6.6.1.2. Instrumentation

- (1) Breast phantom (see Section 7.6.1 for details).
- (2) Cassette (it is recommended that the same cassette always be used) and film.
- (3) Magnifier: 2× or greater. It is recommended that the magnifier contain a second area that provides 4× or 5× magnification.
- (4) Viewbox.

6.6.1.3. Methodology

- (1) Load the cassette with a film from the box that is in use for patients. Wait the amount of time recommended in order to have good screen film contact. (If it is the first film taken from a new box, note the emulsion number in the data collection sheet.)
- (2) Place the loaded cassette in the Bucky.
- (3) Place the phantom on the breast support aligned with the chest wall and centred laterally.
- (4) Place the PMMA disc on the phantom in such a way that it does not interfere with any phantom detail or the AEC sensor. It is recommended that it be located slightly below and between the first and the second fibre.
- (5) Lower the compression paddle until it touches the phantom.
- (6) Confirm that the AEC sensor is under the phantom.
- (7) Select the technique factors that are used in clinical practice: target, filter, kVp, grid, density control position, operation mode (semi-automatic or automatic).
- (8) Make an exposure.
- (9) Record the mAs on the corresponding data collection sheet. Check that the kVp and target/filter are as established at the baseline. If not, comment on the data collection sheet.
- (10) Process the film in the processor normally used for mammography, introducing it always with the same orientation.
- (11) View the film with the magnifier on the viewbox used clinically, placing a mask over the area of the viewbox that is not being used, and with a very low level of environmental illumination in the reading room.
- (12) Evaluate the image according to the evaluation method (see Section 7.6.1). Note the results on the data collection sheet.
- (13) Observe the image and determine if there are significant artefacts. Record in the comments section of the data collection sheet. Investigate the causes of any artefacts.

6.6.1.4. Interpretation of results and conclusions

Tolerances:

- (1) Image contrast index:
 - Achievable: ≥ 0.55 ;
 - Acceptable: ≥ 0.40 for the ODD between PMMA disc and adjacent point;

- (2) Image quality:
Acceptable: fibres: ≥ 4 ;
Microcalcifications: ≥ 3 ;
Masses: ≥ 3 .

6.6.1.5. *Recommendations and corrective actions*

- (1) If the image contrast index or image quality deteriorates over time, it will be necessary to carry out other tests (e.g. kVp, AEC, processor) to determine the source of the change.¹⁰
- (2) If there are several different types of film in use, this test should be carried out for each type.

¹⁰ Due to the subjective component associated with the observer, it is recommended that the test always be performed by the same person, using the same criteria, viewbox, magnifier lens and illumination conditions.

6.7. FILM REJECTION RATE

6.7.1. Film rejection rate

6.7.1.1. Scope

- Objective: To identify the causes of the film rejection and the most urgent needs for action in order to improve the operation of the service;
- References [3, 5, 6, 11, 22];
- Minimum frequency: Quarterly.

6.7.1.2. Instrumentation

Rejected films.

6.7.1.3. Methodology

- (1) Collect the rejected films during a period of 15 days or a month (the actual collection time will depend on the rejection rate).
- (2) Note on the data collection sheet the number of total studies conducted in the same period.
- (3) Note on the data collection sheet the person responsible (radiologist or technologist) who decided that the image or images be repeated.
- (4) Identify in each case the explicit cause or causes of the rejection and standardize them with a code. Take note of them on the data collection sheet together with the projection type (e.g. CC, MLO, Lat).
- (5) Note on the data collection sheet the incidence of each rejection cause.
- (6) Obtain the total number of films used in the same period based on a reliable source (e.g. data from the film warehouse, from the exposure or processed film counters, from the service statistics). If possible, the number of films used for QC testing and other non-patient exposures should be excluded. This will yield a more realistic rejection rate.
- (7) Calculate the rejection rate as the quotient between the number of rejected films and the total number of films used during a certain time period, expressed as a percentage.

6.7.1.4. Interpretation of results and conclusions

Tolerances:

Achievable: <3%;

Acceptable: <8%.

6.7.1.5. Recommendations and corrective actions

Analyse the reasons for the rejection of images so that remedial action may be implemented with regard to the most serious causes. However, on many occasions it will be impossible to judge whether the image system is working adequately. The film rejection rate is a global indicator of the necessity for monitoring the QC programme.

Note: Low rejection rates may in fact indicate that radiologists are accepting films that should be repeated, rather than indicating good practice.

7. MEDICAL PHYSICIST'S TESTS

The tests to be performed by the medical physicist are listed in Table 4. The order in which tests are performed does not necessarily have to be the same as that in which they appear in this publication. The preferred order will depend on various factors relating to the mammography facility and the medical physicist's preferences. It should always be considered that some test results are required before other tests can be performed effectively. In some cases, the medical physicist may also be responsible for performing some of the tests described in the radiographer's section. Data collection sheets for the medical physics tests are found in Annex II and are available in electronic format. All annual tests would also be done at acceptance. Additionally, at commissioning the baseline values for the radiographer's tests would be established. This requirement is indicated in the test descriptions themselves and is also indicated in the column labelled 'suggested frequency' of Table 4.

TABLE 4. LIST OF MEDICAL PHYSICIST'S TESTS

Test	Priority (E ^a , D ^b)	Suggested frequency	Tolerances
<i>Unit assembly evaluation</i>			
Unit assembly evaluation	E	Annually	See Section 7.1
<i>Sensitometry and darkroom</i>			
Sensitometry and darkroom	E	At commissioning and annually	See Section 7.2.1
Darkroom radiation level	D	As required	Acceptable: <20 µGy/week
<i>Radiological equipment</i>			
Radiation leakage	D	At acceptance and after changes	Acceptable: ≤1 mGy/h at 1 m
Accuracy and repeatability of the tube kVp	E	Annually	Acceptable: Accuracy: ±5% Repeatability: difference ≤5% or COV ^c ≤2%
Half-value layer (HVL)	E	Annually	See Section 7.3.3.5
Output: repeatability and linearity	E	Annually	Acceptable: Repeatability: difference ≤5% or COV ≤5%; Linearity: ±10%
Normalized output value	D	Annually	Acceptable: >30 µGy/mAs at 1 m, 28 kVp, Mo/Mo
<i>Compression</i>			
Compression force and thickness	E	Annually	See Section 7.3.5.4
<i>Automatic exposure control</i>			
Repeatability of the automatic exposure control	E	Annually	Acceptable: COV in mAs: ≤5%
Constancy of OD ^d with baseline value	E	Annually	Acceptable: OD = OD _{target} ± 0.20
Exposure time for 45 mm slab	E	Annually	Contact mammography: Achievable: t ≤ 1.5 s Acceptable: t ≤ 2 s Magnification mammography: Achievable: t ≤ 2 s Acceptable: t ≤ 3 s

TABLE 4. LIST OF MEDICAL PHYSICIST'S TESTS (cont.)

Test	Priority (E, D)	Suggested frequency	Tolerances
Compensation of the AEC ^c for different thickness and beam quality	E	Annually	Achievable: $OD = OD_{\text{target}} \pm 0.15$ Acceptable: $OD = OD_{\text{target}} \pm 0.20$
Increase of OD for each step of the density control	E	Annually	Acceptable: $\Delta^f OD = 0.1-0.2$
<i>Collimation system</i>			
Light field/radiation field coincidence	D	Annually	Achievable: $\leq 1\%$ of FFD ^g all edges
Radiation field/image receptor coincidence	E	Annually	Achievable: completely irradiate the image receptor, but does not extend beyond the shielded breast support except at the chest wall where it may extend by ≤ 5 mm (see Section 7.4.1.4) Acceptable: as above for the chest wall and within the breast support by $\leq 2\%$ FFD for the other edges
Compression paddle/breast support alignment	E	Annually	Acceptable: paddle not visible in image and edge of paddle $\leq 1\%$ of FFD beyond chest wall edge of image receptor
<i>Image viewing conditions</i>			
Luminance of the viewboxes	E	Annually	>3000 cd/m ² (nit)
Viewboxes homogeneity and colour	E	Annually	Acceptable: $<30\%$ for each viewbox and $<15\%$ between panels in a viewbox
Ambient interpretation room illumination	E	Annually	Achievable: ≤ 10 lx Acceptable: ≤ 50 lx
<i>Image quality</i>			
Target background density	E	Annually	Acceptable: $OD = OD_{\text{target}} \pm 0.20$
Optical density difference between disc and background	E	Annually	Achievable: ≥ 0.55 OD Acceptable: ≥ 0.40 OD
Phantom image quality evaluation (ACR ^h)	E	Annually	Acceptable: Fibre score ≥ 4 Speck score ≥ 3 Mass score ≥ 3
System spatial resolution	E	Annually	Achievable: ≥ 15 lp/mm Acceptable: ≥ 11 lp/mm
<i>Dosimetry</i>			
Mean glandular dose (D_G)	E	Annually	Achievable: $D_G \leq 2$ mGy Acceptable: $D_G \leq 2.5$ mGy

^a E: Essential, basic requirement.

^b D: Desirable, recommended.

^c COV: Coefficient of variation.

^d OD: Optical density.

^e AEC: Automatic exposure control.

^f Δ : Change in parameter.

^g FFD: Focus film distance.

^h American College of Radiology.

Note: The ACR phantom has been taken as an example because it is probably the one most commonly used.

7.1. UNIT ASSEMBLY

7.1.1. Unit assembly evaluation

The medical physicist should carry out an inspection of the mammography facility to ensure that the mechanical features of the unit are functioning correctly and that basic safety features are in place. The following guidelines should be followed:

- Verify that the free-standing unit is mechanically stable under normal operating conditions.
- Check that indicator lights are working properly.
- Confirm that all moving parts move smoothly, without obstructions to motion, and no obstructions hinder the full range of motions within these limits.
- Check that all locks and detents work properly by setting them and confirming that the relevant mechanical motion is inhibited.
- Check that the angulation indicators function properly.
- Confirm that the image receptor and holder assembly is free from wobble or vibration during normal operation (exposure).
- Establish that the image receptor slides smoothly into the holder assembly and is held securely regardless of gantry orientation.
- Inspect the compression paddles to confirm they are in good condition.
- Verify by direct measurements that the compression breast thickness scale (analog or digital) is accurate to within ± 5 mm and reproducible to within ± 2 mm under conditions of moderate compression. Measurements should be performed for at least the small and large image receptors and compression paddles if applicable (see also Section 7.3.5).
- Test that the automatic compression release following exposure functions correctly and can be overridden to maintain compression if required.
- Confirm, by switching the power OFF, that the compression can be released manually in the event of power failure.
- Verify that the operator shielding is adequate during exposure at the operator's position.
- Inspect the unit to confirm that neither the patient nor the operator is exposed to sharp or rough edges or other hazards.
- Confirm that appropriate, current operator technique control charts are posted.

7.1.1.1. Scope

- Objective: To evaluate mechanical functionality and safety of the mammography unit;
- References [3, 4];
- Minimum frequency: At acceptance, annually, and after changes that can alter the mechanical function of all or part of the unit.

7.1.1.2. Instrumentation

Tape measure or ruler.

7.1.1.3. Methodology

Record the results of the inspection on the data collection sheet.

7.1.1.4. Recommendations and corrective actions

Items that are hazardous, inoperative or fail to operate correctly should be repaired by the appropriate equipment manufacturer.

7.2. SENSITOMETRY AND DARKROOM

7.2.1. Sensitometry and darkroom conditions

It is essential that the medical physicist review the site sensitometric QC records, and desirable that he/she measure the full sensitometric curve of the film in clinical usage. The procedure is described in Section 6.3.5, however, all optical density steps should be measured and the sensitometric curve plotted (OD versus log exposure). Information on the exposure increments will be available in the sensitometer's manual. In addition, the gradient of the film should be calculated between net densities of 1.0 and 2.0 (see Appendix VII). Alternatively, this could be done using an automatic reading densitometer [14].

This is particularly important during commissioning when assisting in the establishment of the site QC programme. Annual test results should be compared with those from previous visits or from other sites using the same film. The purpose of this measurement is to help standardize performance among multiple facilities.

Similarly, it is essential that the physicist review the darkroom conditions and the associated radiographer tests. This involves checking temperature, humidity, ventilation conditions and radiation level survey, if required, as described in Appendix III. This assessment is essential at acceptance testing and ideally would follow previous consultation in the design phase of the facility. Annual tests would include the flutter test (see Section 6.2.1) and any other actions dictated by a review of conditions.

7.3. RADIOLOGICAL EQUIPMENT

7.3.1. Radiation leakage

7.3.1.1. Scope

- Objective: To evaluate radiation leakage of the X ray tube;
- References [3, 4];
- Minimum frequency: At acceptance, and after changes in the X ray housing that could influence radiation levels.

7.3.1.2. Instrumentation

- (1) Lead (sheet, shielded apron with lead, etc.);
- (2) Loaded cassette or radiographic films for direct exposure;
- (3) Appropriate radiation detector for measurement of leakage radiation.

7.3.1.3. Methodology

- (1) Place a lead sheet in the collimator in order to ensure that no primary radiation can escape.
- (2) Cover the head of the X ray tube with loaded cassettes (or films for direct exposure), marking them to allow subsequent identification of their position.
- (3) Make an exposure, selecting the maximum value of kVp and a reasonable value of mAs that ensures damage to the tube will be avoided.
- (4) Process the films and identify those in which the points with the greatest degree of blackening are observed. These points indicate the presence of radiation leakage.
- (5) If leakage is detected, it should be quantified as follows:
 - Measure the air kerma at the point or points at which radiation leakage has been noticed. Place an appropriate ionization chamber at a known distance close enough to obtain a reliable measurement, but no closer than 20 cm from the focal spot.
 - Note these values on the data collection sheet.

7.3.1.4. Calculation procedure

- (1) Transform the readings obtained to mGy/h at a distance of 1 m from the focus and at the maximum tube loading for 1 h, applying the following formula:

$$\dot{K}_{RF} = \frac{\text{mAs}_{\text{max}}}{\text{mAs}_0} \times K_0 \times d^2 \quad (1)$$

where \dot{K}_{RF} is the radiation leakage rate (kerma), d is the distance in metres from the focal spot to the point where the ionization chamber is located, mAs_{max} is the maximum tube loading per hour of the X ray tube specified by the manufacturer, mAs_0 is the actual tube loading used and K_0 is the reading from the radiation detector.

- (2) Record the results on the data collection sheet.

7.3.1.5. *Interpretation of results and conclusions*

Tolerance:

Acceptable: $\dot{K}_{RF} < 1$ mGy in 1 h, at a distance of 1 m from the focus, and at the maximum output rate (power), as averaged over an area that does not exceed 100 cm².

7.3.1.6. *Recommendations and corrective actions*

If values exceed tolerance, it will be necessary to contact the equipment manufacturer to correct the problem.

7.3.2. Accuracy and repeatability of kVp

7.3.2.1. Scope

- Objective: To verify the accuracy and repeatability of the kVp;
- References [3, 15, 20];
- Minimum frequency: At acceptance, annually (essential), or semi-annually (desirable), and after changes that can alter the kVp of the X ray tube.

7.3.2.2. Instrumentation

A kVp meter (non-invasive type).

7.3.2.3. Methodology

- (1) Check that the line voltage of the electricity supply connected to the X ray equipment is stable, if applicable. If not, contact the service engineer.
- (2) Select a number of kVp settings commonly used in the clinical practice range (e.g. 25–35 kVp).
- (3) Note on the data collection sheet the selected values and size of the focuses used.
- (4) Normally, the compression paddle can be left in place.¹¹
- (5) Position the measuring instrument on the breast support and centre it in the radiation field.
- (6) If it is possible, adjust the field size to the one specified for the measuring instrument.
- (7) In manual mode (without AEC), carry out two exposures at 28 kVp (a value of mAs in the range of 30–40 can be used; or the value recommended by the manufacturer of the instrument) and note the measured values and the mAs used.
- (8) Use the calculation procedure to calculate the percentage difference between the two measurements. If the two measurements do not differ by more than 5%, go to step (10).
- (9) If the two measurements differ by more than 5%, take another three measurements at the same setting.
- (10) Carry out a single exposure for the other kVp settings selected and note the values measured on the data collection sheet.

7.3.2.4. Calculation procedure

Repeatability:

- (1) Calculate the measurement percentage difference using the following formula:

Difference (%) = (maximum measurement – minimum measurement)/minimum measurement

- (2) If five measurements are made for the same kVp (28 kVp), calculate the coefficient of variation (COV) of tube voltage.

¹¹ Consult the specifications of the kVp meter to check if it is calibrated for measuring with the breast compression paddle interposed in the field.

Accuracy:

- (1) For each selected kVp, determine the percentage deviation between the nominal value and the measured kVp (for 28 kVp, take the first value or the average value if more than two readings are made) in accordance with the following equation:

$$\text{Deviation (\%)} = 100 \frac{kVp_{\text{nom}} - kVp_{\text{measured}}}{kVp_{\text{nom}}} \quad (2)$$

where kVp_{nom} is the value indicated on the equipment and kVp_{measured} is the measured value. This percentage deviation may be taken as a measure of the accuracy.

- (2) Note on the data collection sheet the values obtained.

7.3.2.5. *Interpretation of results and conclusions*

Tolerances:

- (1) Acceptable:
Accuracy: $\pm 5\%$;
Repeatability: difference $\leq 5\%$ or COV $\leq 2\%$;
- (2) If there is one or more of the measured kVp values outside the tolerance, note it as unacceptable in the summary data collection sheet.

7.3.2.6. *Recommendations and corrective actions*

In case of non-conformance, the equipment should be reviewed by the person in charge of maintenance, within a period of 30 days.

7.3.3. Half-value layer (HVL)

7.3.3.1. Scope

- Objective: To confirm that the total filtration of the X ray beam is in agreement with the minimum requirements of the national and international standards;
- References [3, 4, 15, 20, 23];
- Minimum frequency: At acceptance, annually, and after changes or maintenance in the housing and/or collimation system.

7.3.3.2. Instrumentation

- (1) Appropriate detector and electrometer for mammography (see Appendix IV);
- (2) Aluminium filters;
- (3) Measuring tape.

7.3.3.3. Methodology

- (1) Select the manual mode of operation, a clinically used kVp and target filter combination and an exposure that provides a reliable dosimeter reading (40–80 mAs).
- (2) Place the dosimeter¹² at a height of 45 mm over the breast support, laterally centred and 40 mm from the chest wall, so that the sensitive volume of the chamber remains completely within the radiation field (Fig. 45).




FIG. 45. Set-up for HVL measurement.

¹² The ideal geometry for measuring the HVL should be such that the beam is narrow, in order to reduce the quantity of scattered photons that affect the dosimeter. This is obtained with a lead plate put between the aluminium filters and the dosimeter. The lead plate should have a hole to allow the radiation to pass; this ensures that the beam is narrow and that it fully covers the dosimeter.

- (3) If it is possible, collimate the radiation field to cover the sensitive volume of the chamber.
- (4) Place the compression paddle¹³ approximately halfway between the focus and the dosimeter.
- (5) Carry out an exposure and record the reading on the data collection sheet.
- (6) Put 0.3 mm of aluminium on the compression paddle, totally covering the active volume of the chamber, and make an exposure with the same parameters. Check that the reading is higher than half of the reading without the filter. If it is not, use a thinner aluminium thickness.
- (7) Add 0.1 mm of Al (total 0.4 mm Al) and repeat the previous step. Check that the reading is smaller than half of the reading without the filter. Otherwise, add more Al until the reading falls below half of the reading without the filter.
- (8) Remove all the filters and repeat the exposure again. Take note of the reading.
- (9) Repeat this procedure for other clinically relevant filter/anode combinations.

7.3.3.4. Procedure of calculation

- (1) Calculate the value of the HVL, based on the following expression:

$$HVL = \frac{t_2 \ln[2M_1/M_0] - t_1 \ln[2M_2/M_0]}{\ln[M_1/M_2]} \quad (3)$$

where:

t_1 and t_2 are the thicknesses (in mm) of the filters used;

M_1 and M_2 are the average values of the readings measured in steps (6) and (7);

M_0 is the average value of the reading measured without any added filter.

- (2) Note the calculated HVL on the data collection sheet.

7.3.3.5. Interpretation of results and conclusions

Tolerance:

Acceptable: $kVp/100 + 0.03 \leq HVL \leq kVp/100 + C$

where:

$C = 0.12$ for Mo/Mo

0.19 for Mo/Rh

0.22 for Rh/Rh

0.30 for W/Rh

and kVp is the measured value for the nominal kVp selected.

Record on the data collection sheet whether the value of the HVL is acceptable.¹⁴

7.3.3.6. Recommendations and corrective actions

If the HVL is very low or very high, suspend use of the unit immediately, until the problem is solved.

¹³ For the calculation of mean glandular dose it is necessary to measure the HVL with the compression paddle in the beam. If it is desired to compare HVL with the manufacturer's specification, which follows IEC standards, the measurement of HVL would be done without the compression paddle.

¹⁴ If this tolerance is met, it is generally the case that regulatory requirements for total beam filtration without compression paddle will also be satisfied.

7.3.4. Output repeatability and linearity

7.3.4.1. Scope

- Objective: To evaluate the repeatability of the air kerma for a given mAs, the linearity with the mAs and the normalized output value ($\mu\text{Gy}/\text{mAs}$ at 1 m);
- References [3, 4, 15, 20, 23];
- Minimum frequency: At acceptance, annually, and after changes.

7.3.4.2. Instrumentation

- (1) Appropriate dosimeter and electrometer for mammography.
- (2) Measuring tape.
- (3) Thermometer (if an ionization chamber requiring temperature correction is used).
- (4) Barometer (if an ionization chamber requiring pressure correction is used).

7.3.4.3. Methodology

- (1) Place the dosimeter at a known height above the breast support (such as 45 mm or 100 mm), laterally centred and 40 mm from the chest wall, so that the sensitive volume of the dosimeter remains completely irradiated. Measure and note the focus–dosimeter distance.
- (2) Select manual mode, the Mo/Mo target filter combination, and 28 kVp.
- (3) Select three values of mAs from among the ones commonly used in the clinical practice.
- (4) Fix the first selected value of mAs and carry out two exposures. Record the readings on the data collection sheet.
- (5) Calculate the percentage difference between the two measurements. If the two measurements do not differ by more than 5%, go to step (7).
- (6) If the two measurements differ by more than 5%, take another three measurements at the same setting.
- (7) Take a single exposure for each of the other values of mAs.
- (8) If there are other filter target combinations, the output is also measured for those combinations, selecting the highest kVp most commonly used in the clinical practice. In this case, it is necessary to check that the dosimeter is calibrated for those beam qualities.

7.3.4.4. Calculation procedure

Repeatability of the kerma (mGy):

- (1) Calculate the measurement percentage difference using the following formula:

$$\text{Difference (\%)} = (\text{maximum measurement} - \text{minimum measurement}) / \text{minimum measurement}$$

- (2) If five measurements are made for the same output factors, calculate the coefficient of variation (COV) of the air kerma.

Linearity of the output ($\mu\text{Gy}/\text{mAs}$):

- (1) For each mAs selected, calculate the average value of the obtained readings of air kerma and record it on the data collection sheet.
- (2) Calculate the output, Y, by dividing each average air kerma value obtained by the corresponding mAs and record the results.
- (3) Take two output values (Y_1 and Y_2) at consecutive mAs settings and calculate the linearity (L) using $L = 100 * (Y_1 - Y_2) / (Y_1 + Y_2)$.
- (4) Note the result on the data collection sheet.

Normalized output ($\mu\text{Gy/mAs}$ at 1 m):

- (1) Calculate the average value of the output for the different values of mAs.
- (2) Apply an inverse square law correction to the result to obtain the output at a distance of 1.0 m from the focus. Multiply by the pressure and temperature correction factor if necessary. Record this value in $\mu\text{Gy/mAs}$ at 1 m.
- (3) Note the results on the data collection sheet.

7.3.4.5. *Interpretation of results and conclusions*

Tolerances:

- (1) Acceptable:
 - Repeatability: difference $\leq 5\%$ or COV $\leq 5\%$;
 - Linearity: $< 10\%$;
 - Normalized output: $> 30 \mu\text{Gy/mAs}$ at 1 m (28 kVp, Mo/Mo);
- (2) Indicate in the data collection sheet whether the tolerances are met.

7.3.4.6. *Recommendations and corrective actions*

If the system is operating outside the tolerance range, contact the person in charge of the maintenance service.

7.3.5. Compression

7.3.5.1. Scope

- Objectives: To check that the mammography system provides an adequate compression in manual and automatic mode; to check the accuracy (or deviation) of the indicator of the compression force when it is present on the equipment; to check the accuracy of the compression thickness indicator;
- References [3, 4, 15];
- Frequency: At acceptance, annually, and if a reduction in breast compression is observed.

7.3.5.2. Instrumentation

- (1) Platform scales (conventional, analog type, non-digital).
- (2) Bath towels (cloths), blocks of rubber foam or tennis ball.
- (3) Slabs of PMMA used for AEC testing.

7.3.5.3. Methodology

Power compression mode:

- (1) Place a bath towel on the Bucky and the platform scale over it. Centre the scale directly under the compression paddle.
- (2) Place one or more towels (or block of rubber foam) on the scale in order to protect the compression paddle but such that it does not obscure the reading on the scale.
- (3) Activate the compression paddle so that it operates and stops at the maximum available force.
- (4) Read and write the value of the compression force on the data collection sheet.
- (5) Release the compression.

Manual mode:

- (1) Using the manual compression mode, move the compression paddle until it stops.
- (2) Read and record the compression force on the data collection sheet.
- (3) Release the compression.

Compression thickness:

- (1) Align a 20 mm thickness of PMMA with the chest wall edge of the breast support platform. Ideally, 18 × 24 cm slabs of PMMA should be used to prevent deformation of the plate and reduce measurement inaccuracies on the indication of thickness due to the tilt angle of the plate. If semicircular slabs are used, they may be aligned perpendicular to the chest wall to minimize both of these effects.
- (2) Ensure that the compression plate is in contact with the PMMA.
- (3) Apply a compression force of between 50 N and 100 N. Once the applied force for this measurement has been established, the same force is used for all subsequent measurements.
- (4) Note the distance from the compression plate to the breast support platform on the midline at the chest wall, which is equal to the thickness of the slabs.
- (5) Compare measured and displayed values.
- (6) Repeat procedure for at least 45 mm and 70 mm thicknesses of PMMA.
- (7) If necessary, repeat the measurements in magnification mode.

7.3.5.4. *Interpretation of results and conclusions*

Tolerances:

- (1) Achievable: Displayed thickness within ± 5 mm of phantom thickness.¹⁵
- (2) Acceptable: The maximum compression force for powered compression should be no less than 150 N and no more than 200 N. Displayed value accuracy: ± 20 N. Maximum manual compression force: < 300 N.

7.3.5.5. *Corrective actions*

If the measured values from these tests are outside tolerance, the compression device should be calibrated by a qualified service engineer.

¹⁵ The displayed compressed breast thickness may be used to choose the technique factors so that it is important that this level of accuracy be achieved.

7.3.6. Evaluation of the automatic exposure control

Each facility should establish a target optical density for mammography at commissioning. This is the optical density that the radiologist considers optimal for imaging. It is helpful for the medical physicist to provide advice on this matter, but in any event, the preferred target optical density should be in the range of 1.5–1.9. A service engineer may be required to adjust the AEC to obtain this target value. If a facility is not able to provide viewboxes producing a luminance of at least 3000 cd/m², it will be difficult for the radiologist to view a film exposed to such a density. Under such conditions, it is likely that a lower target optical density will be chosen. If this is done, it must be realized that this will cause the film to be used in a suboptimal manner and there is a high probability that this will reduce diagnostic accuracy.

Important: Prior to conducting this test, the proper operation of the processor should be confirmed by checking the results of the sensitometry test.

7.3.6.1. Scope

- Objectives: To test the repeatability of the AEC; to establish the baseline value for the target optical density under AEC and ensure that this target is maintained with time; to evaluate the exposure time for a 45 mm thick PMMA test slab; to test the compensation of the AEC for different beam qualities and thicknesses; to test the performance of the density control;
- References [3, 4, 15, 20, 21];
- Minimum frequency: At acceptance/commissioning, annually, and after changes in the equipment, in the image receptor and/or in the processor.

7.3.6.2. Instrumentation

- (1) A set of PMMA slabs to give a total thickness of 2.0 cm, 4.5 cm and 7.0 cm.
- (2) Appropriate spacers (e.g. radiolucent U-shaped expanded polystyrene spacers of 8 mm thickness).
- (3) Exposure timer.¹⁶
- (4) Radio-opaque markers for film identification.
- (5) Loaded cassettes.
- (6) Densitometer.

7.3.6.3. Methodology

- (1) Select the clinically relevant automatic exposure mode. Select position '0' or 'normal' or 'N' on the density control selector.
- (2) Use the same cassette to carry out all the tests.
- (3) Place the loaded cassette in the Bucky.
- (4) Place a 45 mm thickness of PMMA aligned with the chest wall margin of the breast support and laterally centred and with the U-shaped polystyrene spacer on top to achieve a total thickness of 53 mm. Confirm that the PMMA covers the AEC sensors completely. It is important that the AEC sensor be placed in the same position every time the test is performed (e.g. at the position closest to the chest wall edge).
- (5) Place a radio-opaque marker to identify the film.

¹⁶ If the system provides a readout of the post-exposure mAs and the mA is known, the exposure time can be calculated from this information without the need of an external exposure timer.

- (6) If the mammography system is fully automatic (i.e. selects automatically the target filter combination and kVp), operate it in this manner. Otherwise, select the kVp and target filter combination used clinically for a breast whose X ray attenuation is equivalent to the 45 mm of PMMA¹⁷, for example 28 kVp; and make three exposures under the same conditions, keeping the same film in the cassette for AEC repeatability assessment. Note the mAs for each exposure.¹⁸
- (7) Replace the film and make a single exposure under the same conditions.
- (8) Process this film and all subsequent films according to the manufacturer's recommendations.
- (9) Measure the OD at a point placed at 40 mm from the chest wall, laterally centred, and record this value on the data collection sheet.
- (10) From the single exposure, determine the exposure time using a suitable exposure timer if needed.
- (11) Compare the measured optical density with the target optical density, OD_{target}, established at commissioning. If it differs by more than 0.2 OD, the density control should be adjusted to provide the correct target optical density. If the target optical density is 1.7 and the measurement is 1.4, for example, a density setting of +1 or +2 is likely to provide the proper OD. The remaining tests are then performed using this setting as the reference. When the next routine servicing is performed on the system, the density control should be recalibrated to provide this OD at the '0' setting.¹⁹
- (12) Select other thicknesses of PMMA (20 mm and 70 mm). Expose a single film at each thickness using the clinically relevant kVp and target filter combination for that thickness. If automatic mode is used clinically, then this mode should be utilized. A technique chart is essential in the absence of automatic mode.
- (13) Make an exposure at those positions of the density control that provide an OD within ± 0.3 of the '0' position using the beam quality settings selected in step (6) and 45 mm of PMMA. Note the mAs.
- (14) Process the films and measure the optical densities at a point placed at 40 mm from the chest wall and laterally centred; record these values on the data collection sheet.
- (15) Repeat the above tests for each receptor size.
- (16) Repeat the above tests in the configuration used for magnification mammography if applicable.

7.3.6.4. Calculation procedure

Repeatability:

- (1) Calculate the mean and standard deviation of the four mAs values.
- (2) Calculate the coefficient of variation (COV) of the four mAs values.

Exposure time for the 45 mm PMMA slab:

If an external exposure timer is not used, divide the post-exposure displayed mAs by the nominal mA and record the calculated time on the data collection sheet.

¹⁷ A 45 mm thick PMMA slab is equivalent to a 53 mm thick breast, and the use of plastic spacers is recommended to ensure the paddle is at 53 mm above the breast support. In some systems the paddle height will affect the automatic selection of exposure technique factors.

¹⁸ If it is noticed that the automatic system causes the kV or target filter combination to change for different exposures, change the thickness of PMMA slightly (by 5 mm) to force the system to select a single, consistent spectrum. Then repeat the test.

¹⁹ If it is found necessary to recalibrate the '0' or 'normal' setting of the AEC control, it is important to understand the underlying reason for this. A shift in sensitivity or X ray attenuation of one of the components of the imaging system can cause a need to employ either more or less radiation to achieve a particular OD. This may be the result of a conscious decision (e.g. replacing the screen film combination with one that is more sensitive or less noisy), in which case this can be acceptable, or it may signal a problem that must be corrected (e.g. poor processor replenishment) rather than masking the problem by readjusting the AEC control.

Constancy of OD with baseline value:

Compare the measured OD in the films obtained with a thickness of 45 mm PMMA and the density control in position '0', 'normal' or 'N' with the baseline OD_{target} established at commissioning (see Section 7.3.6.3, step (9)).

Compensation for different thicknesses and beam quality:

For each film, calculate the difference between the measured OD and the OD_{target}, and record the values on the data collection sheet.

Density control performance:

- (1) Calculate the differences in OD between films obtained for adjacent steps of the density control.
- (2) Note the mAs corresponding to each step of the density control.
- (3) Determine the mAs difference for consecutive steps of the density control, expressed as a percentage of the lower value.

7.3.6.5. Interpretation of results and conclusions

Tolerances:

- (1) Repeatability: Acceptable; COV in mAs $\leq 5\%$.
- (2) Exposure time for 45 mm slab:
Achievable: contact mode ≤ 1.5 s, magnification mode ≤ 2 s ;
Acceptable: contact mode ≤ 2 s, magnification mode ≤ 3 s.
Note: Meeting these tolerances should ensure that the impact of movement artefacts and film reciprocity law failure of the film on the image quality is minimized.
- (3) Constancy of OD with baseline value: Acceptable: $OD = OD_{target} \pm 0.20$.
- (4) Compensation for different thicknesses and beam quality:
Achievable: $OD = OD_{target} \pm 0.15$;
Acceptable: $OD = OD_{target} \pm 0.2$ for both contact and magnification mammography.
- (5) Density control: Change in OD for different steps: Acceptable: $0.1-0.2 OD$.²⁰

7.3.6.6. Recommendations and corrective actions

- (1) The AEC should be adjusted to obtain a target optical density between 1.5–1.9 OD when the density control is in position '0' ('normal' or 'N'). If this is not fulfilled, it should be determined which position of the density control is required to obtain this range of OD. The determined position should be recommended for clinical use until the maintenance person adjusts the AEC.
- (2) Because of the lower mA available with the small focal spot (magnification mode), the exposure times may be expected to be considerably longer than in contact mode. If times are found to be excessive and patient motion is a concern, an increase in kVp should be considered.
- (3) In case of any non-compliance, contact the person in charge of the maintenance service.

²⁰ An alternative approach, which will reduce the amount of film required for testing, is for the physicist to establish the relationship between a change in mAs (near the target operating point) and the change in OD for the film product used, and then to determine that the AEC adjustment provides for an adequately small change in OD to meet this requirement.

7.4. COLLIMATION SYSTEM

7.4.1. Collimation system

7.4.1.1. Scope

- Objective: To evaluate the coincidence between the radiation field and the image receptor; to evaluate the coincidence between the light field and the radiation field; to evaluate the alignment of the compression paddle with the margin of the breast support;
- References [3, 20, 24];
- Minimum frequency: At acceptance, annually, and after changes.

7.4.1.2. Instrumentation

- (1) Six radio-opaque objects (e.g. coins or clips): five of the same size and one larger.
- (2) Two loaded mammography cassettes: one that fits the Bucky and another, larger one, whenever possible.
- (3) Ruler.

7.4.1.3. Methodology

- (1) Place the larger cassette on the breast support, so that the cassette extends 20 mm beyond the margin of the chest wall.²¹
- (2) Place the other cassette in the Bucky.
- (3) Temporarily remove the compression paddle and turn out the room lights to make it easier to see the light field margins.
- (4) Turn on the collimation lights and place four radio-opaque objects on the margins of the light field, but just within the light field. The single large radio-opaque object should be on the margin of the light field on the chest wall.
- (5) Replace the compression paddle and locate it approximately as it would be used clinically (40–60 mm above the breast support).
- (6) Place the sixth radio-opaque object on the lower surface of the compression device, tangent to the inner surface of the lip of the compression device, adjacent to the chest wall. Place it sufficiently off the centre line such that it does not interfere with the sensor of the AEC.
- (7) Make an exposure with the AEC. In the event that the X ray unit does not permit irradiation without compression force, place an object of sufficient thickness under the support arm of the compression device and operate the compression paddle.
- (8) Process the films.
- (9) Determine the coincidence of the light field with the radiation field by measuring the distance between the margins of the X ray field and those of the radio-opaque objects. These measurements should be done on the film that was in the cassette placed on the breast support. Note the result on the data collection sheet.
- (10) Overlap the films of the two cassettes on the viewbox so that the most external margins of the radio-opaque objects are coincident. For each of the four edges, determine the deviation between the radiation field (as determined from the film exposed on the breast support) and the image field (as determined from the film exposed in the Bucky). Note the result on the data collection sheet.
- (11) Determine the coincidence between the edge of the compression paddle and the edge of the image receptor corresponding to the chest wall by measuring the extent to which the radio-opaque object placed on the compression paddle is only partly visible in the image of the film located in the Bucky. Note the result on the data collection sheet.

²¹ If necessary, to avoid overexposing the film, place a developed film between the intensifying screen and the film, or place the film with the emulsion in reverse. If large sized cassettes are not available, two smaller ones placed side by side may be used to ensure complete coverage of the irradiated field.

- (12) Repeat the procedure for all cassette sizes and collimators that are commonly used.
- (13) Repeat the procedure for all X ray tube targets and for magnification mode.

7.4.1.4. *Interpretation of results and conclusions*

Tolerances:

- (1) Coincidence between light field and radiation field: achievable; $\leq 1\%$ of focus-film distance (FFD) on any side.
- (2) Coincidence between film and radiation field: achievable; the beam completely irradiates the image receptor, but does not extend beyond the breast support except at the chest wall where it can extend beyond the breast support to a maximum of 5 mm acceptable; may be inside the breast support by $\leq 2\%$ of the FFD for the three sides other than at the chest wall. For the chest wall, it can extend beyond the breast support to a maximum of 5 mm.²²
- (3) Compression paddle alignment: acceptable; the chest wall edge of the compression paddle is not visible in the image and should not extend beyond the image receptor edge by more than 1% of the FFD.

7.4.1.5. *Recommendations and corrective actions*

If measurements are outside of tolerance, contact the person in charge of the maintenance service.

²² It is preferred that the entire film be exposed to radiation as any unexposed area of film can result in glare when viewing the processed film.

7.5. IMAGE VIEWING CONDITIONS

7.5.1. Luminance and homogeneity of the viewboxes

7.5.1.1. Scope

- Objective: To verify that the viewboxes provide adequate viewing conditions for mammography;
- References [3, 4, 15];
- Minimum frequency: At acceptance, annually, and after changes.

7.5.1.2. Instrumentation

Light meter (photometer) to measure luminance (range: 100–7000 cd/m²).

7.5.1.3. Methodology

- (1) Assess the need for replacing defective fluorescent tubes, noting lack of cleanliness, colour of the tubes and the viewbox screens, vibrations, etc.²³
- (2) Record the information on the data collection sheet.
- (3) Select five measurement points (one point should be centrally located and the other four should be located towards the corners of the viewbox), at least 50 mm away from the edges.
- (4) Place the light meter in contact with the surface of the viewbox at each selected point.
- (5) Measure the luminance at each point. Record the values on the data collection sheet.
- (6) Repeat for each viewbox in use.

7.5.1.4. Calculation procedure

In a single viewbox or single panel of a viewbox bank, select the central luminance value (L_c) and the most discrepant reading recorded (L_{disc}). Apply the following equation:

$$\text{Maximum deviation (\%)} = 100 \left| \frac{L_{disc} - L_c}{L_c} \right| \quad (4)$$

- (1) For viewboxes consisting of several panels or for viewboxes adjacent to one another: the maximum deviation between the central luminance value of any panel, L_{cx} , and the calculated mean of the central luminance of all panels, L_{mean} , is found using the equation:

$$\text{Maximum deviation (\%)} = 100 \left| \frac{L_{cx} - L_{mean}}{L_{mean}} \right| \quad (5)$$

7.5.1.5. Interpretation of results and conclusions

Tolerances:

- (1) Luminance:
Acceptable: Maximum luminance for each panel: >3000 cd/m² (nit);

²³ It is recommended that an inventory of the viewboxes in the institution be maintained, noting their location and age.

- (2) Luminance uniformity:
Acceptable: <30% maximum deviation for different areas of a single viewbox or single panel in a viewbox bank and <15% maximum deviation between central luminance of panels in a viewbox bank or between adjacent viewboxes.

7.5.1.6. Recommendations and corrective actions

- (1) In case of detecting abnormalities during the visual inspection, contact the person in charge of the maintenance service.
- (2) If the luminance level is below the required level or a lack of uniformity exists that surpasses the tolerance, replace the fluorescent tubes in the poorly performing viewboxes as soon as possible.
- (3) Triphosphor or quadphosphor tubes are recommended because of their increased luminance output. If possible, all tubes should be replaced in a viewbox or viewbox bank at the same time to ensure uniformity. It is advisable to purchase fluorescent tubes in one batch to ensure that the colour matches between the tubes.

7.5.2. Interpretation room ambient illumination

7.5.2.1. Scope

- Objective: To verify that the interpretation room has adequate viewing conditions for mammography;
- References [3, 4];
- Minimum frequency: At acceptance, annually, and after changes.

7.5.2.2. Instrumentation

Light meter (photometer) to measure illuminance (range: 1–1000 lx).

7.5.2.3. Methodology

- (1) Place the photometer at a distance of 0.5 m from the previously turned off viewbox²⁴, directed towards the environmental lights, and measure the background light levels.
- (2) Record the value on the data collection sheet.

7.5.2.4. Interpretation of results and conclusions

Tolerances:

Room illuminance:

Achievable: ≤ 10 lx;

Acceptable: ≤ 50 lx.

7.5.2.5. Recommendations and corrective actions

If the area illumination conditions are greater than the recommended values, contact the person in charge of the maintenance service to modify the illumination of the area.

²⁴ If there is more than one viewbox, check that all of them are turned off.

7.6. IMAGE QUALITY

7.6.1. Evaluation of image quality

7.6.1.1. Scope

- Objective: To ensure that the optical density, image contrast and image quality are acceptable;
- Reference [3];
- Minimum frequency: At acceptance, annually, and after changes.

7.6.1.2. Instrumentation

- (1) Breast phantom. In this publication, the ACR mammography accreditation phantom has been taken as a convenient example.²⁵ It contains the following objects:
 - 6 fibres with diameters of 1.56 mm, 1.12 mm, 0.89 mm, 0.75 mm, 0.54 mm and 0.40 mm;
 - 5 groups of simulated microcalcifications with diameters of 0.54 mm, 0.40 mm, 0.32 mm, 0.24 mm and 0.16 mm;
 - 5 masses with thicknesses of 2.0 mm, 1.0 mm, 0.75 mm, 0.5 mm and 0.25 mm (Fig. 46).
- (2) Disc, composed of PMMA, 4 mm in thickness and approximately 10 mm in diameter (included with the ACR phantom).
- (3) Densitometer (range: 0.00–4.00 OD).
- (4) Cassette (it is recommended that the same cassette always be used) and film.
- (5) Magnifier: 2× or greater. It is recommended that the magnifier contain a second area that provides 4× or 5× magnification.
- (6) Viewbox.

7.6.1.3. Methodology

- (1) Load the cassette with a film from the box that is in use for patients. Wait the time recommended in order to have good screen film contact. (If it is the first film taken from a new box, note the emulsion number in the control sheet control.)




FIG. 46. Schematic representation of the ACR phantom — the chest wall edge is at the bottom.

²⁵ Many facilities have a phantom available that is specified by the mammography accreditation phantom of the American College of Radiology. This phantom is equivalent to a thickness of approximately 40 mm of PMMA and simulates a breast that is 42 mm thick and composed of a central region of 50% fibroglandular tissue and 50% fat tissue, surrounded above and below by a layer of fat that is 4 mm thick.

- (2) Place the loaded cassette in the Bucky.
- (3) Place the phantom on the breast support aligned with the chest wall and centred laterally.
- (4) Place the PMMA disc on the phantom in such a way that it does not interfere with any phantom detail or the AEC sensor. It is recommended that it be located slightly below and between the first and the second fibre.
- (5) Lower the compression paddle until it touches the acrylic disc.
- (6) Confirm that the AEC sensor is under the phantom, and the PMMA disc does not interfere with the AEC sensor.
- (7) Select the technique factors that are used in clinical practice: target, filter, kVp, grid, density control position, operation mode (semi-automatic or automatic).
- (8) Make an exposure.
- (9) Record on the corresponding data collection sheet the exposure factors and technique used.
- (10) Process the film in the processor normally used for mammography, always introducing it with the same orientation.
- (11) Determine the optical density at the following points, noting them on the data collection sheet:
 - At the geometric centre of the accreditation phantom image;
 - At the centre of the image of the PMMA disc;
 - Adjacent to the acrylic disc (perpendicular to the direction of the anode–cathode).
- (12) Calculate the difference of densities between the optical density of the PMMA disc and the density of the adjacent point. This is the ‘contrast index’. Note the result on the data collection sheet.
- (13) View the film with the magnifier on the viewbox used clinically, placing a mask in the area of the viewbox that is not being used, and with a very low level of environmental illumination in the reading room.
- (14) Evaluate the image according to the evaluation method outlined in the following. Note the results on the data collection sheet.
- (15) Observe the image and verify if there are artefacts that can be confused with any of the phantom details. With the magnifying lens, examine carefully the image for non-uniform areas, dirt or dust, lines (if the grid is used), processing artefacts or any other type of artefact. Mark them with a circle. Compare the image with previous images.
- (16) Investigate the causes of any artefacts.

7.6.1.4. *Interpretation of results and conclusions*

Fibre score:

- Totally visualized: 1;
- Partially visualized (more than a half-fibre): 0.5;
- Partially visualized (less than a half-fibre): 0.

Each fibre is scored individually, starting with those of greatest thickness and stopping with the fibre that has been assigned a score of 0.5. All the assigned scores should be added. The image is analysed, looking for those artefacts that can be confused with fibres but whose positions are misplaced or in inappropriate orientations. If these fibre-like artefacts are detected, their scores (1 or 0.5) have to be subtracted from the total score.

Microcalcifications group score (it is helpful to view these under 4× or 5× magnification):

- Groups in which 4 or more microcalcifications are visualized: 1;
- Groups in which 2–3 microcalcifications are visualized: 0.5;
- Groups in which less than 2 microcalcifications are visualized: 0.

Each group is scored independently, starting with the microcalcifications group of greatest size. Beginning with the group of largest calcifications, add the scores for the groups, stopping at the first group encountered for which the score is 0 or 0.5 (even if further groups of smaller calcifications are visible). Any ‘apparent’ microcalcifications visualized outside the position indicated by the phantom manufacturer are artefacts. When such

artefacts are observed, their scores should be subtracted from the score given to the last visualized group. For example, if the last evaluated group has received a score of 1 because 4 microcalcifications were visualized, and there is an artefact that can be confused with a microcalcification, then 1 microcalcification is subtracted from the count in the last group. Accordingly, instead of 4 microcalcifications being counted, the group now contains 3 and the group score decreases to 0.5.

Masses score:

- Totally visualized: 1;
- Partially visualized (the form is not generally circular): 0.5.

Each mass is scored independently, starting with the one of greater diameter. The evaluation stops when a score of 0.5 or smaller has been assigned to a mass. Add all the mass scores and note the result on the data collection sheet. Any masses visualized outside the position indicated by the phantom manufacturer are artefacts. When such artefacts are detected, the score given to the last evaluated mass (1 or 0.5) should be subtracted from the total score. For example, if the last evaluated mass has a score of 0.5 and an artefact is detected, 0.5 should be subtracted from the total score.

Tolerances:

- (1) Background optical density: acceptable: $OD = OD_{\text{target}} \pm 0.20$. Note that the OD_{target} should be between 1.5 and 1.9 OD.
- (2) Image contrast index:
Achievable: ≥ 0.55 ;
Acceptable: ≥ 0.40 for the density difference between PMMA disc and adjacent point.
- (3) Image quality:
Acceptable: Fibres: ≥ 4 ;
Microcalcifications: ≥ 3 ;
Masses: ≥ 3 .

From one site visit to the next one, the total score assigned to any of the three objects should not vary by more than 0.5 points. If this value is surpassed, the causes should be investigated.

7.6.1.5. Recommendations and corrective actions

- (1) If the image quality deteriorates over time, it will be necessary to carry out other tests (e.g. kVp, AEC and processor) to determine the source of the change.²⁶
- (2) If there are several different types of film in use, this test should be carried out for each type.

²⁶ Due to the subjective component associated with the observer, it is recommended that the test always be performed by the same person, using the same criteria, viewbox, magnifier lens and illumination conditions.

7.6.2. System spatial resolution

7.6.2.1. Scope

- Objective: To determine the system high contrast resolution;
- References [3, 4, 15];
- Minimum frequency: At acceptance, annually, and after changes.

7.6.2.2. Instrumentation

- (1) Resolution pattern of up to 20 lp/mm.
- (2) PMMA slabs.
- (3) Loaded cassette.
- (4) Magnifier: 2× or greater. It is recommended that the magnifier contain a second area that provides 4× or 5× magnification.

7.6.2.3. Methodology

- (1) Load the cassette with a film from the box in clinical use. Wait the recommended time to ensure that there is good screen film contact. If this is the first image taken from a new box of film, note the emulsion number on the data collection sheet.
- (2) Place the loaded cassette in the Bucky.
- (3) Place the resolution pattern centrally on 40 mm of PMMA, 1 cm from the chest wall.
- (4) Confirm that the AEC sensor is not under the resolution pattern. Make two exposures, on different films, one with the pattern perpendicular to the anode–cathode direction, and the other one with the pattern parallel to the anode–cathode direction. If the equipment does not have AEC, expose the film with the technical factors (kVp, grid, filter, target) clinically used for a compressed breast of 45 mm (equivalent thickness to 40 mm of PMMA).
- (5) Place the films on the viewbox and, using the magnifier lens, note the number of line groups that can be observed clearly, starting with the most easily resolved.
- (6) Note the result on the data collection sheet.

7.6.2.4. Interpretation of results and conclusions

Tolerance:

Achievable: ≥ 15 lp/mm in both directions;

Acceptable: ≥ 11 lp/mm in both directions.

7.6.2.5. Recommendations and corrective actions

If important variations in the resolution are observed, the size of the focal point should be checked.

7.7. DOSIMETRY

7.7.1. Incident air kerma at the entrance surface of the phantom

7.7.1.1. Scope

- Objective: To estimate the incident kerma in air (without backscatter) at the position corresponding to the entrance surface of the phantom;²⁷
- References [3, 4, 25–27];
- Minimum frequency: At acceptance, annually, and after changes.

7.7.1.2. Instrumentation

- (1) Appropriate dosimeter for mammography;
- (2) Phantom simulating the standard breast (45 ± 0.5 mm of PMMA, see Appendix VIII);
- (3) Radiolucent spacer of thickness 8 mm (e.g. U-shaped expanded polystyrene);
- (4) Tape measure;
- (5) Thermometer;
- (6) Barometer.

7.7.1.3. Methodology

The incident air kerma without backscatter is measured using the same exposure factors selected to expose a phantom equivalent to the standard breast (45 mm of PMMA). The method to determine the incident kerma in air at the entrance surface of the phantom consists of two parts:

- (1) Determination of the mAs value for the correct exposure of the phantom in clinical conditions;
- (2) Measurement of the incident kerma in air for the given mAs.

This procedure is proposed since, generally, the exposure is made using the AEC and the mAs may not be known a priori. The procedure is outlined in the following discussion.

- (1) Determination of the mAs value:
 - (i) Select the typical imaging conditions (kVp, anode–filter) for imaging the standard breast, with the grid and with a loaded cassette in the Bucky. For a system with a fully automatic AEC, these factors can be determined by making an exposure of the phantom in this mode.
 - (ii) Place the phantom on the breast support aligned with the corresponding margin of the chest wall and laterally centred. If the AEC is sensitive to the compressed breast thickness, add the 8 mm thick spacer at the outer edges of the phantom. Ensure that the AEC sensor is under the phantom and located at the position closest to the chest wall edge.
 - (iii) Apply a standard compression force (e.g. 100 N) to the phantom so that the total ‘breast’ thickness is 53 mm. Do not use magnification.²⁸
 - (iv) Measure the distance from the surface of the PMMA to the focus of the X ray tube.

EITHER

²⁷ Without the standard breast phantom, because the measurement should be done without scattered radiation.

²⁸ Note that in some systems and in some modes of operation, the compressed breast thickness is utilized to determine the technique factors and this thickness is, in turn, dependent on the degree of compression applied.

For X ray units with AEC and mAs indication at the end of the exposure:

- (v) Make an exposure with the AEC. The density selector should be set as for clinical images of the breast that is being simulated.
- (vi) Record the values of the tube voltage (kVp), anode–filter combination and mAs_{auto} on the data collection sheet.

OR

For X ray units with AEC and without mAs indication at the end of the exposure:

- (v) Place the dosimeter to one side of the phantom at a height so that the centre of its active volume is at the entrance surface of the phantom (i.e. 45 mm above the breast support table).²⁹ Ensure that the AEC sensor is not impeded in any way and the compression paddle is located within the beam without compression.
- (vi) Make an exposure with the AEC, using the same clinical conditions used for a breast similar to the one that is being simulated.
- (vii) Record the reading of the dosimeter (M_{auto}) on the data collection sheet.
- (viii) Select the manual mode with the same set-up and make an exposure at the same kVp and anode filter combination using a clinically relevant mAs value (or exposure time).
- (ix) Record the reading of the dosimeter (M_{manual}) together with the mAs value (mAs_{manual}) on the data collection sheet.
- (x) Deduce the mAs value (mAs_{auto}) for the phantom exposure with the AEC, using the mAs_{manual} value, through the following formula:

$$mAs_{\text{auto}} = (M_{\text{auto}}/M_{\text{manual}}) \times mAs_{\text{manual}} \quad (6)$$

- (xi) Record the kVp and the target–filter combination used on the data collection sheet.
- (2) Measurement of the incident kerma (with the value of mAs obtained):
- (i) Remove the phantom and place the centre of the ionization chamber 45 mm above the breast support, laterally centred and 40 mm from the margin of the chest wall.
 - (ii) Make an exposure in manual mode with the previously determined mAs value, mAs_{auto} . If it is not possible to select the exact mAs value, a numerical interpolation of the dosimeter readings, M_1 and M_2 , can be done for the exposures with mAs values immediately higher and lower than the desired mAs_{auto} , applying the following formula:

$$M_{\text{auto}} = 0.5 \times mAs_{\text{auto}} \left(\frac{M_1}{mAs_1} + \frac{M_2}{mAs_2} \right) \quad (7)$$

- (iii) Obtain the value of the incident kerma K_i from M_{auto} using the following formulas and record the value on the data collection sheet:

$$K_i = M_{\text{auto}} \times N_{\text{mammo}} \times k_{TP} \quad (8)$$

$$k_{TP} = \frac{(273.2 + T)}{(273.2 + T_0)} \times \frac{P_0}{P} \quad (9)$$

²⁹ The dosimeter is placed at 45 mm rather than 53 mm because the inverse square law correction to 53 mm is already incorporated in the conversion factors used in Table 4.

where

k_{TP} is the correction factor for temperature and pressure;

N_{mammo} is the value of the calibration factor for beam quality;

T_0 and P_0 are the values of pressure and temperature at which the dosimeter is calibrated (if applicable).

7.7.2. Determination of the mean glandular dose (D_G)

7.7.2.1. Scope

- Objective: To estimate the mean glandular dose (D_G) for a standard breast;
- References [3, 4, 25–27];
- Minimum frequency: At acceptance, annually, and after changes.

7.7.2.2. Methodology

- (1) The mean glandular dose (D_G) is obtained from the incident kerma in air and relevant conversion coefficients using the following formula:³⁰

$$D_G = g_{53} c_{53} s K_i \quad (10)$$

where

K_i is the entrance air kerma at the surface of the 45 mm thickness of PMMA, measured without backscatter (see Section 7.7.1);

g_{53} is the factor that converts the entrance air kerma to the mean glandular dose for the 53 mm thick standard breast;

c_{53} is the conversion factor which allows for the glandularity of the 53 mm thick standard breast;

s is the factor which gives a correction that depends on the target filter combination.

- (2) Use the g and c factors, which are dependent on the HVL of the spectra utilized, as provided in Table 5. The HVL value is obtained following the method described in Section 7.3.3.
- (3) Apply the values of the s factor for the relevant target filter combination provided in Table 6.
- (4) Record the values of g , c , s , HVL and D_G on the data collection sheet.

TABLE 5. CONVERSION FACTORS g AND c FOR CALCULATING D_G FOR THE STANDARD BREAST (53 mm THICK) FROM MEASUREMENTS WITH A 45 mm THICK PMMA PHANTOM

HVL (mm Al)	g_{53} (mGy/mGy)	c_{53}	Product of g_{53} and c_{53}
0.30	0.155	1.109	0.172
0.35	0.177	1.105	0.196
0.40	0.198	1.102	0.218
0.45	0.220	1.099	0.242
0.50	0.245	1.096	0.269
0.55	0.272	1.091	0.297
0.60	0.295	1.088	0.321

³⁰ This definition uses the United Kingdom and European formula using an expected average glandular content instead of that given previously [28] where the glandularity was assumed to be 50%.

TABLE 6. s FACTORS FOR ANODE–FILTER COMBINATIONS

Target filter combination	s factor
Mo/Mo	1.000
Mo/Rh	1.017
Rh/Rh	1.061
Rh/Al	1.044
W/Rh	1.042

Tolerance:

Achievable: $D_G \leq 2.0$ mGy;

Acceptable: $D_G \leq 2.5$ mGy.

7.7.2.3. Recommendations and corrective actions

- (1) Observe the variation of D_G with time and, if the tolerances are exceeded significantly, investigate the possible causes and take the necessary corrective measures.
- (2) It is recommended that the mean glandular dose also be estimated from actual patient exposures periodically, and values compared to diagnostic reference levels that have been established at the local or national level [25–29].

Appendix I

GUIDELINES FOR THE TRAINING OF A MEDICAL PHYSICIST SPECIALIZING IN MAMMOGRAPHY

The clinically qualified physicist in medical imaging (radiology) specializing in mammography should have at least:

- (1) A university degree in physics, engineering or an equivalent physical science.
- (2) A master's degree in medical physics (or equivalent) in an academic postgraduate programme of at least one year's duration. This requires studies in several areas of medicine (e.g. radiodiagnostics, nuclear medicine and radiotherapy).
- (3) At least two years of full time equivalent comprehensive clinical in-service training in the physics of medical imaging. This physics residence training will be under the supervision of an experienced medical imaging physicist. Where the academic studies include a considerable clinical training component, this should be taken into account in the fulfilment of the time requirement. Preferably, this training should be approved by a suitable professional body, that is, a Board that will issue a professional certification. It is emphasized that the holder of a university degree in medical physics without the required clinical training cannot be considered clinically qualified.
- (4) The responsibilities of the medical imaging physicist cover several major areas: dosimetry, radiation safety, quality control, equipment selection and possibly training.

Appendix II

MAMMOGRAPHY ROOM DESIGN

The requirements concerning the design of a mammography room outlined in the following discussion,³¹ and also in Appendix III, have been found to be useful in practice. The requirements are advisory only, as some of the suggested specifications may be difficult to achieve in Member States. Nevertheless, it is strongly recommended that most of the design recommendations be implemented.

II.1. LAYOUT

The layout of the mammography room should allow the radiographer easy access to the patient from all sides during the procedure. This generally means that there must be at least a distance of 2 m from the centre of rotation of the gantry to the nearest wall. This will allow adequate clearance and access during a mediolateral oblique view. There must be enough room at the operating console for an operator and one observer (e.g. trainee) to stand behind the protective shield.

II.2. CASSETTE STORAGE

There must be a place to store the loaded cassettes in the room, behind the protective shield (cassettes should be stored vertically, not lying flat). It is a good idea to have a number of shelves or slots with clear markings indicating where exposed and unexposed cassettes are to be placed. At least eight cassettes of each size must be accommodated.

II.3. OTHER STORAGE AND WORKING AREAS

There must also be adequate, easily accessible storage for the removable grid assemblies, compression paddles and collimator plates. Convenient storage for the image quality phantom and PMMA blocks should be close to the unit. There should be a work surface at least large enough to accommodate the patient record sheets, and a flat surface for writing.

II.4. VIEWBOX

A viewbox may be mounted in the mammography room for the review of images. It should be placed at a comfortable height for the radiographer, normally with the top edge of the viewing area at approximately 155 cm above the floor. There should be a work ledge below, not to exceed 30 cm in depth, permitting the radiographer to approach the viewbox. The intensity of this viewbox should closely match that of the viewbox used by the radiologist for reading mammograms, and a 'mammoviewer'/magnifier and adequate masking should be available.

II.5. ROOM LIGHTING

In order to maintain a room brightness of less than 50 lx for correct viewing of images for initial quality assurance by the technologist and to allow visualization of the field illuminator, the room lights should be equipped with a dimmer switch. The windows (preferably frosted for privacy) should be fitted with shades and/or drapes.

³¹ This appendix has been adapted from material developed by the Ontario Breast Screening Program, Toronto, Canada. The help of that organization is gratefully acknowledged.

II.6. SINK AND CLEANING PRODUCTS

A sink in the room is not a requirement, however, for infection control it is a good idea for the radiographer to wash his or her hands between patient examinations, especially if a given patient is experiencing discharge from the nipple. Cleaning supplies for the compression paddles and tabletop should be easily accessible, and used frequently (after each patient). Check with the manufacturer of the mammography unit for recommended cleaning products and instructions specific to each unit.

II.7. ELECTRICAL SERVICES REQUIRED

Most current mammography units operate from standard single phase AC circuits with a power requirement of around 10 kVA instantaneous (for up to 6 s), and a standby requirement of 1.5 kVA. Most units are configurable to regional power, and operate at either 50 Hz or 60 Hz. A mains isolator switch or circuit breaker fused at 25 A for 200–240 V is normally required in the room.

II.8. ENVIRONMENTAL SERVICES REQUIRED

Since most mammography rooms are very small, and the door remains closed during the procedure, there are often problems with environmental control. The mammography unit dissipates a considerable amount of heat with the maximum specified being 3 kW. This heating should be considered in addition to that from the two people who are occupying the room for the entire day.

There is a requirement for fresh make-up air as well as a cold air return within the room. The fresh air diffuser should not be directly in the centre of the room (over the patient), as it is sometimes very uncomfortable to have air blown on bare shoulders and back. The room should have independent heating/air conditioning control. There should be adequate humidity regulation so that there is reduced static buildup in the room. This is both for patient comfort (nipple to table discharge) and to prevent static attraction of dust to the compression paddles and static artefacts on the films.

II.9. GENERAL INSTALLATION

As dedicated mammography systems are often of freestanding design, it should be verified that the entire system is installed properly in the examination room and is stable. Floor loading for the gantry can be up to 400 kg, concentrated on the area of the baseplate of approximately 0.3 m². The control unit (up to 150 kg) should be level and there should be convenient access for service and maintenance. There should be no obstruction that hinders the range of mechanical motions provided by the system. It should be ensured that the mechanical stops and magnetic locks are properly set and adjusted to prevent accidental motion or collision against other mechanical parts or the patient.

II.10. RADIATION SHIELDING AND REGULATORY APPROVALS

Before construction of the facility is started, approval for the radiation shielding should be sought, as with any X ray facility. The medical physicist can assist in the calculation of barrier thicknesses and provide advice on documentation requirements. Provision of adequate shielding for mammography is not normally difficult. Shielding is required only to reduce exposure from scattered radiation, since the primary beam is limited to the area of the image receptor support. Typically, a total thickness of 24 mm of gypsum wallboard (two standard sheets) will provide adequate shielding. The radiographer must have a transparent shield (nominally equivalent to 0.3 mm Pb) behind which he or she can stand while making an exposure. The exposure control must be designed so that it cannot be operated from outside this shielded area. The door to the room should be solid core; normally a fire rated door will provide enough shielding.

II.11. RADIATION PROTECTION ACCESSORIES

In modern screening mammography, there is negligible exposure of radiation sensitive sites other than the breast. The main value of radiation protection apparel is psychological. This may be a factor in dealing with patients who are radiationphobic or who have had thyroid disease. If such apparel is to be supplied, it should only be done at the request of the patient, and the apparel should not be kept on display in the examination room. The presence of the aprons and collars in the mammography room might suggest that their use is accepted practice, which is not the case.

Both measurements and calculations show that the amount of radiation reaching the thyroid during mammography is negligible. The amount of radiation reaching the ovaries is even less, due to the attenuation of the breast support and of the overlying tissue. Virtually all primary radiation is stopped by the breast and the breast support, only an extremely low level of scatter reaches other parts of the body. The calculated dose to the thyroid for one four-view examination is less than 0.03 mGy. This is about 1% of the dose that would be received by the breast during the examination, and is equal to the dose that would be received by the thyroid from three days of natural background radiation. In other words, this would be the equivalent of the patient receiving 368 days of natural background radiation per year instead of 365 days of background radiation that would be received without the examination. Natural variations in background radiation from locality to locality are much greater than this.

Appendix III

DARKROOM DESIGN

III.1. GENERAL REQUIREMENTS

The darkroom should be designed to provide the following conditions that are required for mammography:³²

- (1) Appropriate lighting conditions for the handling of photosensitive film.
- (2) Suitable storage for unexposed film in use with particular regard to radiation levels.
- (3) Space for the film processor and accessories.
- (4) Space for the storage of chemical solutions.
- (5) Adequate facilities for the preparation of chemical solutions and for the cleaning of processor rollers (large sink and hose).
- (6) Appropriate ventilation.
- (7) Provision of a silver recovery system.
- (8) Adequate workspace for the loading of films, and a viewing area and workspace for sorting the films.
- (9) Fire alarms and other warning signals should be audible inside the darkroom.

Note that these requirements still apply for daylight, or automatic feeding processors, in order to reduce image artefacts. Ventilation of the processor and handling of chemicals remains the same.

III.2. DARKROOM LOCATION AND RADIATION CONDITIONS

The darkroom should be located adjacent to the mammography room, and there should be a clear path from the mammography room to the entrance of the darkroom. A pass-box directly from the mammography room to the darkroom will permit increased throughput if there is a darkroom technician, and allows for the storage of loaded cassettes.

Care should be taken that there are no excessive radiation levels in the darkroom and film storage areas. This can be checked using TLD dosimeters. The acceptable level is <20 $\mu\text{Gy}/\text{week}$.

III.3. WORK SURFACE

The darkroom area must have a work surface at least 1.3 m long for the loading and unloading of cassettes. This surface should be constructed of a hard, anti-static material that is easily cleaned. A light colour is appropriate for all surfaces in the room, as this reduces the safelighting wattage required and allows for the detection of light leaks. Beige or mid-grey is, however, a poor choice for darkroom counters or floors since film emulsion is this colour, and a film may not be seen if it is placed on the counter. The film can be stored in film bins under the counter; however, there are a number of convenient photographic paper safes available for countertop use, which are easier to use than the standard lead-lined radiographic models. Space should be reserved on the counter for a sensitometer and film identification camera if required. The film-feed tray of the processor should be oriented to allow for good workflow.

III.4. DARKROOM DOOR

The entrance door to the darkroom must be light-tight. This involves careful weather-stripping of the door, including the installation of a threshold sealer (the automatic style is preferred, as there will be less wear and

³² This appendix has been adapted from material developed by the Ontario Breast Screening Program, Toronto, Canada. The help of that organization is gratefully acknowledged.

tear). The seal around the processor is critical and usually requires weather-stripping and perhaps black sealant along the floor. The door should have a lock in the knob that will keep people from unintentionally entering the room while film is being unloaded, however, will permit access in case of emergency. A bathroom or entrance lockset is usually used. In order to evaluate light leak levels, one should allow dark adaptation for at least 5 min. Darkroom fogging tests should be carried out immediately if any fogging is suspected. A darkroom ventilation louvre installed in the door will serve to allow air flow in and out of the darkroom, and to reduce the turbulence generated when closing the door.

III.5. LIGHTING

Two levels of lighting should be available in the darkroom. A strong white light should be provided by surface mounted incandescent fixtures. Fluorescent bulbs tend to emit a long afterglow that may fog films if they are switched off just before the films are unloaded or loaded. Recessed lights (pot lights, downlights) are discouraged because they allow light leaks from above the ceiling area into the room. If a pot light or standard recessed fluorescent fixture is used, there must be a fireproof box built above the ceiling to contain the fixture, with all penetrations sealed.

Safe lighting should be provided, with manufacturer recommended filters in ceiling mounted fixtures. The maximum recommended bulbs are 15 W, for direct safe lighting, but 25 W bulbs may be used if the safe lighting is indirect. *It is important to check current recommended safe light wattages with the film provider*, as some films require much lower power ratings than these levels (e.g. Kodak recommends a 7.5 W bulb with their Min-R 2000 and EV films). There should be a switched duplex receptacle mounted in the ceiling for the safe light. The switches for the safe light and white light should be located at the door, but at different heights above the floor to avoid errors in the dark.

III.6. SAFE LIGHT FILTERS

Photographic safe light filters are of the absorption gelatin type and are designed so that darkrooms have as high a degree of illumination as is consistent with the safe handling of photosensitive materials. Safe lights are designed to transmit in the region of the visible spectrum to which the photographic emulsion is least sensitive. If safe lights are chosen carefully, the total safe illumination in the darkroom can often be sufficiently high for very comfortable vision.

Manufacturers recommend replacement of the safe light filters every two years. At the time of replacement, document the date on the filter with a marking pen. The filter should be inserted in such a way that it will be possible to read the writing and easily determine when the next replacement is due. To ensure proper operation, the safe light test (see Section 6.3.3) should be performed once the filters have been replaced.

III.7. DARKROOM VENTILATION

It is essential that there be adequate filtered and humidified air introduced into the room, and that sufficient exhaust be provided out of the darkroom in order to remove fumes. There should be a positive pressure in the darkroom, drawing air in through a vent in the door. The return duct should have at least four bends between the grill in the darkroom and a grill in another room in order to prevent light fogging of the film due to reflected light. An in-line dust filter mounted at the supply louvres will reduce dust as well as light transmission. Any dust introduced through a supply duct will seriously degrade the image quality and may cause damage to the screens.

As noted previously, a light-tight vent in the darkroom door will provide constant flow-through and turnover of new air at floor level, and reduce pressure variations as the door is opened and closed. Pressure variations stir up dust particles.

The fumes from processor chemicals are often quite obvious at levels well below the listed harmful levels. While this is true, the chemical smell is sometimes unpleasant, and there are developed 'sensitivities' which can result from low level repeated contact with these fumes. This has resulted in an increased level of concern by radiographers about exposure to the fumes. From a practical standpoint, with a well functioning darkroom, there should be no odour.

Proper ventilation will also tend to reduce dust artefacts and improve processor operation. Without adequate 24 hour direct processor exhaust, evaporated solutions from the processor will tend to build up on the inside of the processor and on rollers. The deposits may cause dust artefacts on the films, and the vapours will cause corrosion on the other components of the processor. Adequate exhaust is also required to remove the heat and humidity originating from the dryer of the processor.

Darkroom general supply and return air exchange should be enough to remove any residual fumes from the replenisher tanks, which should be well covered to reduce evaporation.

Darkroom ventilation conditions:

- (1) The processor is to be exhausted to the outside of the building through a special fume exhaust. The outlet separation from the building supply air must meet building code regulations (typically not within 2 m of a window or air intake). This duct should be capable of removing approximately 2.5 m³/min through a 10 cm diameter duct. There should be a 20–50 mm gap in the duct just after it enters the room to ensure that there is less than a 1 mm vacuum at the processor outlet. (See, for example, Kodak service bulletins Nos 101 and 158 [30].)
- (2) The exhaust fan for the processor must be left on at all times, and operate throughout the day and night, even when the processor is turned off.
- (3) It is recommended that facilities have an ‘air proofing switch’ installed in the processor exhaust vent with an indicator lamp visible in the radiographer work area which will light up when the exhaust is working.
- (4) There should be at least 10 air changes per hour in the room to ensure removal of chemical fumes from the area. The supply should be located so that it does not short circuit and feed directly to the exhaust. A darkroom supply of 50% of the exhaust rate is ideal.
- (5) Supply air should be between 15°C and 21°C with a humidity level between 30% and 70%.
- (6) Supply air should be filtered to remove dust and airborne sources of artefact. This may be done by the use of a local filter. Note that visual inspection of fibreglass filters should be undertaken at least every three months and they should be changed at regular intervals, or as required.
- (7) After renovations or if ducts are found to be excessively dusty, duct cleaning must be carried out back to the nearest filter.
- (8) Evaluation of exhaust fan operation. ‘Flutter test’: with the processor turned off, a piece of tissue is to be held near the gap in the exhaust duct. The tissue should be seen to be drawn towards the gap. If there is not noticeable movement in the tissue, the flexible duct is to be removed and the negative pressure measured 25 cm into the opening of the duct, using a manometer. If there is inadequate suction, or if the tissue blows away from the gap, the building manager should be requested to obtain the services of a qualified service person.

III.8. AIR FILTRATION

For local air filtration, HEPA filters can be obtained relatively cheaply combined with an ionizer device. This device charges the air molecules, which henceforth attract dust particles, drawing the dust out of the air. Such filters are often also combined with humidifiers.

III.9. PERSONNEL

Contamination from personnel includes hair, dandruff, skin particles, fibres from clothing, cosmetics and medication. Garments made of monofilament polyester are recommended to avoid the possibility of fibres adding to dust levels. Some garments and footwear also contribute to high static levels. General rules involving cleanliness, such as no eating and smoking in critical areas, must be strictly enforced.

III.10. CEILINGS

The darkroom ceiling should not be of the suspended tile type. A sealed drywall ceiling is required in order to reduce dust and light leaks to a minimum. Suspended ceilings tend to drop dust from the composition tiles,

especially as there tends to be movement in the tiles when the door of the room is opened and closed. Opaque vinyl coated 'sanitary grade' ceiling tiles may be used, if they are clipped down using clips designed to meet fire code specifications and the ceiling track is caulked around the perimeter.

III.11. SINK

A large sink is required, installed at waist height to accommodate the cleaning of the rollers from the processor. This sink should be at least 60 cm wide by 45 cm deep and be provided with both hot and cold water. There should be a laundry tub hose at least long enough to reach into every part of the processor connected to a mixing faucet through a vacuum breaker as required by most local plumbing codes. A large sink is required even for daylight load processors.

III.12. PHYSICAL ARRANGEMENT OF PROCESSOR AND CHEMICALS

There should be an easily accessible area in which to store the replenishment chemistry (normally two tanks, 45 cm diameter) and a provision for the hoses that must reach from the replenisher tanks to the processor. These hoses should be allowed to lie flat on the floor, and not be subject to extreme bending or kinking. A silver recovery unit should be located in this area, and connected to the drain. At least 0.5 m clearance is required on the sides of the processor and 1 m in front and above for service.

III.13. FLOATING LIDS

The use of floating lids on solution storage tanks will extend the storage life of the solution by reducing evaporation and chemical oxidation. Lids are also effective for keeping dirt and dust out of processing solutions. In this way, floating lids reduce the need to discard processing solutions.

III.14. SILVER RECOVERY

The discharge fluid from the processor if allowed to flow directly down the drain will have silver levels which exceed environmental standard limits. Having the fixer enter an electrolytic silver recovery unit, and then allowing the overflow of that unit to pass into a metallic replacement cartridge unit, will ensure that the silver levels are below legal requirements. It is important to regularly (quarterly) replace the units (and harvest the metallic silver).

III.15. ELECTRICAL SERVICES REQUIRED

In the darkroom, there must be an unlighted switch for the white light (1.6 m from the floor) and another unlighted switch for the safe light (1.4 m from floor). The countertop must have at least four outlets, although they can be on the same circuit. Generally, an electrical circuit breaker and switch (isolator or disconnect) is required for the processor. The isolator must be located within 2 m of the processor and can be either inside or outside the darkroom. Connections to the processor should be made through a flexible, liquid-proof cable or conduit.

III.16. PLUMBING SERVICES REQUIRED

Most processors require only a cold water supply, but this should be verified. The temperature and pressure of this supply should be relatively stable, and for some locations, a mixing or tempering valve may be required. If water that is too cold is used in a processor or in an automixer, the film response may be degraded. A 50 μ *in-line water filter*³³ is recommended. This filter should be *changed annually*, or more often as required. A

³³ This filter is not to be confused with the *developer filter* located inside the processor, which should be *replaced monthly*.

drain connection is required for the processor to accommodate chemical and water drainage of about 15 L/min. An oversized drain may prevent overflow due to chemical precipitation and deposits reducing the size of the pipe. Glass, stoneware or chemically inert plastics (PVC) are appropriate materials to use, in order to prevent corrosion from the chemicals. There should also be a floor drain in the area with the processor and silver recovery unit.

III.17. FLOORS

The floors of the darkroom and adjacent area must be non-porous, not slippery when wet, and stain and water-resistant. They must be suitable for frequent cleaning and not collect dust. The adhesive used to secure the flooring must also be waterproof. A static-free mat may be found useful in the area for the loading and unloading of films.

III.18. HUMIDIFIERS: RELATIVE HUMIDITY AND STATIC ELECTRICITY

Static electricity can be the cause of defects in photosensitized materials in two different ways. Firstly, the light produced by the discharge of static electricity can fog the sensitized materials, and secondly, an electrostatic charge on the sensitized materials increases the attraction of dust particles to them. There are several methods for controlling static on film.

Anti-static screen cleaners reduce the amount of charge built up on the screen, cassette and film, but the most effective method for controlling static is to maintain a high relative humidity in the darkroom. A high relative humidity does not affect the conductivity of the air but it does increase the electrical conductivity of nearly all surfaces in any area, for example, the walls, floor and equipment. The higher electrical conductivity helps to drain off static charges as they form.

Corrosion of equipment and tackiness of photographic film generally set an upper limit on the relative humidity that can be used. A relative humidity that is below 30% may cause static attraction of dust to the screens and/or countertops while a relative humidity above 70% may cause the sticking of films.

Steam humidifiers are preferred in a central heating, ventilation and air conditioning (HVAC) system, but evaporative humidifiers work well if they are regularly maintained and have an overflow design. Ultrasonic humidifiers are not recommended for inside or outside darkrooms, as they require constant refilling with distilled water. If tap water is used, fine powder is deposited throughout the facility, as a result of the minerals in the water.

Another method for static reduction is to ionize the air so that static charges can be bled off to the walls, floor and equipment through the air. Air ionizers produce air flow over the needles of a high voltage static bar, thus charging the air molecules. The static bar can be hidden so that the film is not exposed by the associated corona discharge.

III.19. STORAGE OF FILM PROCESSING CHEMICALS

Chemicals may be obtained in either pre-mixed or 'user to mix' format. In order to minimize chemical deterioration associated with ageing and oxidation, especially in the case of pre-mixed chemicals, it is important that large reserve stocks not be kept on hand but rather that supplies be obtained on a frequent and regular basis. The mixed solutions should be stored at room temperature in containers with solid covers to prevent air circulation over the solution.

Appendix IV

TEST EQUIPMENT SPECIFICATIONS

TABLE 7. SPECIFICATIONS OF TEST EQUIPMENT FOR RADIOGRAPHER'S TESTS

Item	Requirements
Breast phantom	Should allow assessment of background optical density, image quality and contrast. Accrediting/regulatory authorities in some jurisdictions may require a specific make and model to be used.
Sensitometer	21 step sensitometer.
Densitometer	Capable of measuring in the range 0–4.0 OD and be accurate to within ± 0.02 OD.
Thermometer (must be non-mercury)	$\pm 0.1^\circ\text{C}$ accuracy (Note: digital thermometers as used in medicine are accurate and inexpensive).
pH meter ^a	Must be temperature compensated and used with calibrated buffer solutions.
Hydrometer	Suitable for measuring specific gravity of processing chemicals.
Silver content estimation paper and fixer retention test kit ^a	Includes hypo-estimator and hypo-test solution.
Screen film contact test tool	40 mesh (or equivalent).

^a Not required for routine QC, but useful for investigating processor problems (see Appendix V).

TABLE 8. SPECIFICATIONS OF TEST EQUIPMENT FOR MEDICAL PHYSICIST'S TESTS

Item	Requirements
Phantom	ACR accreditation phantom or other nationally or internationally recognized breast phantom.
PMMA (AEC)	Able to provide thicknesses of 20–70 mm. Rectangular with dimensions $\geq 150\text{ mm} \times 100\text{ mm}$ or semicircular with radius $\geq 100\text{ mm}$.
PMMA (dosimetry) ^a	Single or multiple pieces to provide a thickness of $45\text{ mm} \pm 0.5\text{ mm}$ (see also Appendix VIII and Glossary). Rectangular with dimensions $\geq 150\text{ mm} \times 100\text{ mm}$ or semicircular with radius $\geq 100\text{ mm}$.
Dosimeter(s) ^b	Mammography energy response of within $\pm 5\%$; accuracy $\pm 5\%$; repeatability $\pm 5\%$. Capable of measuring down to $0.1\ \mu\text{Gy}$ for leakage measurements.
Thermometer and barometer	To be used with dosimeter when automatic temperature and pressure correction is not available.
Lead sheets with apertures of various diameters (optional)	To create narrow beam geometry for HVL measurement.
Aluminium filters	At least 99.9% pure Al and providing thicknesses of 0.1–0.6 mm in increments of 0.1 mm.
kVp meter ^b	Capable of measurement down to 25 kVp with accuracy of $\pm 1\text{ kVp}$ and repeatability of $\pm 0.5\text{ kVp}$.
Resolution pattern	Allowing measurement down to 20 lp/mm.
Radiation timer	Capable of measuring exposure times between 0.1 and 5 s with precision of 5 ms.
Densitometer	Capable of measuring in the range 0–4.0 OD and be accurate to within ± 0.02 OD.
Optical density calibration strip	Current and traceable to an accepted standard.
Magnifier	Capable of at least $2\times$ magnification. It is recommended that the magnifier contain a second area that provides $4\times$ or $5\times$ magnification.
Calibrated photometer	Designed to measure both luminance and illuminance.

^a If the AEC utilizes the thickness measurement from the compression device, it is also useful to have a radiolucent spacer (e.g. U-shaped expanded polystyrene) that is 8 mm thick to be placed between the 45 mm PMMA slab and the compression paddle. In this way, the compression thickness device senses a thickness of 53 mm (standard breast) when the 45 mm thick slab is being imaged in the dosimetry and AEC tests.

^b Regular calibration of dosimeters and kVp meters is essential [28].

Appendix V

AUTOMATIC PROCESSOR DEVELOPER TIME, SPECIFIC GRAVITY, REPLENISHMENT AND pH

V.1. SCOPE

- Objective: To verify the developer chemicals' pH, the developing time, specific gravity and replenishment rate;
- References [3–5, 12, 14, 31];
- Frequency: These tests should be conducted only if problems are detected in routine sensitometry and it is necessary to determine their cause. They are best done by specialized technicians.

V.2. INSTRUMENTATION

- (1) Calibrated, temperature compensated pH meter with required buffer solutions. Note that colour corrected pH paper is of little help in processor QC due to its limited accuracy (about 15%). Please note that silver content estimating paper can be very useful to determine if the fixer replenishment is suitable. Similarly, a hypo-retention test on a sample film will also measure the quality of the wash cycle in a processor and hence the film's archival quality.
- (2) Stopwatch.
- (3) Graduated measuring cylinder 150 ml.
- (4) Plastic containers large enough to contain the replenishment volume.
- (5) Hydrometers.

V.3. METHODOLOGY

pH, specific gravity and development time measurements:

- (1) Confirm that chemical levels are adequate.
- (2) Measure the specific gravity of the developer and fixer in the mixing tank, which is often the replenishment tank, using the hydrometer.
- (3) Measure the pH of the developer and fixer in both the replenisher tanks and in the processor. Record these values on the data collection sheet.
- (4) Place a 'green' film³⁴ in the processor feed tray.
- (5) Activate the stopwatch when the film is dragged by the processor.
- (6) Stop the stopwatch when the film starts leaving the processor through the exit rollers. This is the film system process time. Record the time read on the data collection sheet. **HAZARD WARNING: *special care needed.*** In order to measure the development time, it is usually necessary to remove the processor lid, but to engage the lid microswitch so the processor will still operate (see the following). The development time is the time from when the leading edge of the film enters the development processing rack to when the leading edge then leaves the rack. This test may be done in conjunction with the replenishment rate test below.

Replenishment rate:

- (1) Refer to the process instruction manual for detailed instructions for particular units. If it is not available [31], proceed with care as described in the following.
- (2) Remove the processor lid that covers the developer and fixer sections.

³⁴ This is a film that has not been processed, however, it can be old stock or light damaged film. While previously processed film is sometimes used, there is a danger of contaminating the developer with the use of this film if the wash process has not been complete in removing the fixer during the initial film processing.

- (3) HAZARD WARNING: *special care needed*. Block — or use a magnet to override — the relevant microswitch so that the processor can be used with the cover open.
- (4) Attach one end of the hoses to the replenishment pump outlets and place the other end into the disposable plastic glasses. Note: both developer and fixer rates should be measured. This may be done separately.
- (5) Activate replenishment pumps by developing a film or other appropriate method.³⁵
- (6) Collect each liquid in a plastic glass and measure the volume collected from each with the test tube.
- (7) Record the measurements on the data collection sheet.
- (8) Repeat the test three times to confirm that the same volume ($\pm 10\%$) is measured.
- (9) Compare with the value provided by the manufacturer using the following relation:

$$\text{Volume variation (\%)} = \frac{\text{manufacturer value} - \text{measured value}}{\text{manufacturer value}} \times 100 \quad (11)$$

V.4. INTERPRETATION OF RESULTS AND CONCLUSIONS

Tolerances:

- (1) Total development time: $\pm 3\%$ of the value indicated by the film manufacturer.
- (2) Specific gravity: Manufacturer's specification, typically ± 0.004 . Note that an error of 0.007 represents a dilution error of 10% [12]. Note also that used developer has an increased density compared with that in the replenishment tanks, as a result of evaporation.
- (3) pH: ± 0.5 value indicated by the manufacturer.
- (4) Replenishment rate: $\pm 10\%$ of manufacturer's specification.

V.5. RECOMMENDATIONS AND CORRECTIVE ACTIONS

- (1) If the developing time is found to be outside of tolerances, contact the maintenance service.
- (2) If the specific gravity is outside of tolerances, check the mixing ratios of concentrate and water used for the developer or fixer solutions. An error may be caused by changes of volume in the mixing tank due to distortion of the tank walls. In some cases, the chemical concentrate may be at fault. Contact the maintenance service.
- (3) Note: the pH measurement alone is not a good indicator of processor performance which is best seen from sensitometry charts. If pH levels are outside of tolerances, contact the maintenance service.
- (4) If the replenishing rate is outside of tolerances, contact the relevant service provider immediately. Note that after examination of the monthly sensitometric chart, the maintenance engineer may change the replenishment rate in order to stabilize the speed and/or density difference. Also, as pumps wear, adjustments need to be made to maintain replenishment rates, however, pumps should be replaced with high quality pumps when indicated.

³⁵ Please note that replenishment rates are defined differently for different processing units. They can be defined simply as a rate (how much replenisher is pumped per minute) or as the amount of replenisher pumped per length of film or per area of film. Please note that the rates may be different for the developer and the fixer.

Appendix VI

EXAMPLE OF A CALCULATION FOR FILM CROSSOVER PROCEDURE

Before a film emulsion batch is finished, it is important for the new film emulsion batch to be used for densitometry in parallel with the remaining five films of the old emulsion batch (Tables 9 and 10). In order to keep the same density of the clinical films, the sensitometry operating levels need to be varied by a calculated amount as shown in Table 11.

TABLE 9. EXAMPLE OF OLD AND NEW EMULSION BATCHES

1 Jan. 2006	Batch identification
Brand, model and number of old emulsion	Kodak abc
Brand, model and number of new emulsion	Kodak xyz

TABLE 10. EXAMPLE OF SENSITOMETRIC DATA FOR BOTH THE OLD AND NEW EMULSION BATCHES

Day	Old emulsion						New emulsion						Difference between average values
	1	2	3	4	5	Average	1	2	3	4	5	Average	
S ^a	1.50	1.49	1.51	1.52	1.48	1.50	1.51	1.52	1.49	1.53	1.52	1.51	0.01
HD ^b	2.60	2.50	2.64	2.52	2.65	2.58	2.70	2.68	2.64	2.70	2.60	2.66	0.08
DD ^c	1.10	1.03	1.13	1.00	1.17	1.08	1.21	1.16	1.15	1.17	1.08	1.15	0.07
B+F ^d	0.22	0.19	0.21	0.18	0.18	0.20	0.21	0.20	0.17	0.18	0.20	0.19	-0.01

^a S: Speed index.

^b HD: High density.

^c DD: Density difference.

^d B+F: Base plus fog level.

TABLE 11. EXAMPLE OF DETERMINATION OF NEW OPERATING LEVELS WITH THE NEW EMULSION BATCH

	S ^a	DD ^b	B+F ^c
Old initial operating level ^d	1.52	1.10	0.21
Difference between old and new ^e	0.01	0.07	-0.01
New initial operating level	1.53	1.17	0.20

^a S: Speed index.

^b DD: High density.

^c B+F: Base plus fog level.

^d Hypothetical example.

^e See Table 10.

Appendix VII

CALCULATION OF AVERAGE FILM GRADIENT

The film gradient between two densities, D_1 and D_2 , is given by the following equation:

$$\text{Film gradient} = \frac{D_2 - D_1}{\log E_2 - \log E_1} \quad (12)$$

where D_1 and D_2 are densities that are usually chosen to be 1.0 and 2.0 above the base plus fog level (B+F). E_1 and E_2 are the corresponding exposure levels.

Most modern sensitometers have a step size that is 0.15 in terms of log exposure units. The calculation is illustrated (Fig. 47 and Table 12) using the sensitometric data (Table 13). The values of the log exposure can be approximated by linear interpolation between the data points provided.




FIG. 47. Sensitometric data used for estimation of film gradient.

TABLE 12. CALCULATION OF FILM GRADIENT

	OD ^a	Log exposure
B+F ^b	0.22	—
D ₁ ^c	1.22	Log E ₁ ^d = 1.63
D ₂	2.22	Log E ₂ = 1.81
Film gradient = (2.22-1.22)/(1.81-1.63)		
= 5.55		

^a OD: Optical density.

^b B+F: Base plus fog level.

^c D: Density.

^d E: Exposure level.

TABLE 13. SENSITOMETRIC DATA

Step	Log exposure	OD ^a
1	0.00	0.22
2	0.15	0.22
3	0.30	0.22
4	0.45	0.22
5	0.60	0.22
6	0.75	0.23
7	0.90	0.25
8	1.05	0.29
9	1.20	0.38
10	1.35	0.51
11	1.50	0.77
12	1.65	1.30
13	1.80	2.17
14	1.95	3.10
15	2.10	3.82
16	2.25	4.21
17	2.40	4.21
18	2.55	4.06
19	2.70	4.13
20	2.85	4.18
21	3.00	4.20

^a OD: Optical density.

Appendix VIII

THE STANDARD BREAST

For testing purposes, the 'standard breast' is simulated by PMMA of thickness 45 ± 0.5 mm with a semicircular section with a radius of at least 100 mm or rectangular with dimensions of at least 150 mm \times 100 mm [4, 28]. The attenuation is equivalent to a compressed breast of 53 mm thickness of which the central thickness is composed of 29% fibroglandular tissue and 71% fat tissue and this region is surrounded above and below by a layer of fat that is 5 mm thick [25].

REFERENCES

- [1] WORLD HEALTH ORGANIZATION, IARC Handbooks of Cancer Prevention, Vol. 7, Breast Cancer Screening (INTERNATIONAL AGENCY FOR RESEARCH ON CANCER, Ed.), IARC Press, Lyons (2002).
- [2] TABAR, L., et al., Mammography service screening and mortality in breast cancer patients: 20 year follow up before and after introduction of screening, *Lancet* **361** (2003) 1405–1410.
- [3] AMERICAN COLLEGE OF RADIOLOGY, Mammography Quality Control Manual, American College of Radiology, Reston, VA (1999).
- [4] EUROPEAN COMMISSION, European Guidelines for Quality Assurance in Mammography Screening, 4th edn, Europe against Cancer, Office for Official Publications of the European Communities, Luxembourg (2006).
- [5] INTERNATIONAL ATOMIC ENERGY AGENCY, Control de Calidad en Mamografía, IAEA-TECDOC-1517, IAEA, Vienna (2006).
- [6] SOCIEDAD ESPAÑOLA DE FÍSICA MÉDICA–SOCIEDAD ESPAÑOLA DE PROTECCIÓN RADIO-LÓGICA, Protocolo Español de Control de Calidad en Radiodiagnóstico (Aspectos Técnicos), 2nd edn, Edicomplet, Madrid (2002).
- [7] PISANO, E.T., et al., Diagnostic Performance of Digital versus Film Mammography for Breast-Cancer Screening, *N. Engl. J. Med.* **353** (2005) 1–11.
- [8] INTERNATIONAL ELECTROTECHNICAL COMMISSION, Medical Electrical Equipment: X Ray Tube Assemblies for Medical Diagnosis — Characteristics of Focal Spots, Rep. IEC 60336, IEC, Geneva (2005).
- [9] WORLD HEALTH ORGANIZATION, Quality Assurance in Diagnostic Radiology, WHO, Geneva (1982).
- [10] INTERNATIONAL ATOMIC ENERGY AGENCY, Applying Radiation Safety Standards in Diagnostic Radiology and Interventional Procedures Using X Rays, Safety Reports Series No. 39, IAEA, Vienna (2006).
- [11] MOORES, B.M., HENNSHAW, E.T., WATKINSON, S.A., PEARCY, B.J., Practical Guide to Quality Assurance in Medical Imaging, John Wiley & Sons, Chichester (1987).
- [12] MCKINNEY, W.E.J., Radiographic Processing and Quality Control, JB Lippincott Co., Philadelphia, PA (1988).
- [13] INTERNATIONAL ELECTROTECHNICAL COMMISSION, Evaluation and Routine Testing in Medical Imaging Departments, Part 2–1: Consistency Tests — Film Processors, Rep. IEC 61223-2-1, IEC, Geneva (1993).
- [14] HAUS, A.G. (Ed.), Advances in Film Processing Systems Technology and Quality Control in Medical Imaging, Medical Physics Publishing, Madison, WI (2001).
- [15] AMERICAN ASSOCIATION OF PHYSICISTS IN MEDICINE, Equipment Requirements and Quality Control for Mammography, Rep. 29, AAPM, New York (1990).
- [16] GRAY, J.E., “Mammographic quality control for the technologists and the medical physicist as consultant to the technologists”, Syllabus: A Categorical Course in Physics — Technical Aspects of Breast Imaging (HAUS, A.G., YAFFE, M.J., Eds), Radiological Society of North America (1992) 131–160.
- [17] INTERNATIONAL ELECTROTECHNICAL COMMISSION, Evaluation and Routine Testing in Medical Imaging Departments, Part 2. 10 Constancy tests — X ray equipment for mammography, Rep. IEC 61223-2-10, IEC, Geneva (1999).
- [18] GRAY, J.E., “Mammographic quality control for the technologists and the medical physicist as consultant to the technologists”, Syllabus: A Categorical Course in Physics — Technical Aspects of Breast Imaging (HAUS, A.G., YAFFE, M.J., Eds), Radiological Society of North America (1994) 153–183.
- [19] WIDMER, J.H., LILLIE, K.R., JASHULSKI, S.M., HAUS, A.G. (Eds), Identifying and Correcting Processing Artefacts, Technical and Scientific Monograph No. 4, Kodak, Health Sciences, Rochester, NY (1994).
- [20] INTERNATIONAL ELECTROTECHNICAL COMMISSION, Evaluation and Routine Testing in Medical Imaging Departments, Part 3–2: Acceptance Testing — Imaging Performance of Mammographic X-ray Equipment, Rep. 61223-3-2, edn 2.0, IEC, Geneva (2004).
- [21] YOUNG, K.C., WALLIS, M.G., BLANKS, R.G., MOSS, S.M., Influence of number of views and mammographic film density on the detection of invasive cancers: Results from the NHS Breast Screening Programme, *Br. J. Radiol.* **70** (1997) 482–488.
- [22] BRITISH INSTITUTE OF RADIOLOGY, Assurance of Quality in the Diagnostic X Ray Department, QA Working Group of the Diagnostic Methods Committee of the BIR, London (1988).
- [23] WAGNER, L.K., ARCHER, B.R., CERRA, F., On the measurement of half-value layer in film-screen mammography, *Med. Phys.* **17** (1990) 989–997.
- [24] INTERNATIONAL ELECTROTECHNICAL COMMISSION, Medical Electrical Equipment, Part 2–45: Particular Requirements for the Safety of Mammographic X-ray Equipment and Mammographic Stereotactic Devices, Rep. IEC 60601-2-45, IEC, Geneva (1998).
- [25] DANCE, D.R., SKINNER, C.L., YOUNG, K.C., BECKETT, J., KOTRE, C.J., Additional factors for the estimation of mean glandular breast dose using the UK mammographic dosimetry protocol, *Phys. Med. Biol.* **45** (2000) 3225–3240.

- [26] WU, X., BARNES, G.T., TUCKER, D.M., Spectral dependence of glandular tissue dose in screen-film mammography, *Radiology* **179** (1991) 143–148.
- [27] WU, X., GINGOLD, E.L., BARNES, G.T., TUCKER, D.M., Normalised average glandular tissue dose in molybdenum target – rhodium filter and rhodium – target rhodium filter mammography, *Radiology* **193** (1994) 83–89.
- [28] INTERNATIONAL ATOMIC ENERGY AGENCY, *Dosimetry in Diagnostic Radiology: An International Code of Practice*, Technical Reports Series No. 457, IAEA, Vienna (2007).
- [29] INTERNATIONAL COMMISSION ON RADIATION UNITS AND MEASUREMENTS, *Patient Dosimetry for X Ray Used in Medical Imaging*, Rep. No. 74, ICRU, Bethesda, MD (2006).
- [30] EASTMAN KODAK COMPANY, *KODAK X-OMAT Processor*, Service Bulletin No. 158, Eastman Kodak Company/Customer & Technical Services/Health Sciences Markets Division/Rochester, NY (1993).
- [31] EASTMAN KODAK COMPANY, *Processing Recommendations*, Service Bulletin No. 30, Eastman Kodak Company/Health Group, Rochester, New York (2007),
<http://www.carestreamhealth.com/en/global/plugins/acrobat/en/health/support/tib3107.pdf>

BIBLIOGRAPHY

- AMERICAN ASSOCIATION OF PHYSICISTS IN MEDICINE, Equipment Requirements and Quality Control for Mammography, AAPM Rep. 29, AAPM, New York (1990).
- AMERICAN COLLEGE OF RADIOLOGY, Mammography Quality Control Manual, ACR, Reston, VA (1999).
- BRITISH INSTITUTE OF RADIOLOGY, Assurance of Quality in the Diagnostic X Ray Department, 2nd edn, QA Working Group of the Radiation Protection Committee of the BIR, London (2001).
- DANCE, D.R., SKINNER, C.L., YOUNG, K.C., BECKETT, J., KOTRE, C.J., Additional factors for the estimation of mean glandular breast dose using the UK mammographic dosimetry protocol, *Phys. Med. Biol.* **45** (2000) 3225–3240.
- EUROPEAN COMMISSION, European Guidelines for Quality Assurance in Mammography Screening, 4th edn, Europe against Cancer, Office for Official Publications of the European Communities, Luxembourg (2006).
- GUNN, C., *Radiographic Imaging: A Practical Approach*, 3rd edn, Churchill Livingstone, Edinburgh (2002).
- HAUS, A.G. (Ed.), *Advances in Film Processing Systems Technology and Quality Control in Medical Imaging*, Medical Physics Publishing, Madison, WI (2001).
- INSTITUTE OF PHYSICS AND ENGINEERING IN MEDICINE, *The Commissioning and Routine Testing of Mammographic X Ray Systems*, IPEM Rep. 89, IPEM, York (2005).
- LLOYD, P.J., *Quality Assurance Workbook for Radiographers and Radiological Technologists*, World Health Organization, Geneva (2001).
- McCLELLAND, I.R., *X Ray Equipment Maintenance and Repairs Book for Radiographers and Radiological Technologists*, World Health Organization, Geneva (2004).
- WIDMER, J.H., LILLIE, K.R., JASHULSKI, S.M., HAUS, A.G. (Eds), *Identifying and Correcting Processing Artefacts*, Technical and Scientific Monograph No. 4, Kodak, Health Sciences, Rochester, NY (1994).

Annex I

RADIOGRAPHER'S DATA COLLECTION SHEETS

List of data collection sheets in this annex:

- 6 Unit Information
 - 6.1 Visual inspection checklist
 - 6.2 Film storage
 - 6.3.1 Darkroom cleanliness
 - 6.3.2 Darkroom temperature, humidity & ventilation conditions
 - 6.3.3 White light leakage and safe lights
 - 6.3.4 Automatic processor developer temperature
 - 6.3.5 Sensitometry (2 sheets)
 - 6.3.5 X ray processing control chart
 - 6.3.6 Artefact detection
 - 6.3.7 Transition between emulsion numbers
 - 6.4.1 Cleaning of intensifying screens
 - 6.4.2 Screen film contact
 - 6.4.3 Uniformity of cassette sensitivity, attenuation and artefacts
 - 6.5.1 Automatic exposure control consistency
 - 6.5.1 X ray system AEC consistency chart
 - 6.5.2 AEC thickness compensation
 - 6.5.2 X ray system AEC thickness compensation chart
 - 6.6 Evaluation of image quality
 - 6.7 Film rejection rate
- Appendix III Automatic processor developer time, pH and replenishment rate, etc.

The electronic version of these data collection sheets can be found by accessing the IAEA web site at http://www-pub.iaea.org/MTCD/publications/PDF/Pub1381Files/Annex1_DataCollectionSheets_Radiographer.xls

Unit Information

Name of Unit:	
Unit no.:	
Date:	

Mammography unit manufacturer:		
Mammography unit model:		
Anode material (indicate all possible options)	Anode 1	
	Anode 2	
Filter material (indicate all possible types and thicknesses):	Filter 1	
	Filter 2	
Grid (lines/cm):		
Automatic exposure control (Yes/No):		

Processor manufacturer:	
Processor model:	

Screen manufacturer:	
Screen type:	

Film manufacturer:	
Film type:	

6.1 Visual inspection checklist

		Month	Year	
Parameter			Yes	No
1	Free standing unit is mechanically stable.			
2	Indicator lights working properly.			
3	All moving parts move smoothly, without obstructions to motion.			
4	All locks and detents work properly.			
5	Angulation indicators function properly.			
6	Image receptor and holder is free from vibrations during exposure.			
7	Image receptor slides smoothly into holder assembly.			
8	Image receptor is held securely by assembly in any orientation.			
9	The compression plate is in good condition.			
10	The compression breast thickness scale (analog or digital) is accurate and reproducible.			
11	The automatic compression release following exposure functions correctly.			
12	The manual release of compression is possible when power fails.			
13	The compression release override works properly.			
14	The radiation shield for the operator is adequate.			
15	There are no sharp edges on the breast support or compression paddle.			
16	The face guard is in place.			

6.3.1 Darkroom cleanliness

Month			Year	
Day	Comment	Floor	Surfaces	Processor
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				
16				
17				
18				
19				
20				
21				
22				
23				
24				
25				
26				
27				
28				
29				
30				
31				

Comments:

6.3.4. Automatic processor developer temperature

Month		Year	
Manufacturer value (°C)			
Processor display temperature (°C)			
Tolerance (°C)	±1.0°C		

Day	Measured value (°C)	Acceptable
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
15		
16		
17		
18		
19		
20		
21		
22		
23		
24		
25		
26		
27		
28		
29		
30		
31		

6.3.5 Sensitometry

Establishing initial operating levels (IOLs)						
Step	Optical density (OD)					Average
	Day 1	Day 2	Day 3	Day 4	Day 5	
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18						
19						
20						
21						
Initial Operating Levels						
	S (OD closest to but not <1.2)	HD (OD closest to but not <2.2)	DD	B+F		
Step				1		
OD						
	Tolerance (B+F)		Acceptable	≤0.25		
			Achievable	≤0.21		
			Acceptable			

Comments:

Parameter	IOL	Tolerance	
		Acceptable	Achievable
B+F		$\leq \text{IOL} + 0.03$	$\leq \text{IOL} + 0.02$
S		$\text{IOL} \pm 0.15$	$\text{IOL} \pm 0.10$
DD		$\text{IOL} \pm 0.15$	$\text{IOL} \pm 0.10$

Daily sensitometry							
Day	Optical density (OD)						
	B+F		S		HD	DD	
	Measured	Acceptable	Measured	Acceptable	Measured	Calculated	Acceptable
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							
17							
18							
19							
20							
21							
22							
23							
24							
25							
26							
27							
28							
29							
30							
31							

6.3.7 Transition between emulsion numbers

Date:

	Old emulsion	New emulsion
Brand		
Model		
Emulsion number		

Day	Old emulsion						New emulsion					
	1	2	3	4	5	average	1	2	3	4	5	average
S												
HD												
DD												
B+F												

	S	DD	B+F
Old initial operation levels			
Difference between new and old			
New initial operation levels			

Comments:

6.5.1 Automatic exposure control consistency						
Density control setting			AEC detector position			
Target			Filter			
kVp			Phantom		45 mm PMMA	
Cassette ID			Film type			
100			Baseline mAs			
Year						
	OD		mAs		Artefacts	
Date	Measured	Acceptable	Displayed	Acceptable	Y/N	Acceptable 1
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
12						
13						
14						
15						
16						
17						
18						
19						
20						
21						
22						
23						
24						
25						
26						
27						
28						
29						
30						
31						

Tolerance		
	Achievable	Acceptable
OD _{target} ± 0.10	OD _{target} ± 0.20	
mAs	±10% of baseline mAs	



Comments:

--

X-RAY SYSTEM AEC CONSISTENCY CHART
Department of Diagnostic Radiology

_____ Casette: _____ Month: _____ Year: _____
 _____ Density setting: _____ AEC detector position: _____
 _____ Anode/Filter: _____ Phantom: _____

Day: _____
 Month: _____
 Year: _____



Notes

6.5.2 AEC thickness compensation			
Density control setting		AEC detector position	
Cassette ID		Film type	
Month		Year	

Phantom thickness 25 mm				
Target		Filter		
kVp				
OD _{target}		Baseline mAs		
	OD		mAs	
Week	Measured	Acceptable	Displayed	Acceptable
1				
2				
3				
4				

Phantom thickness 45 mm				
Target		Filter		
kVp				
Target optical density		Baseline mAs		
	OD		mAs	
Week	Measured	Acceptable	Displayed	Acceptable
1				
2				
3				
4				

Phantom thickness 70 mm				
Target		Filter		
kVp				
Target optical density		Baseline mAs		
	OD		mAs	
Week	Measured	Acceptable	Displayed	Acceptable
1				
2				
3				
4				

Tolerance		
	Achievable	Acceptable
OD	OD _{target} ± 0.15	OD _{target} ± 0.20
mAs	±10% of baseline mAs for each thickness	




Comments:

--

X-RAY SYSTEM AEC THICKNESS COMPENSATION CHART
Department of Diagnostic Radiology

Room: _____ Cassette: _____ Year: _____
 AEC Mode: _____ Density setting: _____ AEC detector position: _____

Day: _____
 Month: _____
 Initials: _____



Comments

Date	Comment	Date	Comment

6.6 Evaluation of image quality

Phantom used			
Cassette size			
Cassette ID			
	Previous film	Current film	Change
Date			
Operation mode (manual, semiautomatic, automatic)			
Target/filter			
kVp setting			
AEC detector position:			
Density control			
Target optical density (1.5–1.9)			
mAs			
Background OD (central location)			
Target OD - Background OD			
OD outside disc (OD ₁)			
OD inside disc (OD ₂)			
OD₁ - OD₂ (image contrast)			
Number of fibres seen			
Number of fibre-like artefacts			
Fibres seen after deduction			
Number of microcalc. groups seen			
Number of micro-like artefacts			
Speck groups after deduction			
Number of masses seen			
Number of masses-like artefacts			
Masses seen after deduction			
Tolerance			
Image contrast index: Achievable: ≥ 0.55; Acceptable: ≥ 0.40			
Pass?			
Fibres: ≥ 4; microcalcifications: ≥ 3; masses: ≥ 3			
Pass (Y/N)?			
Change in fibres, microcalcifications, masses ≤ 0.5			
Pass (Y/N)?			

Comments:

6.7 Film rejection rate

Person responsible							
Period of analysis		from				to	
Total of films used for the period				Total studies			
Cause	Frequency of rejection cause						
	CC		MLO		Other		Total
	Left	Right	Left	Right	Left	Right	
Technical factors:							
Positioning							
Motion							
Incorrect labelling							
Technique factors							
Artefacts							
Darkroom and processor							
Static							
Fog							
Light and dark films							
Artefacts							
Others:							
Total of rejected films							
QC films:		Rejection rate (%)*:					

Tolerance			
Acceptable	≤8%	Acceptable	
Achievable	≤3%		

***Reject rate (%) = 100 * Total reject films / (Total of films used for period - QC films)**

Comments:

Appendix III Automatic processor developer time, pH and replenishment rate, etc.

Developer time	Manufacturer value	Measured	Tolerance	Acceptable
Total processing time (s)			±3%	
Developer time (s)			±3%	
Specific gravity	Manufacturer value	Measured	Tolerance	Acceptable
Measured value			±0.004	
pH	Manufacturer value	Measured	Tolerance	Acceptable
pH developer			±0.5	
pH fixer			±0.5	
Developer replenishment rate (volume)	Manufacturer value	Measured	Tolerance	Acceptable
Trial 1			±10%	
Trial 2			±10%	
Trial 3			±10%	
Fixer replenishment rate (volume)	Manufacturer value	Measured	Tolerance	Acceptable
Trial 1			±10%	
Trial 2			±10%	
Trial 3			±10%	

$$\text{Volume variation (\%)} = \frac{\text{manufacturer value} - \text{measured value}}{\text{manufacturer value}} \times 100$$

Comments:

Annex II

MEDICAL PHYSICIST'S DATA COLLECTION SHEETS

List of data collection sheets in this annex:

- 7 Mammography equipment evaluation (2 sheets)
- 7 Mammography equipment information
- 7.1 Unit assembly evaluation
- 7.2.1 Sensitometry (2 sheets)
- 7.3.1 Radiation leakage
- 7.3.2 Accuracy and repeatability of kVp
- 7.3.3 Half-value layer (HVL)
- 7.3.4 Output repeatability and linearity
- 7.3.5 Compression
- 7.3.6 Evaluation of the automatic exposure control (3 sheets)
- 7.4 Collimation system
- 7.5 Image viewing conditions
- 7.6.1 Evaluation of image quality
- 7.6.2 System spatial resolution
- 7.7.1 Incident kerma at entrance surface of the phantom
- 7.7.2 Mean glandular dose (D_G)

The electronic version of these data collection sheets can be found by accessing the IAEA web site at http://www-pub.iaea.org/MTCD/publications/PDF/Pub1381Files/Annex2_DataCollectionSheets_MedicalPhysicist.xls

Mammography equipment evaluation

Site	
X ray unit manufacturer	
Date of installation	
Medical physicist	
Radiographer	
Type of inspection*	

(*Annual; semi-annual; monthly)

Report date	<input style="width: 100%;" type="text"/>
Survey date	<input style="width: 100%;" type="text"/>
Model	<input style="width: 100%;" type="text"/>
Room ID	<input style="width: 100%;" type="text"/>
Signature	<input style="width: 100%;" type="text"/>

UNIT ASSEMBLY EVALUATION

Pass

DARKROOM AND SENSITOMETRY

Sensitometry

Darkroom radiation level < 20 μ Gy/week

RADIOLOGICAL EQUIPMENT

Radiation leakage < 1 mGy/h at 1 m

Accuracy and repeatability of the tube kV

Accuracy: $\leq 5\%$

Repeatability: Difference $\leq 5\%$ or COV $\leq 2\%$

Half-value layer

HVL is within acceptable lower and upper limits at all kVp values tested

Output: Repeatability and linearity

Repeatability: Difference $\leq 5\%$ or COV $\leq 5\%$

Linearity: $\pm 10\%$

Normalized output value: >30 μ Gy/mAs at 1 m, 28 kV, Mo/Mo

COMPRESSION

Compression force: Power Comp: 150 N–200 N; Manual Comp: <300 N

Compression force: Accuracy $\leq \pm 20$ N

Compression thickness accuracy: ≤ 5 mm

AUTOMATIC EXPOSURE CONTROL

Density Control in "Normal" or "0" position provides OD_{target} ± 0.20

Repeatability of the automatic exposure control: mAs: COV $\leq 5\%$

Compensation for different kV and thicknesses (and target/filter combination):

Acceptable: Maximum deviation from OD target ≤ 0.2 ;

Achievable: Maximum deviation from OD target ≤ 0.15

Density control setting: difference per step: %mAs change: 12%–15%

Δ OD : 0.1–0.2

Exposure time 45 mm thick PMMA test

Contact mammography:

Achievable: $t \leq 1.5$ s

Acceptable: $t \leq 2$ s

Magnification mammography:

Achievable: $t \leq 2$ s

Acceptable: $t \leq 3$ s

Mammography equipment evaluation (cont.)

COLLIMATION SYSTEM

Light field/radiation field coincidence: $\leq 1\%$ of FFD on any side

Radiation field/image receptor coincidence

Acceptable: Chest side 0 to ≤ 5 mm; $\leq 2\%$ of FFD for other three sides

Achievable: Chest side 0 to ≤ 5 mm; ≥ 0 for other three sides

Compression paddle/breast support alignment: $\leq +1\%$ of FFD

IMAGE VIEWING CONDITIONS

Luminance of the viewboxes: > 3000 cd/m² (nit)

Viewboxes homogeneity and colour

Single viewbox: $< 30\%$

Between viewboxes: $< 15\%$

Ambient interpretation room illumination: Achievable: ≤ 10 lux; Acceptable: ≤ 50 lux

IMAGE QUALITY

Phantom image quality evaluation (ACR)

Background optical density: $OD_{\text{target}} \pm 0.2$

Optical density difference between disk and background

Acceptable: ≥ 0.40 OD

Achievable: ≥ 0.55 OD

Phantom score

Acceptable: Fibre score: ≥ 4

Speck score: ≥ 3

Mass score: ≥ 3

Change in score for each group: ≤ 0.5

System spatial resolution

Achievable: ≥ 15 lp/mm

Acceptable: ≥ 11 lp/mm

DOSIMETRY

Mean glandular dose (D_G)

Achievable: $D_G \leq 2.0$ mGy;

Acceptable: $D_G \leq 2.5$ mGy

Mammography equipment information

Site	
Radiographer	

Equipment

Room ID X-ray unit manufacturer: X-ray tube manufacturer: Processor manufacturer: Screen manufacturer: Film manufacturer: Focus size: Contact Targets: AEC sensor positions: AEC sensor for large cassette:		Date: Model: Serial#: Model: Type: Type: Magnif.: Filters: AEC modes:	
---	--	---	--

Mammography phototimer technique chart

Compressed breast thickness (cm)	Fatty breast			50% Fatty - 50% Dense			Dense breast		
	Target & filter	kVp	Density	Target & filter	kVp	Density	Target & filter	kVp	Density
<3									
3 to 5									
5 to 7									
>7									

7.1 Unit assembly evaluation

Parameter	Pass (Y/N/NA)
1 Free standing unit is mechanically stable.	
2 Indicator lights work properly.	
3 All moving parts move smoothly, without obstructions to motion.	
4 All locks and detents work properly.	
5 Angulation indicators function properly.	
6 Image receptor and holder is free from vibrations during exposure.	
7 Image receptor slides smoothly into holder assembly.	
8 Image receptor is held securely by assembly in any orientation.	
9 The compression plate is in good condition.	
10 The compression breast thickness scale (analog or digital) is accurate and reproducible.	
11 The automatic compression release following exposure functions correctly.	
12 The manual release of compression is possible when power fails.	
13 The compression release override works properly.	
14 The radiation shield for the operator is adequate	
15 There are no sharp edges on the breast support or compression paddle.	
16 The face guard is in place.	
17 Confirm that appropriate, current operator technique control charts are posted.	

Unit assembly acceptable (Y/N)?

Comments:

7.2.1 Sensitometry

Establishing initial operating levels (IOLs)							
Step	log E	Optical density (OD)					
		Day 1	Day 2	Day 3	Day 4	Day 5	Average
1	0.00						
2	0.15						
3	0.30						
4	0.45						
5	0.60						
6	0.75						
7	0.90						
8	1.05						
9	1.20						
10	1.35						
11	1.50						
12	1.65						
13	1.80						
14	1.95						
15	2.10						
16	2.25						
17	2.40						
18	2.55						
19	2.70						
20	2.85						
21	3.00						
Initial Operating Levels							
		S (OD closest to but not < 1.2)	HD (OD closest to but not < 2.2)	DD	B+F		
Step							
OD							
		Tolerance (B+F)		Acceptable	≤0.25		
				Achievable	≤0.21		
				Acceptable (Y/N)?			
Gradient estimation							
	$D_1 = 1.0 + (B+F)$	$D_2 = 2.0 + (B+F)$	log E ₁	log E ₂	Gradient		

Log E₁ and log E₂ is obtained by linear interpolation of data points near relevant density steps

Comments:

7.2.1 Sensitometry (cont.)

Annual Site Visit		
Step	log E	OD
1	0.00	
2	0.15	
3	0.30	
4	0.45	
5	0.60	
6	0.75	
7	0.90	
8	1.05	
9	1.20	
10	1.35	
11	1.50	
12	1.65	
13	1.80	
14	1.95	
15	2.10	
16	2.25	
17	2.40	
18	2.55	
19	2.70	
20	2.85	
21	3.00	

Operating levels				
	S (OD closest to but not < 1.2)	HD (OD closest to but not < 2.2)	DD	B+F
Step				
OD				
IOL				
Acceptable				

Parameter	Tolerance	
	Acceptable	Achievable
B+F	$\leq \text{IOL} + 0.03$	$\leq \text{IOL} + 0.02$
S	$\text{IOL} \pm 0.15$	$\text{IOL} \pm 0.10$
DD	$\text{IOL} \pm 0.15$	$\text{IOL} \pm 0.10$

Gradient Estimation				
$D_1 = 1.0 + (B+F)$	$D_2 = 2.0 + (B+F)$	$\log E_1$	$\log E_2$	Gradient

Sensitometry acceptable (Y/N)?

7.3.1 Radiation leakage

Dosimetry system used		Energy correction factor (N_{mammo})		1.00
Pressure (mbar)	Temperature (°C)	Pressure and temperature correction factor (k_{TP})		1.00
mAs ₀ used:	50	mAs _{max}	14400	
Measurements				
Point	Focal spot-chamber distance, d (m)	Reading (mGy)	K ₀ (mGy)	kerma rate at 1m (mGy/h)
1				
2				
3				
4				
Tolerance				
\dot{K}_{RF}		<1 mGy/h at 1 m		
Pass (Y/N)?				

$$\dot{K}_{\text{RF}} = \frac{\text{mAs}_{\text{max}}}{\text{mAs}_0} \times K_0 \times d^2$$

Comments:

7.3.2 Accuracy and repeatability of kVp

kVp meter used:					
Meter setting:					
Nominal kVp setting					
Focal spot (large/small)					
mA setting					
mAs setting					
Measured kVp values					
kVp1					
kVp2					
Repeatability: difference (%)					
Additional measurements (if needed)					
kVp3					
kVp4					
kVp5					
Mean kVp <kVp>					
Standard deviation (SD)					
Repeatability: COV (%)					
Repeatability: acceptable?					
Nominal kVp - <kVp>					
Accuracy (%)					
Tolerance					
Pass (Y/N)?	Repeatability: Difference \leq 5% or			Accuracy: \pm5%	

Comments:

7.3.3 Half-value layer

Dosimetry system used					
Units (mGy, mR)					
Nominal kVp setting					
Anode					
Filter					
Parameter, C					
mAs setting:					
Air kerma or exposure measurements:					
No aluminium filtration, M ₀					
0.2 mm of added aluminium, M ₁					
0.3 mm of added aluminium, M ₂					
0.4 mm of added aluminium, M ₃					
0.5 mm of added aluminium, M ₄					
0.6 mm of added aluminium, M ₅					
Repeat no aluminium filtration, M ₀					
Average no aluminium filtration, M ₀					
Record thicknesses (t _a <t _b) and air kerma or exposure values that bracket k ₀ /2: (K _a >K _b)					
t _a					
t _b					
M _a					
M _b					
Calculated HVL (mm Al)					
Tolerance					
Minimum allowed HVL (mm Al)					
Maximum allowed HVL (mm Al)					
Acceptable (HVL)					
All HVLs acceptable					

$$\text{Calculated HVL} = \frac{t_b \ln[2M_a/M_0] - t_a \ln[2M_b/M_0]}{\ln[K_a/K_b]}$$

Minimum allowed HVL (mm Al):	kVp/100 + 0.03 (in mm Al)
Maximum allowed HVL (mm Al):	(kVp/100) + C (in mm Al)
where C = 0.12 for Mo/Mo, 0.19 for Mo/Rh, 0.22 for Rh/Rh, and 0.30 for W/Rh	

Comments:

7.3.4 Output repeatability and linearity						
Dosimetry system used				Energy correction factor (N _{mammo})		
Units (enter "mGy" or "mR")				Focus-detector distance (cm)		
P(mbar)		Temp (°C)		T ₀ (°C)		P ₀ (mbar)
Pressure and temperature correction factor (k _{TP})* (for auto correction = 1)						1.00
Focus size		Large				
Anode		Mo				
Filter		Mo				
Nominal kVp setting		28				
mAs ₁ =	40	R1				
		R2				
		R3				
		R4				
		R5				
Repeatability: Difference (%)						
Average value						
Standard deviation						
Repeatability: COV (%)						
Output (Y ₁)						
mAs ₂ =	80	R1				
		R2				
Average value						
Output (Y ₂)						
mAs ₃ =	120	R1				
		R2				
Average value						
Output (Y ₃)						
Linearity		L ₁				
		L ₂				
Normalized output (μGy/mAs at 1m)						

Tolerance				
Repeatability @ 28 kVp: Difference ≤ 5% or COV ≤ 5%				
Linearity: <±10%				
Normalized output @ 28 kV with Mo/Mo: >30 μGy/mAs				

$$k_{TP} = \frac{(273.2 + T) P_0}{(273.2 + T_0) P}$$

Comments:

--

7.3.5 Compression

Compression mode		Powered				
Imaging mode	Contact	Magnification				
Maximum force measured						
Displayed force value						
Displayed value accuracy						
Compression mode		Manual				
Imaging mode	Contact	Magnification				
Maximum force measured						
Displayed force value						
Displayed value accuracy						
Thickness measurements						
Compression mode (powered or manual)						
PMMA thickness (mm)	20	45	70	20	45	70
Measured thickness (mm)						
Displayed thickness (mm)						
Displayed accuracy (mm)						
Tolerance						
Maximum force measured: Power compression: 150–200 N; Manual compression: <300 N						
Acceptable (Y/N)?						
Displayed force accuracy		±20 N				
Acceptable (Y/N)?						
Displayed thickness accuracy		≤5 mm				
Acceptable (Y/N)?						
Does breast thickness indication depend on degree of compression (Y/N)?						

Comments:

7.3.6 Evaluation of automatic exposure control (AEC)			
Automatic exposure mode (at clinical conditions)			
Density control	0 (N)	OD _{target} (1.5–1.9)	
Small Cassette ID		Large Cassette ID	
AEC sensor position			
Different sensor for large and small cassette (Y/N)?			
AEC Repeatability (45 mm PMMA)			
	Contact mode		Magnification
	Small cassette	Large cassette	
Target filter:			
kVp:			
mAs1			
mAs2			
mAs3			
mAs4			
Mean (mAs)			
COV (%)			
OD (last image)			
Tolerances			
OD = OD_{target} ± 0.2			
COV ≤ 5%			
Exposure time 45 mm thick test			
	Contact mode		Magnification
	Small cassette	Large cassette	
mAs			
mA			
Time (s) (measured)			
Time (s) (calculated)			
Tolerances			
Acceptable	2.0 s		3.0 s
Achievable	1.5 s		2.0 s
Acceptable (Y/N)?			
Achievable (Y/N)?			

Comments:

--

7.3.6 Evaluation of automatic exposure control (AEC) (cont.)

Compensation for different thicknesses (and kVp and target/filter combination)							
Contact mode				Small cassette			
OD_{target}							
Density control				AEC Sensor position			
PMMA thickness (mm)	Image no.	Target filter	kV	mAs	OD	Deviation	
20							
40							
45							
60							
70							
Contact mode				Large cassette			
(NB: If AEC sensor is common to both image receptor sizes then 45 mm phantom thickness only is required)							
Density control				AEC sensor position			
PMMA thickness (mm)	Image no.	Target filter	kV	mAs	OD	Deviation	
20							
40							
45							
60							
70							
Magnification mode				Small cassette			
Density control				AEC Sensor position			
PMMA thickness (mm)	Image no.	Target filter	kV	mAs	OD	Deviation	
20							
40							
45							
60							
70							
Tolerance							
Acceptable: Maximum deviation from $OD_{target} \leq 0.2$							
Achievable: Maximum deviation from $OD_{target} \leq 0.15$							
	Small cassette	Large cassette	Magnification				
Acceptable							
Achievable							

Comments:

7.4 Collimation system				
Focus-film distance (FFD):			mm	
Coincidence between X ray field and light field: (Note: +ve deviation indicates X ray field extends beyond light field)				
Anode material				
Collimator (cm x cm)				
Left edge deviation (mm)				
Deviation as % of FFD				
Right edge deviation (mm)				
Deviation as % of FFD				
Nipple edge deviation (mm)				
Deviation as % of FFD				
Chest wall edge deviation (mm)				
Deviation as % of FFD				
Tolerance				
Achievable: $\leq 1\%$ of FFD on any side				
Pass (Y/N)?				
Light field alignment (Pass (Y/N)?)				
Coincidence between X ray field and image receptor: (Note: +ve deviation indicates X ray field extends beyond image receptor)				
Left edge deviation (mm)				
% of FFD (retain sign)				
Right edge deviation (mm)				
% of FFD (retain sign)				
Anterior edge deviation (mm)				
% of FFD (retain sign)				
Chest edge deviation (mm)				
Tolerance				
Acceptable: Chest edge 0 to ≤ 5 mm; $\leq 2\%$ of FFD for other three edges Achievable: Chest edge 0 to ≤ 5 mm; > 0 for other three edges				
Pass?	Chest edge			
	Other edges			
X-ray field alignment (Pass (Y/N)?)				
Alignment of chest wall edge of compression paddle and image receptor: (Note: +ve deviation indicates paddle extends beyond image receptor)				
Diff. paddle edge & film (mm)				
Differences as % of FFD (mm)				
Tolerance				
Acceptable: Not visible at the chest wall & between 0 to $\leq 1\%$ of FFD				
Pass (Y/N)?				
Compression paddle alignment (Pass (Y/N)?)				
Comments:				

7.5 Image viewing conditions						
Viewbox information		Location				
Adjacent viewbox panel identification						
Fluorescent tube age						
Fluorescent tube type						
Dirt & marks present						
Panel colour difference?						
Functioning masks						
Luminance homogeneity						
Top left corner						
Bottom left corner						
Top right corner						
Bottom right corner						
Centre (L_c)						
Maximum deviation in luminance (%)	$100 \frac{L_{disc} - L_c}{L_c}$					
Average viewbox central luminance for all panels $L_{mean} =$						
Maximum deviation between panels in viewbox or adjacent panels (%)	$100 \frac{L_c - L_{mean}}{L_{mean}}$					
Tolerances						
Luminance uniformity within each panel in viewbox < 30%						
Pass (Y/N)?						
Acceptable (Y/N)?						
Viewbox central luminance for each panel > 3000 cd/m ² (nit)						
Pass (Y/N)?						
Acceptable (Y/N)?						
Luminance uniformity between panels in viewbox <15%						
Acceptable (Y/N)?						
Ambient illumination (viewboxes turned off)						
Viewbox location						
Illumination (lux)						
Tolerance						
Luminance uniformity		Achievable: ≤10 lux; Acceptable: ≤50 lux				
Acceptable (Y/N)?						
Achievable (Y/N)?						
Acceptable (Y/N)?						

Comments:

7.6.1 Evaluation of image quality			
Phantom used:			
Cassette size:			
Cassette ID:			
	Previous Film	Current Film	Change
Date			
Operation mode (manual, semiaut. automatic)			
Target/filter			
kVp setting			
AEC detector position:			
Density control			
OD _{target} (1.5-1.9)			
mAs			
Background OD (central location)			
OD_{target} - Background OD			
OD outside disc (OD ₁)			
OD inside disc (OD ₂)			
OD₁ - OD₂ (contrast)			
Number of fibres seen			
Number of fiber-like artefacts			
Fibres seen after deduction			
Number of micro-calc. groups seen			
Number of micro-like artefacts			
Speck groups after deduction			
Number of masses seen			
Number of masses-like artefacts			
Masses seen after deduction			
Tolerance			
Contrast: Achievable: ≥ 0.55; Acceptable: ≥ 0.40			
Acceptable (Y/N)?			
Achievable (Y/N)?			
Fibres: ≥ 4; Micro-calcifications: ≥ 3; Masses: ≥ 3			
Acceptable (Y/N)?			
Change in fibres, micro, masses ≤ 0.5			
Acceptable (Y/N)?			
Background optical density: OD_{target} ± 0.20			
Acceptable (Y/N)?			

7.6.2 System spatial resolution

X ray tube manufacturer		Model no.			
Emulsion no. (if the film box is new)					
Nominal focal spot size, f_{nom}					
Target material					
Nominal kVp setting					
Nominal mA setting					
Density control setting					
PMMA (mm)					
mAs					
Magnification factor	Contact				
Limiting resolution	bars // to A-C axis				
	bars \perp to A-C axis				
Tolerance					
Achievable: ≥ 15 lp/mm in both directions; Acceptable: ≥ 11 lp/mm in both directions					
bars // to A-C axis	Acceptable				
	Achievable				
bars \perp to A-C axis	Acceptable				
	Achievable				
Acceptable (Y/N)?					
Achievable (Y/N)?					

Comments:

7.7.1 Incident kerma at entrance surface of the phantom

Imaging mode		SID (mm)	
Screen type		Phantom thickness (mm)	
Film type		Source-Bucky distance (mm)	
Cassette size (cm)		Source-phantom distance (mm)	
Collimator		Compression force applied	
Dosimetry system used		Energy correction factor (N_{mammo})	
Pressure and temperature correction factor (k_{TP})			1.00

Determination of the mAs value for AEC

Exposure factors (with AEC)			
Nominal kVp setting			
Target/filter			
AEC mode			
AEC density setting			
Systems with post-exposure mAs display			
mAs_{auto}			
Systems without post-exposure mAs display			
<i>In AEC mode</i>			
Dosimeter reading (K_{auto})			
<i>In manual mode</i>			
mAs_{manual}			
Dosimeter reading (K_{manual})			
mAs_{auto}^1			
Measurement of the incident kerma			
Select manual mode and the above technique factors (target filter, kVp) and mAs_{auto}			
Dosimeter reading (mGy) M_{auto}			
If it is not possible to select the exact mAs_{auto} value			
$mAs_1 > mAs_{\text{auto}}$			
Dosimeter reading (mGy) M_1			
$mAs_2 < mAs_{\text{auto}}$			
Dosimeter reading (mGy) M_2			
Interpolated M_{auto}^2			
Incident kerma (mGy) K_i			

$$^1 \quad mAs_{\text{auto}} = \frac{M_{\text{meas}}}{M_{\text{man}}} \times mAs_{\text{man}}$$

$$^2 \quad M_{\text{auto}} = 0.5 \times mAs_{\text{auto}} \left(\frac{M_1}{mAs_1} + \frac{M_2}{mAs_2} \right)$$

7.7.2 Mean glandular dose (D_G)

Incident kerma (mGy) K_i		
Target/filter		
Measured HVL (mm Al):		
Product of g_{53} and c_{53}		
s factor		
D_G (mGy)		
Tolerance		
Achievable: $D_G \leq 2.0$ mGy; Acceptable: $D_G \leq 2.5$ mGy		
Acceptable (Y/N)?		
Achievable (Y/N)?		

HVL	Product of g_{53} and c_{53}
0.30	0.172
0.35	0.196
0.40	0.218
0.45	0.242
0.50	0.269
0.55	0.297
0.60	0.321

Target/filter	s factors
Mo/Mo	1.000
Mo/Rh	1.017
Rh/Al	1.044
Rh/Rh	1.061
W/Rh	1.042

Comments:

GLOSSARY

absorbed dose. The amount of energy per unit mass deposited in a medium. The SI unit of absorbed dose is the gray (Gy) representing 1 J/kg.

air kerma. The energy deposited per unit mass in air. The unit used to measure air kerma is the gray (Gy). For X rays with energies less than 300 kiloelectronvolts (keV), the magnitude of air kerma and absorbed dose in air are equivalent.

artefact. Any structure or pattern visible in the image that is not part of the object being imaged.

automatic exposure control (AEC) systems.

Fully automatic AEC. A device designed to determine the spectrum (target material, filtration material and kVp) and/or the exposure (mAs) needed to produce an adequately penetrated X ray image. This is typically done by sampling the X ray intensity after it passes through the patient and image receptor.

Automatic exposure time (phototiming). The beam quality parameters are selected manually and the exposure time is controlled automatically.

cassette. A light tight case usually made of thin, low X ray absorption plastic, that houses a film whose emulsion is, or whose emulsions are, in close contact with an intensifying screen(s).

coefficient of variation (COV). $COV = \sigma_x/\bar{x}$, where \bar{x} is the mean of a set of measurements x_i , and σ_x is their standard deviation.

compression device. A plastic paddle used to flatten and immobilize the breast. Compression helps reduce motion blurring in the breast, separates structures within the breast, and decreases the thickness of breast tissue. This minimizes the amount of radiation required and the amount of scattered radiation reaching the image receptor. Ideally, the compression device is made of rigid, thin plastic and has a flat bottom surface that is parallel to the plane of the image receptor and with edges perpendicular to the plane of the image receptor to assist in moving breast tissue away from the chest wall and into the field of view.

control chart. A graphical means of displaying data in which the variable of interest is plotted on the vertical axis as a function of time on the horizontal axis. The control chart allows for easy and rapid review of the data to determine whether the process is within the desired control limits ('in control').

control limit. The upper and lower values indicating that the process is 'out of control' and requiring that corrective action be taken. It is prudent to immediately repeat the measurement to verify that the system is 'out of control' before taking corrective action. If the repeated measurement is 'out of control', then corrective action is required immediately (or in some cases within 30 days). Synonym for 'action limit'.

craniocaudal (CC) view. One of two routine views for mammography. The image receptor is placed caudal to (below) the breast and the vertical X ray beam is directed from cranial to caudal (downward) through the breast.

D_{max} , D_{min} . The maximum and minimum optical density on a film. D_{max} is the darkest area of the film, where the highest exposure of X rays occurs. D_{min} is the base plus fog optical density of the film.

densitometer. An instrument for measuring the optical density or degree of blackening of film.

density difference (DD). The difference in optical density between the high density and the speed index. This provides an index of contrast in the mammography film.

detents. Mechanical settings that limit or prevent the motion, rotation or exposure of an X ray tube, cassette assembly or image receptor system.

diagnostic mammography. Mammography performed on patients who, by virtue of symptoms, physical findings or a prior screening mammography examination, are considered to have some likelihood of having breast cancer.

dose. See absorbed dose.

dosimeter. Equipment which uses ionization chambers and/or semiconductor detectors for the measurement of air kerma, in the beam of an X ray machine used for mammography examinations.

exposure. The act of initiating and producing X radiation from an X ray unit.

film reciprocity law failure. Refers to the loss of film speed as exposure time increases. This results in a loss of expected film density and is particularly noticeable in mammography for exposure times in excess of 2 s.

focal film distance (FFD). The distance from the focal spot to the film.

focal spot. The area of the target or anode that is bombarded by electrons from the cathode of the X ray tube to produce X rays. The smaller the focal spot, the better the limiting spatial resolution of the X ray system, especially in magnification mammography.

fog. The unwanted signal added to an image by the exposure of the image receptor to light, radiation or heat between patient exposures.

FSL. Fog produced on the film due to the safety light.

grid. A set of thin, closely spaced strips of highly attenuating material, such as lead, interspaced by a radiolucent support material, such as carbon fibre. In mammography, the grid is placed between the breast and image receptor to reduce scattered radiation reaching the image receptor. Scattered radiation reduces image contrast in mammography and limits the detection of low contrast structures such as fibres and masses. Grids improve the contrast of radiographic images at the price of increased dose to the patient.

half-value layer (HVL). The thickness of a specified substance which, when introduced into the path of a beam of radiation, reduces the exposure rate by one half. HVL is a measure of beam quality and is usually specified in millimetres of aluminium for diagnostic X ray equipment. The higher the HVL, the more penetrating the X ray beam.

illuminance. A photometric quantity describing the light intensity per unit area falling on a surface. The SI unit for illuminance is lux (lx), candela–steradians per square metre.

image receptor. A device that detects and records the distribution of X rays to form an image.

IOL. Initial operating level.

kilovoltage, peak (kVp). The maximum value of the potential difference (kV) between anode and cathode in an X ray tube.

linearity. System response where the output increases in direct proportion to the input signal.

luminance. A photometric quantity describing the light power per unit area per unit solid angle emitted by a light source. The SI unit for luminance is candelas per square metre (cd/m^2) (also known as 'nit').

mean glandular dose (MGD). The energy deposited per unit mass of fibroglandular tissue (the most radiosensitive tissue in the breast) averaged over all the fibroglandular tissue in the breast. The MGD is calculated from values of entrance air kerma, the X ray beam quality (half-value layer), and compressed breast thickness.

mediolateral oblique view (MLO). Now one of the standard two views of the breast. The image receptor is angled 30–60° from horizontal so that the cassette assembly is parallel to the pectoral muscle and the corner of the cassette holder fits comfortably into the axilla. The X ray beam is directed from the superomedial to the inferolateral aspect of the breast.

operating level. The central value about which we expect day to day measurements to fluctuate, for example, the empirically determined speed index on a sensitometric film.

phantom. A test object that simulates some aspect of human anatomy. A breast phantom simulates a typical breast in terms of size, composition and X ray attenuation, and may contain test objects that simulate anatomy in the breast.

PMMA. See *polymethyl methacrylate*.

polymethyl methacrylate (PMMA). Also known by the generic name acrylic, and trade names Plexiglas, Acrylate, Lucite and Perspex.

processor artefact. Any unwanted or artificial image feature appearing on an image due to a malfunction in the imaging chain or of the film processor.

resolution. The ability to determine small discrete objects with an imaging system. This may be measured using a bar pattern and is frequently described in terms of the number of line pairs visible per millimetre.

resolution pattern. A tool for determining the limiting spatial resolution of an imaging system. It is composed of groups of highly X ray attenuating strips spaced by an equal length of non-attenuating material. Each group of strips is smaller than the previous. The last group in which each bar is visible as a distinct line indicates the limiting spatial resolution of the system.

screen film mammography. Radiographic images of the breast performed with high detail intensifying screen(s) that are in close contact with the film in the cassette.

screening mammography. X ray breast examination of asymptomatic women in an attempt to detect breast cancer when it is small, non-palpable and confined to the breast.

speed index (S). In film sensitometry, the optical density of a specified step on a processed sensitometric image is chosen to reflect the operating point of the film.

standard breast. A 53 mm thick compressed breast consisting of 29% fibroglandular and 71% adipose tissue used as the ‘average’ breast for dosimetry calculations. This may be represented by 45 mm ± 0.5 mm of PMMA, a thickness which attenuates approximately the same amount as the standard breast. Rectangular with dimensions ≥150 mm × 100 mm or semicircular with radius ≥100 mm.

target optical density (OD_{target}). The optical density that the unit is aiming for. It should be selected by taking into account local factors (such as the type of film in use) and should be in the range of 1.5–1.9 OD.

TLD. A radiation measurement device typically comprising a chip or powder material that absorbs radiation and subsequently, upon heating, produces light whose intensity is proportional to the amount of radiation absorbed.

tolerance values. The range over which the parameter is allowed to vary before the item is no longer considered to be operating within limits. These ranges are classified into two categories: 'achievable' and 'acceptable'.

viewbox. A device providing a relatively uniform surface luminance for viewing mammographic films. Mammographic viewboxes should have a luminance level of at least 3000 cd/m² (nit).

CONTRIBUTORS TO DRAFTING AND REVIEW

Brandan, M.-E.	Instituto de Fisica, UNAM, Mexico
Burch, A.	Breast Test Wales, United Kingdom
Chevalier, M.	Complutense University of Madrid, Spain
Heggie, J.	St Vincents Hospital, Australia
Kanal, K.	University of Washington, United States of America
McLean, I.D.	International Atomic Energy Agency
Mora, P.	University of Costa Rica, Costa Rica
Morán, P.	Complutense University of Madrid, Spain
Ortiz López, P.	International Atomic Energy Agency
Pernička, F.	National Radiation Protection Institute, Czech Republic
Thomson, F.	Auckland Hospital, New Zealand
Yaffe, M.	Sunnybrook Health Sciences Centre, University of Toronto, Canada
Young, K.	Royal Surrey County Hospital, United Kingdom

Consultants Meetings

Vienna, Austria: 5–9 December 2005, 11–16 December 2006



Where to order IAEA publications

In the following countries IAEA publications may be purchased from the sources listed below, or from major local booksellers. Payment may be made in local currency or with UNESCO coupons.

Australia

DA Information Services, 648 Whitehorse Road, Mitcham Victoria 3132
Telephone: +61 3 9210 7777 • Fax: +61 3 9210 7788
Email: service@dadirect.com.au • Web site: <http://www.dadirect.com.au>

Belgium

Jean de Lannoy, avenue du Roi 202, B-1190 Brussels
Telephone: +32 2 538 43 08 • Fax: +32 2 538 08 41
Email: jean.de.lannoy@infoboard.be • Web site: <http://www.jean-de-lannoy.be>

Canada

Bernan Associates, 4611-F Assembly Drive, Lanham, MD 20706-4391, USA
Telephone: 1-800-865-3457 • Fax: 1-800-865-3450
Email: order@bernan.com • Web site: <http://www.bernan.com>

Renouf Publishing Company Ltd., 1-5369 Canotek Rd., Ottawa, Ontario, K1J 9J3
Telephone: +613 745 2665 • Fax: +613 745 7660
Email: order.dept@renoufbooks.com • Web site: <http://www.renoufbooks.com>

China

IAEA Publications in Chinese: China Nuclear Energy Industry Corporation, Translation Section, P.O. Box 2103, Beijing

Czech Republic

Suweco CZ, S.R.O. Klecakova 347, 180 21 Praha 9
Telephone: +420 26603 5364 • Fax: +420 28482 1646
Email: nakup@suweco.cz • Web site: <http://www.suweco.cz>

Finland

Akateeminen Kirjakauppa, PL 128 (Keskuskatu 1), FIN-00101 Helsinki
Telephone: +358 9 121 41 • Fax: +358 9 121 4450
Email: akatilaus@akateeminen.com • Web site: <http://www.akateeminen.com>

France

Form-Edit, 5, rue Janssen, P.O. Box 25, F-75921 Paris Cedex 19
Telephone: +33 1 42 01 49 49 • Fax: +33 1 42 01 90 90 • Email: formedit@formedit.fr

Lavoisier SAS, 14 rue de Provigny, 94236 Cachan Cedex
Telephone: + 33 1 47 40 67 00 • Fax +33 1 47 40 67 02
Email: livres@lavoisier.fr • Web site: <http://www.lavoisier.fr>

Germany

UNO-Verlag, Vertriebs- und Verlags GmbH, August-Bebel-Allee 6, D-53175 Bonn
Telephone: +49 02 28 949 02-0 • Fax: +49 02 28 949 02-22
Email: info@uno-verlag.de • Web site: <http://www.uno-verlag.de>

Hungary

Librotrade Ltd., Book Import, P.O. Box 126, H-1656 Budapest
Telephone: +36 1 257 7777 • Fax: +36 1 257 7472 • Email: books@librotrade.hu

India

Allied Publishers Group, 1st Floor, Dubash House, 15, J. N. Heredia Marg, Ballard Estate, Mumbai 400 001,
Telephone: +91 22 22617926/27 • Fax: +91 22 22617928
Email: alliedpl@vsnl.com • Web site: <http://www.alliedpublishers.com>

Bookwell, 24/4800, Ansari Road, Darya Ganj, New Delhi 110002
Telephone: +91 11 23268786, +91 11 23257264 • Fax: +91 11 23281315
Email: bookwell@vsnl.net • Web site: <http://www.bookwellindia.com>

Italy

Libreria Scientifica Dott. Lucio di Biasio "AEIOU", Via Coronelli 6, I-20146 Milan
Telephone: +39 02 48 95 45 52 or 48 95 45 62 • Fax: +39 02 48 95 45 48

Japan

Maruzen Company, Ltd., 13-6 Nihonbashi, 3 chome, Chuo-ku, Tokyo 103-0027
Telephone: +81 3 3275 8582 • Fax: +81 3 3275 9072
Email: journal@maruzen.co.jp • Web site: <http://www.maruzen.co.jp>

Korea, Republic of

KINS Inc., Information Business Dept. Samho Bldg. 2nd Floor, 275-1 Yang Jae-dong SeoCho-G, Seoul 137-130
Telephone: +02 589 1740 • Fax: +02 589 1746
Email: sj8142@kins.co.kr • Web site: <http://www.kins.co.kr>

Netherlands

Martinus Nijhoff International, Koraalrood 50, P.O. Box 1853, 2700 CZ Zoetermeer
Telephone: +31 793 684 400 • Fax: +31 793 615 698 • Email: info@nijhoff.nl • Web site: <http://www.nijhoff.nl>

Swets and Zeitlinger b.v., P.O. Box 830, 2160 SZ Lisse
Telephone: +31 252 435 111 • Fax: +31 252 415 888 • Email: infoho@swets.nl • Web site: <http://www.swets.nl>

New Zealand

DA Information Services, 648 Whitehorse Road, MITCHAM 3132, Australia
Telephone: +61 3 9210 7777 • Fax: +61 3 9210 7788
Email: service@dadirect.com.au • Web site: <http://www.dadirect.com.au>

Slovenia

Cankarjeva Založba d.d., Kopitarjeva 2, SI-1512 Ljubljana
Telephone: +386 1 432 31 44 • Fax: +386 1 230 14 35
Email: import.books@cankarjeva-z.si • Web site: <http://www.cankarjeva-z.si/uvoz>

Spain

Díaz de Santos, S.A., c/ Juan Bravo, 3A, E-28006 Madrid
Telephone: +34 91 781 94 80 • Fax: +34 91 575 55 63 • Email: compras@diazdesantos.es
carmela@diazdesantos.es • barcelona@diazdesantos.es • julio@diazdesantos.es
Web site: <http://www.diazdesantos.es>

United Kingdom

The Stationery Office Ltd, International Sales Agency, PO Box 29, Norwich, NR3 1 GN
Telephone (orders): +44 870 600 5552 • (enquiries): +44 207 873 8372 • Fax: +44 207 873 8203
Email (orders): book.orders@tso.co.uk • (enquiries): book.enquiries@tso.co.uk • Web site: <http://www.tso.co.uk>

On-line orders:

DELTA Int. Book Wholesalers Ltd., 39 Alexandra Road, Addlestone, Surrey, KT15 2PQ
Email: info@profbooks.com • Web site: <http://www.profbooks.com>

Books on the Environment:

Earthprint Ltd., P.O. Box 119, Stevenage SG1 4TP
Telephone: +44 1438748111 • Fax: +44 1438748844
Email: orders@earthprint.com • Web site: <http://www.earthprint.com>

United Nations (UN)

Dept. 1004, Room DC2-0853, First Avenue at 46th Street, New York, N.Y. 10017, USA
Telephone: +800 253-9646 or +212 963-8302 • Fax: +212 963-3489
Email: publications@un.org • Web site: <http://www.un.org>

United States of America

Bernan Associates, 4611-F Assembly Drive, Lanham, MD 20706-4391
Telephone: 1-800-865-3457 • Fax: 1-800-865-3450
Email: order@bernan.com • Web site: <http://www.bernan.com>

Renouf Publishing Company Ltd., 812 Proctor Ave., Ogdensburg, NY, 13669
Telephone: +888 551 7470 (toll-free) • Fax: +888 568 8546 (toll-free)
Email: order.dept@renoufbooks.com • Web site: <http://www.renoufbooks.com>

Orders and requests for information may also be addressed directly to:

Sales and Promotion Unit, International Atomic Energy Agency

Vienna International Centre, PO Box 100, 1400 Vienna, Austria
Telephone: +43 1 2600 22529 (or 22530) • Fax: +43 1 2600 29302
Email: sales.publications@iaea.org • Web site: <http://www.iaea.org/books>

This manual is intended to provide a standardized framework for quality control (QC) for mammography which can be used in Member States. It is intended to provide practical tests, using the simplest test equipment possible, to help ensure high quality screen film mammography. In order to maximize the benefit to those undergoing mammography where resources may be limited, this framework is supplemented with additional instructional material including a detailed section on clinical image quality.

IAEA HUMAN HEALTH SERIES