ACCELERATOR MASS SPECTROMETRY PROGRAMME AT MUMBAI PELLETRON ACCELERATOR FACILITY

P. Surendran, A. Shrivastava, A. K. Gupta, J. P.Nair, M. L. Yadav, M. Hemalatha, H. Sparrow, K. Mahata,
R. G. Thomas, P. V. Bhagwat, S. Kailas, Suman Sharma, A. S. Deodhar, U. Saravana Kumar & R. M. Kale

> AMS Programme based on <sup>36</sup> Cl - Environment and water samples

#### **INTRODUCTION**

#### AMS STABLE ISOTOPES OF VERY LOW ABUNDANCE AND LONG LIVED RADIOACTIVE ATOMS IN LOW CONCENTRATIONS


FROM 1977NEARLY60AMS FACILITIES


WHY AMS : Conventional MS

**Interference from molecular ions, isobars** 

**Radioactive decay counting** 


**Too long a counting time for small sample** 





#### **BEAM CHOPPER**

- A Chopper is designed to reduce beam current by a factor of 20 while transporting <sup>35</sup> Cl and <sup>37</sup> Cl.
- Chopper will be off for <sup>36</sup>Cl.
- **Duty Cycle = 1/20**
- Beam on time =  $50 \mu S$
- Beam off time =  $950 \mu S$




## **MULTI CATHODE ION SOURCE**

- A multi cathode sputter ion source
- 30 samples mounting cathode wheel undergoing testing at our ion source test bench set up.
- <u>The design and development of this source is</u> <u>fully indigenous in all respects.</u>
- Initial test results are encouraging. It is at an advanced stage of completion

# Important Applications of <sup>36</sup>Cl

- •
- measure the ages of ground water
- measure terrestrial ages of meteorites
- trace the movement of ground water
- trace the leakage of nuclear waste
- calibrate the cosmic ray flux using ice
- cores
- measure the neutron flux from the
- nuclear bombs
- dating of exposure ages of rocks




Accelerator based Mass Spectrometry (AMS) : ultra-sensitive means of counting individual atoms of long half life.

**Background :** Stable isotope of sulfur – <sup>36</sup>S


Detector : Multi-anode gas detector.
Length of ∆E anodes, matched according to the energy loss curves of <sup>36</sup>Cl and <sup>36</sup>S

## Energy Loss Curves <sup>36</sup>Cl, <sup>36</sup>S




In  $\Delta E1$  &  $\Delta E2$ , signal height of  ${}^{36}Cl > {}^{36}S$ In  $\Delta E3$ , Eres, signal height of  ${}^{36}S > {}^{36}Cl$  ${}^{36}S$  background can be reduced by Factor 10<sup>4</sup>


## **Hybrid Ionisation Chamber**



 $V_{anode} = +150 V$   $V_{grid} = +63 V$   $V_{cathode} = -150 V$ Gas : P10 ( 52 torr flow mode) Window : Mylar (3.5 $\mu$ m) Si -Detector : 300 $\mu$ m

**Fig.** 1





#### **Experiment**

To test and calibrate the hybrid Detector

<sup>32</sup>S, <sup>34</sup>S, <sup>36</sup>S, <sup>35</sup>Cl, <sup>37</sup>Cl, <sup>36</sup>Cl ions from 14UD Pelletron in low intensity, 10<sup>3</sup> pps were <u>detected directly</u> in the hybrid-gas detector

#### **Procedure Followed**

Tune beam – Maximise current in FC in front of detector Ensure that <sup>35</sup>Cl / <sup>37</sup>Cl ratio is as per natural abundance Optimise with chopper on and chopper off

Adjust IM – select isotopes of interest Keep AM field, quadrupoles ...all magnetic elements fixed Terminal voltage scaled As per mass - constant ME/q<sup>2</sup>

Calibrate detector with different ions - directly into the detector - Reduce filament current to do this


Detector optimised – signals from segmented detector <sup>32,34,36</sup>S and <sup>35,37</sup> Cl ions used

First measurement with BLANK. Later with STANDARD Finally unknown samples


#### **Table: Energy loss at different stages**


ion	VT MV	Ebeam MeV	Eres(after Mylar)MeV	A1 MeV	A2 MeV	A3 MeV
32S	7.873	62.99	55.6	15.3±0.5	17.0±0.8	9.7±0.4
34S	7.410	59.28	51.7	15.8±0.6	17.7±0.7	9.3±0.4
36S	6.999	55.99	48.0	16.5±0.6	17.8±0.7	7.9±0.4
<sup>35</sup> Cl	7.198	57.59	48.3	17.2±0.6	18.5±0.6	7.7±0.4
<sup>36</sup> Cl	6.998	55.99	47.5	17.4±0.8	18.4±0.7	7.1±0.4
<sup>37</sup> Cl	6.810	54.48	45.8	17.7±0.7	18.2±0.5	6.2±0.3

## Energy loss curves for <sup>32,34,36</sup>S and <sup>35,36,37</sup>Cl



### **Results**





## **RESULTS Old Water Samples**

Sample	Ratio <sup>36</sup> Cl/Cl	Sample	Ratio <sup>36</sup> Cl/Cl
Blank sample	~ 7×10 <sup>-14</sup>	Ground Water-2	4.28±0.27 ×10 <sup>-12</sup>
Standard Sample	<b>4.16×10</b> <sup>-11</sup>	Ground Water-3	5.00±0.50 ×10 <sup>-12</sup>
Ground Water-1	3.83±0.40 ×10 <sup>-12</sup>	Ground Water-4	2.80±0.27 ×10 <sup>-12</sup>

### CONCLUSION

AMS facility is ready for regular <sup>36</sup>Cl/ Cl ratio measurements

<sup>36</sup>Cl / Cl ratio of water samples range between 2 to 5 x 10<sup>-12</sup>

Controlled experiments planned to get the age of these Water samples ( expected to be over 35,000 years – 1- 4 pMC)

<sup>129</sup>I is planned for future

<sup>14</sup> C at IOP

<sup>10</sup> Be at IUAC